

Laboratory Chemicals

Price List - 2024



Scan to download

Part-I

Product Range

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---|------------|----------|--------|-----|
| 20287SG100 | 441-38-3 | α-Benzoin Oxime AR | 2942.00.90 | 100 gm. | 2728 | 18% |
| 10000SG500 | 9000-01-5 | Acacia powder extrapure | 1301.20.00 | 500 gm. | 963 | 5% |
| 10003SG500 | 60-35-5 | Acetamide for synthesis | 2924.29.90 | 500 gm. | 1262 | 18% |
| 20012LM500 | 64-19-7 | Acetic acid glacial AR DRY SOLVENT | 2915.21.00 | 500 ml. | 811 | 18% |
| 70010LM500 | 64-19-7 | Acetic acid glacial AR/ACS | 2915.21.00 | 500 ml. | 655 | 18% |
| 70010LC250 | | | | 2.5 ltr. | 2121 | 18% |
| 70010LL025 | | | | 25 ltr. | 12061 | 18% |
| 10010LM500 | 64-19-7 | Acetic acid glacial extrapure | 2915.21.00 | 500 ml. | 457 | 18% |
| 10010LC250 | | | | 2.5 ltr. | 1635 | 18% |
| 10010LL025 | | | | 25 ltr. | 11641 | 18% |
| 10010LL200 | | | | 200 ltr. | 93030 | 18% |
| 10010LK200 | | | | 200 kg. | 92400 | 18% |
| 30010LC100 | 64-19-7 | Acetic acid glacial HPLC & Spectroscopy | 2915.21.00 | 1 ltr. | 1518 | 18% |
| 20021LM500 | 67-64-1 | Acetone AR DRY SOLVENT | 2914.11.00 | 500 ml. | 713 | 18% |
| 70020LC250 | 67-64-1 | Acetone AR/ACS | 2914.11.00 | 2.5 ltr. | 2209 | 18% |
| 70020LL025 | | | | 25 ltr. | 15988 | 18% |
| 70020LM500 | 67-64-1 | Acetone AR/ACS | 2914.11.00 | 500 ml. | 591 | 18% |
| 10020LM500 | 67-64-1 | Acetone extrapure | 2914.11.00 | 500 ml. | 482 | 18% |
| 10020LC250 | | | | 2.5 ltr. | 1936 | 18% |
| 10020LL025 | | | | 25 ltr. | 16800 | 18% |
| 10020GL025 | 67-64-1 | Acetone extrapure - GI Drum | 2914.11.00 | 25 ltr. | 20188 | 18% |
| 10020LL200 | 67-64-1 | Acetone extrapure - MS Drum | 2914.11.00 | 200 ltr. | 128800 | 18% |
| 10020PL200 | 67-64-1 | Acetone extrapure - Plastic Drum | 2914.11.00 | 200 ltr. | 133000 | 18% |
| 66020LM250 | 67-64-1 | Acetone GC Grade | 2914.11.00 | 250 ml. | 1276 | 18% |
| 30020LC100 | 67-64-1 | Acetone HPLC & Spectroscopy | 2914.11.00 | 1 ltr. | 1532 | 18% |
| 30020LC250 | | | | 2.5 ltr. | 3523 | 18% |
| 20031LM500 | 75-5-8 | Acetonitrile AR DRY SOLVENT | 2926.90.00 | 500 ml. | 1035 | 18% |
| 20031LC250 | | | | 2.5 ltr. | 3015 | 18% |
| 70030LM500 | 75-05-8 | Acetonitrile AR/ACS | 2926.90.00 | 500 ml. | 760 | 18% |
| 70030LC250 | | | | 2.5 ltr. | 3032 | 18% |
| 70030LL025 | | | | 25 ltr. | 25083 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---|------------|----------|--------|-----|
| 10030LM500 | 75-05-8 | Acetonitrile extrapure | 2926.90.00 | 500 ml. | 720 | 18% |
| 10030LC250 | | | | 2.5 ltr. | 2972 | 18% |
| 10030LL025 | | | | 25 ltr. | 25375 | 18% |
| 10030LL200 | | | | 200 ltr. | 179200 | 18% |
| 66030LM250 | 75-05-8 | Acetonitrile GC Grade | 2926.90.00 | 250 ml. | 1366 | 18% |
| 40030LC250 | 75-05-8 | Acetonitrile Gradient Grade | 2926.90.01 | 2.5 ltr. | 3459 | 18% |
| 40030LC400 | | | | 4 lt. | 6371 | 18% |
| 30030LC250 | 75-05-8 | Acetonitrile HPLC & Spectroscopy | 2926.90.01 | 2.5 ltr. | 3400 | 18% |
| 10032LM500 | 98-86-2 | Acetophenone extrapure | 2914.39.10 | 500 ml. | 1212 | 18% |
| 10032LC250 | | | | 2.5 ltr. | 3876 | 18% |
| 10028LM500 | 75-36-5 | Acetyl chloride for synthesis | 2915.90.10 | 500 ml. | 905 | 18% |
| 10044S5K50 | 7440-44-0 | Activated charcoal granular (2.0 - 5.0 mm) (5X10) | 3802.10.00 | 50 kg. | 75572 | 18% |
| 10044SG500 | 7440-44-0 | Activated charcoal granular (2.0 - 5.0 mm) | 3802.10.00 | 500 gm. | 833 | 18% |
| 10044SR500 | 7440-44-0 | Activated charcoal granular | 3802.10.00 | 5 kg. | 6102 | 18% |
| 10040S5K25 | 7440-44-0 | Activated charcoal powder extrapure (5X5) | 3802.10.00 | 25 kg. | 29715 | 18% |
| 10043SG500 | 124-04-9 | Adipic acid extrapure | 2917.12.00 | 500 gm. | 994 | 18% |
| 10050SG250 | 9002-18-0 | Agar agar (for bacteriology) | 1302.31.00 | 250 gm. | 3214 | 18% |
| 10050SG500 | | | | 500 gm. | 6326 | 18% |
| 10054SG025 | 130-22-3 | Alizarin red S for microscopy | 2817.00.10 | 25 gm. | 1123 | 18% |
| 10060SG500 | 7784-26-1 | Aluminium ammonium sulphate dodecahydrate extrapure | 2833.22.90 | 500 gm. | 491 | 18% |
| 10062SG500 | 7446-70-0 | Aluminium chloride anhydrous extrapure | 2827.32.00 | 500 gm. | 773 | 18% |
| 20063S5K25 | 7784-13-6 | Aluminium chloride hexahydrate AR (5X5) | | 25 kg. | 46690 | 18% |
| 20063SG500 | 7784-13-6 | Aluminium chloride hexahydrate AR | 2827.32.00 | 500 gm. | 3530 | 18% |
| 70064SG500 | 7784-27-2 | Aluminium nitrate nonahydrate AR/ACS | 2834.29.90 | 500 gm. | 683 | 18% |
| 10064SG500 | 7784-27-2 | Aluminium nitrate Nonahydrate extrapure | 2834.29.90 | 500 gm. | 610 | 18% |
| 10071SG500 | 1344-28-1 | Aluminium oxide acidic | 2818.20.00 | 500 gm. | 745 | 18% |
| 10072S5K50 | 1344-28-1 | Aluminium oxide basic (5X10) | 2818.20.00 | 50 kg. | 60235 | 18% |
| 10072SG500 | 1344-28-1 | Aluminium oxide basic | 2818.20.00 | 500 gm. | 735 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 10073S5K50 | 1344-28-1 | Aluminium oxide neutral (5X10) | 2818.20.00 | 50 kg. | 56739 | 18% |
| 10073SG500 | 1344-28-1 | Aluminium oxide neutral | 2818.20.00 | 500 gm. | 665 | 18% |
| 20065SG500 | 7784-24-9 | Aluminium potassium sulphate dodecahydrate AR | 2833.22.90 | 500 gm. | 509 | 18% |
| 20065SR500 | | | | 5 kg. | 4865 | 18% |
| 10065SG500 | 7784-24-9 | Aluminium potassium sulphate dodecahydrate extrapure | 2833.22.90 | 500 gm. | 381 | 18% |
| 12067SK050 | 16828-11-8 | Aluminium sulphate Battery grade | 2833.22.10 | 50 kg. | 22127 | 18% |
| 10067SG500 | 16828-11-8 | Aluminium sulphate extrapure | 2833.22.10 | 500 gm. | 356 | 18% |
| 10066SG025 | 83-07-8 | 4-Aminoantipyrine extrapure | 2933.11.00 | 25 gm. | 1514 | 18% |
| 10066SG100 | | | | 100 gm. | 5883 | 18% |
| 10068SG100 | 150-13-0 | 4-Aminobenzoic acid extrapure | 2922.29.90 | 100 gm. | 1170 | 18% |
| 10068SG500 | | | | 500 gm. | 4227 | 18% |
| 20075SG025 | 116-63-2 | 1-Amino-2-naphthol-4-sulphonic acid AR | 2922.21.90 | 25 gm. | 815 | 18% |
| 20075SG100 | | | | 100 gm. | 2466 | 18% |
| 10075SG025 | 116-63-2 | 1-Amino-2-naphthol-4-sulphonic acid extrapure | | 25 gm. | 680 | 18% |
| 10069LM500 | 1336-21-6 | Ammonia buffer solution | 3822.00.90 | 500 ml. | 313 | 12% |
| 20070LM500 | 1336-21-6 | Ammonia solution 0.9ld AR-Glass bottle | 2814.20.00 | 500 ml. | 302 | 18% |
| 20070LC250 | | | | 2.5 ltr. | 902 | 18% |
| 10070LM500 | 1336-21-6 | Ammonia solution 0.9ld extrapure | 2814.20.00 | 500 ml. | 278 | 18% |
| 10070LC250 | | | | 2.5 ltr. | 798 | 18% |
| 10070LC500 | | | | 5 ltr. | 1131 | 18% |
| 10070LL025 | | | | 25 ltr. | 7525 | 18% |
| 10070TP250 | 1336-21-6 | Ammonia solution 0.9ld thermocol packing extrapure | | 2.5 ltr. | 2330 | 18% |
| 70080SG500 | 631-61-8 | Ammonium acetate AR/ACS | 2915.29.90 | 500 gm. | 719 | 18% |
| 10080SG500 | 631-61-8 | Ammonium acetate extrapure | 2915.29.90 | 500 gm. | 588 | 18% |
| 30080SG500 | 631-61-8 | Ammonium acetate HPLC | 2915.29.90 | 500 gm. | 1385 | 18% |
| 10082S5K25 | 1066-33-7 | Ammonium bicarbonate extrapure (5X5) | 2836.99.90 | 25 kg. | 9557 | 18% |
| 10082SG500 | 1066-33-7 | Ammonium bicarbonate extrapure | 2836.99.90 | 500 gm. | 341 | 18% |
| 10085S5K50 | 1341-49-7 | Ammonium bifluoride extrapure (5X10) | 2827.39.90 | 50 kg. | 92750 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|---------|-------|-----|
| 10085SG500 | 1341-49-7 | Ammonium bifluoride extrapure | 2827.39.90 | 500 gm. | 877 | 18% |
| 10087SG500 | 12124-97-9 | Ammonium bromide extrapure | 2827.59.90 | 500 gm. | 1463 | 18% |
| 20090SG500 | 10361-29-2 | Ammonium carbonate AR | 2836.20.90 | 500 gm. | 751 | 18% |
| 10090SG500 | 10361-29-2 | Ammonium carbonate extrapure | 2836.20.90 | 500 gm. | 569 | 18% |
| 70092SG100 | 16774-21-3 | Ammonium ceric nitrate AR/ACS | 3102.60.00 | 100 gm. | 3762 | 18% |
| 10092SG100 | 16774-21-3 | Ammonium ceric nitrate extrapure | 3102.60.00 | 100 gm. | 2932 | 18% |
| 70094SG100 | 10378-47-9 | Ammonium ceric sulphate AR/ACS | 3102.29.90 | 100 gm. | 3800 | 18% |
| 10094SG100 | 10378-47-9 | Ammonium ceric sulphate extrapure | 3102.29.90 | 100 gm. | 3747 | 18% |
| 70095SK025 | 12125-2-9 | Ammonium chloride AR/ACS | 2827.10.00 | 25 kg. | 17192 | 18% |
| 70095SG500 | 12125-02-9 | Ammonium chloride AR/ACS | 2827.10.00 | 500 gm. | 637 | 18% |
| 10095S5K50 | 12125-02-9 | Ammonium chloride extrapure (5X10) | 2827.10.00 | 50 kg. | 33205 | 18% |
| 10095SG500 | 12125-2-9 | Ammonium chloride extrapure | 2827.10.00 | 500 gm. | 463 | 18% |
| 10095SR500 | | | | 5 kg. | 2993 | 18% |
| 10095SK025 | | | | 25 kg. | 18235 | 18% |
| 20118SG500 | 7789-09-5 | Ammonium dichromate AR | 2841.50.90 | 500 gm. | 1515 | 18% |
| 10118SG500 | 7789-09-5 | Ammonium dichromate extrapure | 2841.50.90 | 500 gm. | 1482 | 18% |
| 70120SG500 | 7722-76-1 | Ammonium di-hydrogen ortho-phosphate AR/ACS | 3105.40.00 | 500 gm. | 951 | 18% |
| 10120S5K50 | 7722-76-1 | Ammonium di-hydrogen ortho-phosphate extrapure (5X10) | 3105.40.00 | 50 kg. | 89135 | 18% |
| 10120SG500 | 7722-76-1 | Ammonium di-hydrogen ortho-phosphate extrapure | 3105.40.00 | 500 gm. | 846 | 18% |
| 30120SG500 | 7722-76-1 | Ammonium di-hydrogen ortho-phosphate HPLC | 3105.40.00 | 500 gm. | 2938 | 18% |
| 70130SG500 | 7783-85-9 | Ammonium ferrous sulphate hexahydrate AR/ACS | 2833.40.00 | 500 gm. | 651 | 18% |
| 10130S5K50 | 7783-85-9 | Ammonium ferrous sulphate hexahydrate extrapure (5X10) | | 50 kg. | 28560 | 18% |
| 10130SG500 | 7783-85-9 | Ammonium ferrous sulphate hexahydrate extrapure | 2833.40.00 | 500 gm. | 408 | 18% |
| 10133SG500 | 12125-01-8 | Ammonium fluoride extrapure | 2826.11.00 | 500 gm. | 1306 | 18% |
| 10136SG500 | 540-69-2 | Ammonium formate extrapure | 2915.12.10 | 500 gm. | 724 | 18% |
| 20145SG100 | 7803-55-6 | Ammonium metavanadate AR | 2825.30.90 | 100 gm. | 2347 | 18% |
| 20145SG500 | | | | 500 gm. | 9445 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|----------|-------|-----|
| 10145SG100 | 7803-55-6 | Ammonium metavanadate extrapure | 2825.30.90 | 100 gm. | 1715 | 18% |
| 10145SG500 | | | | 500 gm. | 8127 | 18% |
| 70150SG100 | 12054-85-2 | Ammonium molybdate tetrahydrate AR/ACS | 2841.70.90 | 100 gm. | 2533 | 18% |
| 70150SG250 | | | | 250 gm. | 4620 | 18% |
| 70150SG500 | | | | 500 gm. | 10436 | 18% |
| 10150SG100 | 12054-85-2 | Ammonium molybdate tetrahydrate extrapure | 2841.70.90 | 100 gm. | 2029 | 18% |
| 10150SG250 | | | | 250 gm. | 4463 | 18% |
| 10150SG500 | | | | 500 gm. | 8594 | 18% |
| 70170SG500 | 6009-70-7 | Ammonium oxalate monohydrate AR/ACS | 2917.11.90 | 500 gm. | 800 | 18% |
| 10170SG500 | 6009-70-7 | Ammonium oxalate monohydrate extrapure | 2917.11.90 | 500 gm. | 595 | 18% |
| 70180SG500 | 7727-54-0 | Ammonium persulphate AR/ACS | 2833.40.00 | 500 gm. | 697 | 18% |
| 10180S5K50 | 7727-54-0 | Ammonium persulphate extrapure (5X10) | 2833.40.00 | 50 kg. | 39410 | 18% |
| 10180SG500 | 7727-54-0 | Ammonium persulphate extrapure | 2833.40.00 | 500 gm. | 608 | 18% |
| 10180SR500 | | | | 5 kg. | 4088 | 18% |
| 20185SG025 | 13573-16-5 | Ammonium reineckate 96% AR | 2842.90.90 | 25 gm.* | 3855 | 18% |
| 70190SG500 | 7783-20-2 | Ammonium sulphate AR/ACS | 3102.21.00 | 500 gm. | 418 | 18% |
| 10190SG500 | 7783-20-2 | Ammonium sulphate extrapure | 3102.21.00 | 500 gm. | 338 | 18% |
| 10190SR500 | | | | 5 kg. | 2038 | 18% |
| 10190S5K50 | | | | 50 kg. | 31150 | 18% |
| 70200SG500 | 1762-95-4 | Ammonium thiocyanate AR/ACS | 2833.30.10 | 500 gm. | 1164 | 18% |
| 10200SG500 | 1762-95-4 | Ammonium thiocyanate extrapure | 2833.30.10 | 500 gm. | 944 | 18% |
| 11005SG100 | 90-15-3 | α-Naphthol extrapure | 2907.15.10 | 100 gm. | 737 | 18% |
| 11005SG500 | | | | 500 gm. | 2334 | 18% |
| 70204LM500 | 62-53-3 | Aniline AR/ACS | 2921.41.10 | 500 ml. | 1043 | 18% |
| 10204LM500 | 62-53-3 | Aniline extrapure | 2921.41.10 | 500 ml. | 930 | 18% |
| 10204LC250 | | | | 2.5 ltr. | 3442 | 18% |
| 10212SG500 | 1309-64-4 | Antimony trioxide extrapure | 2825.80.00 | 500 gm. | 3114 | 18% |
| 10219SG500 | 1327-53-3 | Arsenic trioxide tech | 2811.29.10 | 500 gm. | 942 | 18% |
| 20221SG010 | 5941-07-1 | Azomethine H AR | 2925.29.90 | 10 gm. | 37081 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|--------|-----|
| 10222SG500 | 543-80-6 | Barium acetate extrapure | 2915.29.90 | 500 gm. | 1225 | 18% |
| 10224SG500 | 7791-28-8 | Barium bromide extrapure | 2827.59.90 | 500 gm. | 1901 | 18% |
| 10225SG500 | 513-77-9 | Barium carbonate extrapure | 2836.60.00 | 500 gm. | 873 | 18% |
| 10231LM500 | 10326-27-9 | Barium chloride 10% w/v solution | 2827.39.90 | 500 ml. | 370 | 18% |
| 70230SG500 | 10326-27-9 | Barium chloride dihydrate AR/ACS | 2827.39.90 | 500 gm. | 565 | 18% |
| 10230SG500 | 10326-27-9 | Barium chloride dihydrate extrapure | 2827.39.90 | 500 gm. | 502 | 18% |
| 10230SR500 | | | | 5 kg. | 3371 | 18% |
| 20233SG005 | 6211-24-1 | Barium diphenylamine sulphonate AR | 2921.11.10 | 5 gm. | 612 | 18% |
| 20233SG025 | | | | 25 gm. | 1799 | 18% |
| 70240SG500 | 10022-31-8 | Barium nitrate AR/ACS | 2834.29.30 | 500 gm. | 778 | 18% |
| 10240SG500 | 10022-31-8 | Barium nitrate extrapure | 2834.29.30 | 500 gm. | 657 | 18% |
| 10245S5K25 | 7727-43-7 | Barium sulphate (suitable for radiology) extrapure (5X5) | 2833.27.00 | 25 kg. | 44636 | 18% |
| 10245SG500 | 7727-43-7 | Barium sulphate (suitable for radiology) extrapure | 2833.27.00 | 500 gm. | 538 | 18% |
| 12245SK050 | 7727-43-7 | Barium sulphate Battery grade | 2833.27.00 | 50 kg. | 958608 | 18% |
| 10247LM500 | 63126-89-6 | Benedict's reagent (Qualitative) | 3822.00.90 | 500 ml. | 413 | 12% |
| 10247LC500 | | | | 5 ltr. | 2795 | 12% |
| 10248SG500 | 1302-78-9 | Bentonite powder | 2508.10.90 | 500 gm. | 742 | 5% |
| 20248S5K25 | 1302-78-9 | Bentonite powder AR (5X5) | 2509.10.90 | 25 kg. | 20090 | 5% |
| 10250LM500 | 100-52-7 | Benzaldehyde extrapure | 2912.21.00 | 500 ml. | 985 | 18% |
| 10250LC250 | | | | 2.5 ltr. | 3838 | 18% |
| 10261LC500 | 8001-54-5 | Benzalkonium chloride 5% soln | 2923.90.00 | 5 ltr. | 1928 | 18% |
| 10260LM500 | 8001-54-5 | Benzalkonium chloride 50% soln | 2923.90.00 | 500 ml. | 911 | 18% |
| 10260LC500 | | | | 5 ltr. | 4987 | 18% |
| 10265SG500 | 55-21-0 | Benzamide | 2924.29.90 | 500 gm. | 2179 | 18% |
| 20270LM500 | 71-43-2 | Benzene AR | 2902.20.00 | 500 ml. | 704 | 18% |
| 20270LC250 | | | | 2.5 ltr. | 2687 | 18% |
| 10270LM500 | 71-43-2 | Benzene extrapure | 2902.20.00 | 500 ml. | 561 | 18% |
| 10270LC250 | | | | 2.5 ltr. | 2061 | 18% |
| 10270LL025 | | | | 25 ltr. | 22078 | 18% |
| 30270LC100 | 71-43-2 | Benzene HPLC & Spectroscopy | 2902.20.00 | 1 ltr. | 1642 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---------------------------------------|------------|----------|-------|-----|
| 10275SG500 | 98-11-3 | Benzene sulphonic acid EP 500 GM | 2902.20.00 | 500 gm. | 3175 | 18% |
| 10280SG025 | 121-54-0 | Benzethonium chloride (Hyamine 1622) | 2923.10.00 | 25 gm. | 2709 | 18% |
| 10279LC100 | 121-54-0 | Benzethonium chloride 0.004M solution | 2923.10.00 | 1 ltr. | 2856 | 18% |
| 70285SG500 | 65-85-0 | Benzoic acid AR/ACS | 2916.31.10 | 500 gm. | 2674 | 18% |
| 10285S5K25 | 65-85-0 | Benzoic acid extrapure (5X5) | 2916.31.10 | 25 kg. | 44266 | 18% |
| 10285SG500 | 65-85-0 | Benzoic acid extrapure | | 500 gm. | 690 | 18% |
| 10282SG100 | 119-53-9 | Benzoin for synthesis | 2942.00.90 | 100 gm. | 989 | 18% |
| 10282SG250 | | | | 250 gm. | 2509 | 18% |
| 10284SG500 | 119-61-9 | Benzophenone, min 99% for synthesis | 2914.30.90 | 500 gm. | 1377 | 18% |
| 10288LM500 | 98-88-4 | Benzoyl chloride for synthesis | 2916.32.00 | 500 ml. | 1184 | 18% |
| 10288LC250 | | | | 2.5 ltr. | 4946 | 18% |
| 70290LM500 | 100-51-6 | Benzyl alcohol AR/ACS | 2906.21.00 | 500 ml. | 1154 | 18% |
| 70290LC250 | | | | 2.5 ltr. | 4451 | 18% |
| 10290LM500 | 100-51-6 | Benzyl alcohol extrapure | 2906.21.00 | 500 ml. | 1015 | 18% |
| 10290LC250 | | | | 2.5 ltr. | 3645 | 18% |
| 10290LL025 | | | | 25 ltr. | 72891 | 18% |
| 30290LC100 | 100-51-6 | Benzyl alcohol HPLC & Spectroscopy | 2906.21.00 | 1 ltr. | 5351 | 18% |
| 10293LM500 | 120-51-4 | Benzyl benzoate extrapure | 2833.29.90 | 500 ml. | 821 | 18% |
| 10295LM500 | 100-44-7 | Benzyl chloride extrapure | 2903.69.40 | 500 ml. | 708 | 18% |
| 10295LC250 | | | | 2.5 ltr. | 2916 | 18% |
| 10292LM500 | 100-46-9 | Benzylamine,min.98% for synthesis | 2918.99.00 | 500 ml. | 1910 | 18% |
| 10301SG100 | 592-10-4 | Bismuth carbonate extrapure | 2916.19.20 | 100 gm.* | 1539 | 18% |
| 10301SG500 | | | | 500 gm.* | 6841 | 18% |
| 10302SG100 | 10035-06-0 | Bismuth nitrate extrapure | 2916.19.20 | 100 gm. | 1321 | 18% |
| 10302SG500 | | | | 500 gm. | 5446 | 18% |
| 10304SG100 | 10361-46-3 | Bismuth subnitrate extrapure | 2916.19.20 | 100 gm. | 2725 | 18% |
| 10304SG500 | | | | 500 gm. | 13062 | 18% |
| 10306SG250 | 7787-68-0 | Bismuth sulphate | 2833.29.90 | 250 gm.* | 5559 | 18% |
| 11006SG500 | 135-19-3 | b-Naphthol extrapure | 2907.15.20 | 500 gm. | 1365 | 18% |
| 70310SG500 | 10043-35-3 | Boric acid crystals AR/ACS | 2810.00.20 | 500 gm. | 1225 | 18% |
| 70310SR500 | | | | 5 kg. | 8775 | 18% |
| 70310S5K50 | | | | 50 kg. | 56494 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|----------|-------|-----|
| 10310SG500 | 10043-35-3 | Boric acid crystals extrapure | 2810.00.20 | 500 gm. | 773 | 18% |
| 10310SR500 | | | | 5 kg. | 4851 | 18% |
| 10310S5K50 | | | | 50 kg. | 56494 | 18% |
| 10311SG500 | 10043-35-3 | Boric acid powder extrapure | 2810.00.20 | 500 gm. | 819 | 18% |
| 10311SR500 | | | | 5 kg. | 6051 | 18% |
| 10311S5K50 | | | | 50 kg. | 40250 | 18% |
| 10312LM500 | 109-63-7 | Boron trifluoride ethyl etherate Extrapure | 2942.00.90 | 500 ml. | 4841 | 18% |
| 40314LM125 | 81029-05-2 | Brilliant cresyl blue indicator solution | 3822.00.90 | 125 ml. | 619 | 12% |
| 10313LM020 | 7726-95-6 | Bromine ampules 5X20ML | 2801.30.20 | 5x20 ml. | 2311 | 18% |
| 20315SG005 | 76-60-8 | Bromocresol green AR | 2934.99.00 | 5 gm. | 1279 | 18% |
| 20315SG025 | | | | 25 gm. | 4472 | 18% |
| 40315LM125 | 76-60-8 | Bromocresol green indicator solution | 2934.99.00 | 125 ml. | 249 | 18% |
| 20316SG005 | 115-40-2 | Bromocresol purple AR | 2934.99.00 | 5 gm. | 402 | 18% |
| 20316SG025 | | | | 25 gm. | 1510 | 18% |
| 40316LM125 | 115-40-2 | Bromocresol purple indicator solution | 2934.99.00 | 125 ml. | 246 | 18% |
| 20318SG005 | 115-39-9 | Bromophenol blue AR | 2934.99.00 | 5 gm. | 418 | 18% |
| 20318SG025 | | | | 25 gm. | 1434 | 18% |
| 40318LM125 | 115-39-9 | Bromophenol blue indicator solution | 2934.99.00 | 125 ml. | 250 | 18% |
| 20320SG005 | 76-59-5 | Bromothymol blue AR | 2934.99.00 | 5 gm. | 479 | 18% |
| 20320SG025 | | | | 25 gm. | 1534 | 18% |
| 40320LM125 | 76-59-5 | Bromothymol blue indicator solution | 2934.99.00 | 125 ml. | 236 | 18% |
| 20323SG025 | 357-57-3 | Brucine AR | 2934.99.00 | 25 gm. | 2759 | 18% |
| 10325CP010 | NA | Buffer capsules pH 4.00±0.05 10CP | 3822.00.90 | 10 caps. | 323 | 12% |
| 10326CP010 | NA | Buffer capsules pH 7.00 ± 0.05 10CP | 3822.00.90 | 10 caps. | 323 | 12% |
| 10327CP010 | NA | Buffer capsules pH 9.20±0.05 10CP | 3822.00.90 | 10 caps. | 323 | 12% |
| 10336LM500 | NA | Buffer Solution pH 4.0 | 3822.00.90 | 500 ml. | 607 | 12% |
| 10337LM500 | NA | Buffer solution pH 7.0 | 3822.00.90 | 500 ml. | 607 | 12% |
| 10338LM500 | NA | Buffer solution pH 9.2 | 3822.00.90 | 500 ml. | 607 | 12% |
| 70330LM500 | 71-36-3 | Butan-1-ol AR/ACS (n-Butanol) | 2905.13.00 | 500 ml. | 660 | 18% |
| 70330LC250 | 71-36-3 | Butan-1-ol AR/ACS (n-Butanol) | 2905.13.00 | 2.5 ltr. | 2488 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|--------|-----|
| 10330LM500 | 71-36-3 | Butan-1-ol extrapure | 2905.13.00 | 500 ml. | 563 | 18% |
| 10330LC250 | | | | 2.5 ltr. | 2204 | 18% |
| 10330LL025 | | | | 25 ltr. | 21840 | 18% |
| 30330LC100 | 71-36-3 | Butan-1-ol HPLC & Spectroscopy | 2905.13.00 | 1 ltr. | 1700 | 18% |
| 20340LM500 | 78-93-3 | Butan-2-one AR | 2914.12.00 | 500 ml. | 1152 | 18% |
| 20340LC250 | | | | 2.5 ltr. | 4398 | 18% |
| 10340LM500 | 78-93-3 | Butan-2-one extrapure | 2914.12.00 | 500 ml. | 1119 | 18% |
| 10340LC250 | | | | 2.5 ltr. | 4012 | 18% |
| 10340LL025 | | | | 25 ltr. | 22108 | 18% |
| 10340LL200 | | | | 200 ltr. | 151326 | 18% |
| 30340LC100 | 78-93-3 | Butan-2-one HPLC & Spectroscopy | 2914.12.00 | 1 ltr. | 2307 | 18% |
| 30350SG025 | 2386-54-1 | 1-Butane sulphonic acid sodium salt anhydrous HPLC | 2904.10.90 | 25 gm. | 1519 | 18% |
| 32350SG025 | 2386-54-1 | 1-Butane sulphonic acid sodium salt monohydrous for HPLC | 2904.10.90 | 25 gm. | 1519 | 18% |
| 10375LM500 | 111-76-2 | Butyl cellosolve extrapure | 2909.43.00 | 500 ml. | 759 | 18% |
| 10375LC250 | | | | 2.5 ltr. | 2822 | 18% |
| 10375LL025 | | | | 25 ltr. | 26873 | 18% |
| 10370SG100 | 25013-16-5 | Butylated hydroxyanisole extrapure | 2909.30.19 | 100 gm. | 1724 | 18% |
| 10370SG500 | | | | 500 gm. | 4750 | 18% |
| 10380SG500 | 128-37-0 | Butylated hydroxytoluene extrapure | 2909.30.19 | 500 gm. | 1330 | 18% |
| 10380SR500 | | | | 5 kg. | 9829 | 18% |
| 10380S5K50 | | | | 50 kg. | 213679 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|---------|--------|-----|
| 20382SG100 | 35658-65-2 | Cadmium chloride monohydrate AR | 2827.39.90 | 100 gm. | 1259 | 18% |
| 20382SG500 | | | | 500 gm. | 4677 | 18% |
| 10382SG100 | 35658-65-2 | Cadmium chloride monohydrate extrapure | 2827.39.90 | 100 gm. | 942 | 18% |
| 10382SG500 | | | | 500 gm. | 4016 | 18% |
| 10388SG100 | 7790-80-9 | Cadmium iodide | 2827.60.00 | 100 gm. | 3385 | 18% |
| 20383SG100 | 10022-68-1 | Cadmium nitrate tetrahydrate AR | 2834.29.90 | 100 gm. | 741 | 18% |
| 10383SG100 | 10022-68-1 | Cadmium nitrate tetrahydrate extrapure | 2834.29.90 | 100 gm. | 713 | 18% |
| 10383SG500 | | | | 500 gm. | 1883 | 18% |
| 10384SG100 | 1306-19-0 | Cadmium oxide extrapure | 2825.90.20 | 100 gm. | 1088 | 18% |
| 20385SG100 | 7790-84-3 | Cadmium sulphate AR | 2833.29.90 | 100 gm. | 1217 | 18% |
| 10385SG100 | 7790-84-3 | Cadmium sulphate extrapure | 2833.29.90 | 100 gm. | 830 | 18% |
| 10385SG500 | | | | 500 gm. | 3062 | 18% |
| 10405S5K25 | 62-54-4 | Calcium acetate extrapure (5X5) | 2915.29.10 | 25 kg. | 49193 | 18% |
| 10405SG500 | 62-54-4 | Calcium acetate extrapure | 2915.29.10 | 500 gm. | 780 | 18% |
| 70410S5K50 | 471-34-1 | Calcium carbonate AR/ACS (5X10) | | 50 kg. | 159429 | 18% |
| 70410SG500 | 471-34-1 | Calcium carbonate AR/ACS | 2836.50.00 | 500 gm. | 666 | 18% |
| 10410S5K50 | 471-34-1 | Calcium carbonate extrapure (5X10) | 2836.50.00 | 50 kg. | 18813 | 18% |
| 10410SG500 | 471-34-1 | Calcium carbonate extrapure | 2836.50.00 | 500 gm. | 403 | 18% |
| 10410SR500 | | | | 5 kg. | 1833 | 18% |
| 70430SG500 | 10035-04-8 | Calcium chloride dihydrate AR/ACS | 2827.20.00 | 500 gm. | 1046 | 18% |
| 10430S5K50 | 10035-04-8 | Calcium chloride dihydrate extrapure (5X10) | | 50 kg. | 22904 | 18% |
| 10430SG500 | 10035-04-8 | Calcium chloride dihydrate extrapure | 2827.20.00 | 500 gm. | 408 | 18% |
| 10430SR500 | | | | 5 kg. | 2370 | 18% |
| 10420S1K25 | 10043-52-4 | Calcium chloride fused extrapure (1X25) | | 25 kg. | 8540 | 18% |
| 10420S5K25 | 10043-52-4 | Calcium chloride fused extrapure (5X5) | | 25 kg. | 8061 | 18% |
| 10420SG500 | 10043-52-4 | Calcium chloride fused extrapure | 2827.20.00 | 500 gm. | 331 | 18% |
| 10440S5K50 | 1305-62-0 | Calcium hydroxide extrapure (5X10) | 2825.90.40 | 50 kg. | 17997 | 18% |
| 10440SG500 | 1305-62-0 | Calcium hydroxide extrapure | 2825.90.40 | 500 gm. | 364 | 18% |
| 10441SG500 | 7778-54-3 | Calcium Hypochlorite | 2828.10.00 | 500 gm. | 323 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 10442S5K50 | 13477-34-4 | Calcium nitrate tetrahydrate extrapure (5X10) | 2834.29.90 | 50 kg. | 27367 | 18% |
| 10442SG500 | 13477-34-4 | Calcium nitrate tetrahydrate extrapure | 2834.29.90 | 500 gm. | 337 | 18% |
| 10443SG500 | 1305-78-8 | Calcium oxide powder extrapure | 2825.90.90 | 500 gm. | 452 | 18% |
| 10443SR500 | | | | 5 kg. | 3505 | 18% |
| 10448SG500 | 10101-41-4 | Calcium sulphate dihydrate extrapure | 2833.29.90 | 500 gm. | 460 | 18% |
| 10451SG025 | 4197-24-4 | Carbol fuchsin powder (M.S) for microscopy | 3204.13.43 | 25 gm. | 533 | 18% |
| 10451SG100 | | | | 100 gm. | 1648 | 18% |
| 10452LM125 | NA | Carbol fuchsin Stain Solution Dilute for microscopy | 3204.13.43 | 125 ml. | 228 | 18% |
| 10453LM125 | NA | Carbol fuchsin stain solution Strong for microscopy | 3204.13.43 | 125 ml. | 267 | 18% |
| 10453LM500 | | | | 500 ml. | 718 | 18% |
| 10450LM500 | 75-15-0 | Carbon disulphide extrapure | 2813.10.00 | 500 ml. | 1125 | 18% |
| 10469LM125 | 1390-65-4 | Carmine aceto staining solution for microscopy | 3204.19.90 | 125 ml. | 1813 | 18% |
| 10459SG005 | 1390-65-4 | Carmine for microscopy | 3204.19.90 | 5 gm. | 1348 | 18% |
| 10464LM500 | 110-80-5 | Cellosolve extrapure | 2909.43.00 | 500 ml. | 608 | 18% |
| 10464LC250 | | | | 2.5 ltr. | 2374 | 18% |
| 10455SG500 | 9004-34-6 | Cellulose powder (For Column Chromatography) | 2909.43.00 | 500 gm. | 1185 | 18% |
| 10472SG100 | 534-17-8 | Cesium carbonate, min 99% for synthesis | 2836.92.00 | 100 gm. | 7536 | 18% |
| 20466SG025 | 7647-17-8 | Cesium chloride AR | 2827.39.00 | 25 gm. | 2075 | 18% |
| 20466SG100 | | | | 100 gm. | 5869 | 18% |
| 20474SG010 | 10294-54-9 | Cesium sulphate, min 99.5% AR | 2833.29.90 | 10 gm. | 2464 | 18% |
| 10476SG500 | 36653-82-4 | Cetyl alcohol extrapure | 2905.17.00 | 500 gm. | 1082 | 18% |
| 10462LM500 | 79-04-9 | Chloroacetyl chloride for synthesis | 2915.90.00 | 500 ml. | 2380 | 18% |
| 70465LM500 | 108-90-7 | Chlorobenzene AR/ACS | 2903.61.10 | 500 ml. | 703 | 18% |
| 70465LC250 | | | | 2.5 ltr. | 2591 | 18% |
| 10465LM500 | 108-90-7 | Chlorobenzene extrapure | 2903.61.10 | 500 ml. | 664 | 18% |
| 10465LC250 | | | | 2.5 ltr. | 2536 | 18% |
| 20471LM500 | 67-66-3 | Chloroform AR DRY SOLVENT | 2903.13.00 | 500 ml. | 717 | 18% |
| 70470LM500 | 67-66-3 | Chloroform AR/ACS | 2903.13.00 | 500 ml. | 640 | 18% |
| 70470LC250 | | | | 2.5 ltr. | 2504 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 70470PL025 | 67-66-3 | Chloroform AR/ACS-Plastic drum | 2903.13.00 | 25 ltr. | 19623 | 18% |
| 10470LM500 | 67-66-3 | Chloroform extrapure | 2903.13.00 | 500 ml. | 481 | 18% |
| 10470LC250 | | | | 2.5 ltr. | 1949 | 18% |
| 10470PL025 | 67-66-3 | Chloroform extrapure-plastic drum | 2903.13.00 | 25 ltr. | 19623 | 18% |
| 66470LM250 | 67-66-3 | Chloroform GC Grade | 2903.13.00 | 250 ml. | 1681 | 18% |
| 30470LC250 | 67-66-3 | Chloroform HPLC & Spectroscopy | 2903.13.00 | 2.5 ltr. | 3083 | 18% |
| 30470LC100 | 67-66-3 | Chloroform HPLC & Spectroscopy(stabilized) | 2903.13.00 | 1 ltr. | 1628 | 18% |
| 10473SG001 | 16941-12-1 | Chloroplatinic acid | 3815.90.00 | 1 gm. | 14327 | 18% |
| 10480SG500 | 1308-38-9 | Chromium oxide green extrapure | 2819.90.00 | 500 gm. | 1693 | 18% |
| 10485SG500 | 1333-82-0 | Chromium trioxide extrapure | 2819.10.00 | 500 gm. | 1415 | 18% |
| 20486SG010 | 5808-22-0 | Chromotropic acid disodium salt dihydrate AR | 3204.90.00 | 10 gm. | 1529 | 18% |
| 20486SG025 | | | | 25 gm. | 2992 | 18% |
| 70489SG500 | 77-92-9 | Citric acid anhydrous AR/ACS | 2918.14.00 | 500 gm. | 1033 | 18% |
| 70489S5K25 | | | | 25 kg. | 85902 | 18% |
| 10489SG500 | 77-92-9 | Citric acid anhydrous extrapure | 2918.14.00 | 500 gm. | 687 | 18% |
| 10489SR500 | | | | 5 kg. | 3920 | 18% |
| 10489S5K25 | | | | 25 kg. | 16940 | 18% |
| 70490SG500 | 5949-29-1 | Citric acid monohydrate AR/ACS | 2918.14.00 | 500 gm. | 812 | 18% |
| 10490SG500 | 5949-29-1 | Citric acid monohydrate extrapure | 2918.14.00 | 500 gm. | 511 | 18% |
| 10490S5K50 | | | | 50 kg. | 30391 | 18% |
| 10496SG100 | 6147-53-1 | Cobalt acetate extrapure | 2915.29.90 | 100 gm. | 2356 | 18% |
| 10496SG500 | | | | 500 gm. | 5442 | 18% |
| 70500SG100 | 7791-13-1 | Cobalt chloride hexahydrate AR/ACS | 2827.10.00 | 100 gm. | 1951 | 18% |
| 70500SG500 | | | | 500 gm. | 9638 | 18% |
| 10500SG100 | 7791-13-1 | Cobalt chloride hexahydrate extrapure | 2827.10.00 | 100 gm. | 1420 | 18% |
| 10500SG500 | | | | 500 gm. | 6243 | 18% |
| 70501SG250 | 10026-22-9 | Cobalt nitrate hexahydrate AR/ACS | 2834.29.90 | 250 gm. | 4690 | 18% |
| 10501SG100 | 10026-22-9 | Cobalt nitrate hexahydrate extrapure | 2834.29.90 | 100 gm. | 1346 | 18% |
| 10501SG500 | | | | 500 gm. | 5488 | 18% |
| 10502SG100 | 1308-16-1 | Cobalt oxide extrapure | 2822.00.10 | 100 gm. | 5215 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|----------|--------|-----|
| 10503SG100 | 10026-24-1 | Cobalt sulphate heptahydrate extrapure | 2833.29.90 | 100 gm. | 1291 | 18% |
| 10503SG500 | | | | 500 gm. | 6155 | 18% |
| 10503S5K25 | | | | 25 kg. | 252035 | 18% |
| 10528BK010 | NA | Congo red indicator papers | 3204.19.00 | 10 bkt. | 486 | 18% |
| 10506SG005 | 1733-12-6 | Cresol red indicator powder | 2907.12.90 | 5 gm. | 346 | 18% |
| 10506SG025 | | | | 25 gm. | 1171 | 18% |
| 10507SG250 | 6046-93-1 | Cupric acetate extrapure | 2915.39.90 | 250 gm. | 1090 | 18% |
| 10509SG500 | 12069-69-1 | Cupric carbonate extrapure | 2836.99.90 | 500 gm. | 2245 | 18% |
| 10511SG500 | 10125-13-0 | Cupric chloride extrapure | 2827.39.40 | 500 gm. | 1692 | 18% |
| 10508SG500 | 10031-43-3 | Cupric nitrate trihydrate extrapure | 2834.29.90 | 500 gm. | 1388 | 18% |
| 10512SG500 | 1317-38-0 | Cupric oxide extrapure | 2931.00.90 | 500 gm. | 2924 | 18% |
| 12510SG500 | 7758-98-7 | Cupric sulphate anhydrous extrapure | 2833.25.00 | 500 gm. | 1771 | 18% |
| 70510SG500 | 7758-99-8 | Cupric sulphate pentahydrate AR/ACS | 2833.25.00 | 500 gm. | 1429 | 18% |
| 10510SG500 | 7758-99-8 | Cupric sulphate pentahydrate extrapure | 2833.25.00 | 500 gm. | 1078 | 18% |
| 10510SR500 | | | | 5 kg. | 9580 | 18% |
| 10510S5K50 | | | | 50 kg. | 95326 | 18% |
| 10514SG500 | 7758-89-6 | Cuprous Chloride extrapure | 2827.39.40 | 500 gm. | 1752 | 18% |
| 20516LM500 | 110-82-7 | Cyclohexane AR DRY SOLVENT | 2902.11.00 | 500 ml. | 1006 | 18% |
| 70515LM500 | 110-82-7 | Cyclohexane AR/ACS | 2902.11.00 | 500 ml. | 865 | 18% |
| 70515LC250 | | | | 2.5 ltr. | 3196 | 18% |
| 10515LM500 | 110-82-7 | Cyclohexane extrapure | 2902.11.00 | 500 ml. | 667 | 18% |
| 10515LC250 | | | | 2.5 ltr. | 2542 | 18% |
| 10515LL025 | | | | 25 ltr. | 21583 | 18% |
| 66515LM250 | 110-82-7 | Cyclohexane GC Grade | 2902.11.00 | 250 ml. | 1436 | 18% |
| 30515LC100 | 110-82-7 | Cyclohexane HPLC & Spectroscopy | 2902.11.00 | 1 ltr. | 1838 | 18% |
| 30515LC250 | | | | 2.5 ltr. | 2599 | 18% |
| 70517LM500 | 108-94-1 | Cyclohexanone AR/ACS | 2914.39.10 | 500 ml. | 948 | 18% |
| 70517LC250 | | | | 2.5 ltr. | 3984 | 18% |
| 10517LM500 | 108-94-1 | Cyclohexanone extrapure | 2914.39.10 | 500 ml. | 892 | 18% |
| 10517LC250 | | | | 2.5 ltr. | 3902 | 18% |
| 10517LL025 | | | | 25 ltr. | 29495 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 32520SG025 | 13419-61-9 | 1-Decane sulphonic acid sodium salt anhydrous HPLC | 2904.10.90 | 25 gm. | 1519 | 18% |
| 30520SG025 | 13419-61-9 | 1-Decane sulphonic acid sodium salt monohydrate HPLC | 2904.10.90 | 25 gm. | 1519 | 18% |
| 10714SG100 | 57-48-7 | D(-) Fructose extrapure | 1702.50.00 | 100 gm. | 433 | 18% |
| 10714SG500 | | | | 500 gm. | 1452 | 18% |
| 20880SG500 | 69-65-8 | D(-)-Mannitol AR | 2905.43.00 | 500 gm. | 2925 | 18% |
| 10880SG500 | 69-65-8 | D(-)-Mannitol extrapure | 2905.43.00 | 500 gm. | 2530 | 18% |
| 10880SR500 | | | | 5 kg. | 10625 | 18% |
| 70530SG500 | 50-99-7 | Dextrose anhydrous AR/ACS | 2940.00.00 | 500 gm. | 588 | 18% |
| 70530SR500 | | | | 5 kg. | 3812 | 18% |
| 10530S5K50 | 50-99-7 | Dextrose anhydrous extrapure (5X10) | 2940.00.00 | 50 kg. | 31679 | 18% |
| 10530SG500 | 50-99-7 | Dextrose anhydrous extrapure | 2940.00.00 | 500 gm. | 546 | 18% |
| 10530SR500 | | | | 5 kg. | 3321 | 18% |
| 10532LM500 | 123-42-2 | Diacetone alcohol extrapure | 2914.40.00 | 500 ml. | 669 | 18% |
| 10532LC250 | | | | 2.5 ltr. | 3186 | 18% |
| 20526SG025 | 57-71-6 | Diacetyl monoxime AR | 2928.00.10 | 25 gm. | 951 | 18% |
| 20526SG100 | | | | 100 gm. | 3120 | 18% |
| 70140SG500 | 7783-28-0 | di-Ammonium hydrogen ortho-phosphate AR/ACS | 3105.40.00 | 500 gm. | 1155 | 18% |
| 10140S5K50 | 7783-28-0 | di-Ammonium hydrogen orthophosphate extrapure (5X10) | 3105.40.00 | 50 kg. | 77000 | 18% |
| 10140SG500 | 7783-28-0 | di-Ammonium hydrogen orthophosphate extrapure | 3105.40.00 | 500 gm. | 1001 | 18% |
| 10535LM500 | 84-74-2 | Dibutyl phthalate extrapure | 2917.39.10 | 500 ml. | 729 | 18% |
| 10535LC250 | | | | 2.5 ltr. | 3030 | 18% |
| 10445SG500 | 7789-77-7 | di-Calcium phosphate dihydrate extrapure | 2835.25.00 | 500 gm. | 452 | 18% |
| 20538LM500 | 95-50-1 | 1,2-Dichlorobenzene AR | 2903.61.20 | 500 ml. | 922 | 18% |
| 20538LC250 | | | | 2.5 ltr. | 3498 | 18% |
| 10538LM500 | 95-50-1 | 1,2-Dichlorobenzene extrapure | 2903.61.20 | 500 ml. | 741 | 18% |
| 10538LC100 | | | | 1 ltr. | 1404 | 18% |
| 10539SG500 | 106-46-7 | 1,4-Dichlorobenzene for synthesis | 2903.61.30 | 500 gm. | 799 | 18% |
| 70540LM500 | 107-06-2 | 1,2-Dichloroethane AR/ACS | 2903.12.00 | 500 ml. | 667 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|---|------------|----------|--------|-----|
| 10540LM500 | 107-06-2 | 1,2-Dichloroethane extrapure | 2903.12.00 | 500 ml. | 558 | 18% |
| 10540LC250 | | | | 2.5 ltr. | 1787 | 18% |
| 10540LL025 | | | | 25 ltr. | 17248 | 18% |
| 30540LC100 | 107-06-2 | 1,2-Dichloroethane HPLC & Spectroscopy | 2903.12.00 | 1 ltr. | 2932 | 18% |
| 30540LC250 | | | | 2.5 ltr. | 4398 | 18% |
| 20551LM500 | 75-09-2 | Dichloromethane AR DRY SOLVENT | 2903.12.00 | 500 ml. | 805 | 18% |
| 20551LC250 | | | | 2.5 ltr. | 3156 | 18% |
| 70550LL025 | 75-09-2 | Dichloromethane AR/ACS | 2903.12.00 | 25 ltr. | 14380 | 18% |
| 70550LM500 | 75-09-2 | Dichloromethane AR/ACS (Methylene dichloride) | 2903.12.00 | 500 ml. | 719 | 18% |
| 70550LC250 | | | | 2.5 ltr. | 2630 | 18% |
| 70550PL025 | 75-09-2 | Dichloromethane AR/ACS- Plastic Drum | 2903.12.00 | 25 ltr. | 14525 | 18% |
| 40550LL025 | 75-09-2 | Dichloromethane CR Grade | 2903.12.00 | 25 ltr. | 14432 | 18% |
| 40550PL025 | 75-09-2 | Dichloromethane CR Grade- Plastic Drum | 2903.12.00 | 25 ltr. | 23557 | 18% |
| 10550LM500 | 75-09-2 | Dichloromethane extrapure | 2903.12.00 | 500 ml. | 649 | 18% |
| 10550LC250 | | | | 2.5 ltr. | 2273 | 18% |
| 10550LL025 | | | | 25 ltr. | 14380 | 18% |
| 10550LL200 | | | | 200 ltr. | 114800 | 18% |
| 10550LK250 | | | | 250 kg. | 148750 | 18% |
| 10550PL025 | 75-09-2 | Dichloromethane extrapure- Plastic Drum | 2903.12.00 | 25 ltr. | 14700 | 18% |
| 66550LM250 | 75-09-2 | Dichloromethane GC Grade | 2903.12.00 | 250 ml. | 1260 | 18% |
| 30550LC100 | 75-09-2 | Dichloromethane HPLC & Spectroscopy | 2903.12.00 | 1 ltr. | 1436 | 18% |
| 30550LC250 | | | | 2.5 ltr. | 3721 | 18% |
| 20598SGI00 | 100-10-7 | 4-Dimethylamino benzaldehyde AR | 2913.00.90 | 100 gm. | 1534 | 18% |
| 10553LM500 | 111-42-2 | Diethanolamine extrapure | 2921.11.10 | 500 ml. | 749 | 18% |
| 20581LM500 | 60-29-7 | Diethyl ether AR DRY SOLVENT | 2909.11.00 | 500 ml. | 784 | 18% |
| 20581LC250 | | | | 2.5 ltr. | 3119 | 18% |
| 70580LC250 | 60-29-7 | Diethyl ether AR/ACS | 2909.11.00 | 2.5 ltr. | 3023 | 18% |
| 70580LM500 | 60-29-7 | Diethyl ether AR/ACS | 2909.11.00 | 500 ml. | 777 | 18% |
| 10580LM500 | 60-29-7 | Diethyl ether extrapure | 2909.11.00 | 500 ml. | 640 | 18% |
| 10580LC250 | | | | 2.5 ltr. | 2519 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|--|------------|----------|--------|-----|
| 66580LM250 | 60-29-7 | Diethyl ether GC Grade | 2909.11.00 | 250 ml. | 2267 | 18% |
| 30580LC100 | 60-29-7 | Diethyl ether HPLC & Spectroscopy | 2909.11.00 | 1 ltr. | 5038 | 18% |
| 10584LM500 | 84-66-2 | Diethyl phthalate for synthesis | 2917.39.90 | 500 ml. | 783 | 18% |
| 10584LC250 | | | | 2.5 ltr. | 2350 | 18% |
| 70575LM500 | 109-89-7 | Diethylamine AR/ACS | 2921.11.10 | 500 ml. | 871 | 18% |
| 10575LM500 | 109-89-7 | Diethylamine extrapure | 2921.11.10 | 500 ml. | 805 | 18% |
| 10575LC250 | | | | 2.5 ltr. | 2641 | 18% |
| 10575LL025 | | | | 25 ltr. | 27064 | 18% |
| 30575LC100 | 109-89-7 | Diethylamine HPLC & Spectroscopy | 2921.11.10 | 1 ltr. | 3928 | 18% |
| 10576LM500 | 111-46-6 | Diethylene glycol for synthesis | 2909.41.00 | 500 ml. | 624 | 18% |
| 10576LC250 | | | | 2.5 ltr. | 2566 | 18% |
| 10577SG100 | 67-43-6 | Diethylene triamine pentaacetic acid extrapure | 2934.90.99 | 100 gm. | 871 | 18% |
| 10577SG500 | | | | 500 gm. | 3705 | 18% |
| 70586LM500 | 108-20-3 | Di-isopropyl ether AR/ACS | 2909.19.00 | 500 ml. | 1158 | 18% |
| 10586LM500 | 108-20-3 | Di-isopropyl ether extrapure | 2909.19.00 | 500 ml. | 955 | 18% |
| 10586LC250 | | | | 2.5 ltr. | 3257 | 18% |
| 10586LL025 | | | | 25 ltr. | 40276 | 18% |
| 10586LL200 | | | | 200 ltr. | 287112 | 18% |
| 10587LM500 | 110-71-4 | 1,2-Dimethoxyethane | 2909.19.00 | 500 ml. | 5513 | 18% |
| 10587LC250 | | | | 2.5 ltr. | 20056 | 18% |
| 20587LM500 | 110-71-4 | 1,2-Dimethoxyethane AR | 2909.19.00 | 500 ml. | 6899 | 18% |
| 20588LM500 | 110-71-4 | 1,2-Dimethoxyethane AR DRY SOLVENT | 2909.19.00 | 500 ml. | 8484 | 18% |
| 10602LM500 | 616-38-6 | Dimethyl carbonate, min 98.0% for synthesis | 2920.90.99 | 500 ml. | 802 | 18% |
| 10602LC250 | | | | 2.5 ltr. | 2995 | 18% |
| 10601SG100 | 95-45-4 | Dimethyl glyoxime extrapure | 2928.00.90 | 100 gm. | 855 | 18% |
| 20621LM500 | 67-68-5 | Dimethyl sulphoxide AR DRY SOLVENT | 2930.90.60 | 500 ml. | 2835 | 18% |
| 70620LM500 | 67-68-5 | Dimethyl sulphoxide AR/ACS | 2930.90.60 | 500 ml. | 2495 | 18% |
| 70620LC250 | | | | 2.5 ltr. | 9185 | 18% |
| 10620LM500 | 67-68-5 | Dimethyl sulphoxide extrapure | 2930.90.60 | 500 ml. | 2326 | 18% |
| 10620LC250 | | | | 2.5 ltr. | 8968 | 18% |
| 10620LL025 | | | | 25 ltr. | 98791 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|--|------------|----------|--------|-----|
| 90620LM500 | 67-68-5 | Dimethyl Sulphoxide GC-HS | 2930.90.60 | 500 ml. | 4511 | 18% |
| 90620LC100 | | | | 1 ltr. | 7165 | 18% |
| 30620LC100 | 67-68-5 | Dimethyl sulphoxide HPLC & Spectroscopy | 2930.90.60 | 1 ltr. | 4432 | 18% |
| 30620LC250 | | | | 2.5 ltr. | 9777 | 18% |
| 10596SG500 | 506-59-2 | Dimethylammonium chloride | 2825.10.90 | 500 gm. | 1153 | 18% |
| 10624SG025 | 119-26-6 | 2,4-Dinitrophenyl hydrazine for synthesis | 2905.41.00 | 25 gm. | 423 | 18% |
| 10624SG100 | | | | 100 gm. | 1254 | 18% |
| 10625SG025 | 609-99-4 | 3,5-Dinitrosalicylic acid for synthesis | 2918.29.90 | 25 gm. | 2682 | 18% |
| 10625SG100 | | | | 100 gm. | 8510 | 18% |
| 10628SG500 | 577-11-7 | Dioctyl sodium sulphosuccinate 98% (DOSS)(Aerosol OT) | 2833.11.00 | 500 gm. | 3196 | 18% |
| 70630LM500 | 123-91-1 | 1,4-Dioxan (stabilized) AR/ACS | 2932.99.00 | 500 ml. | 984 | 18% |
| 70630LC250 | | | | 2.5 ltr. | 4231 | 18% |
| 20631LM500 | 123-91-1 | 1,4-Dioxan AR DRY SOLVENT | | 500 ml. | 1043 | 18% |
| 10630LM500 | 123-91-1 | 1,4-Dioxan extrapure | 2932.99.00 | 500 ml. | 961 | 18% |
| 10630LC250 | | | | 2.5 ltr. | 3156 | 18% |
| 10630LL025 | | | | 25 ltr. | 49840 | 18% |
| 66630LM250 | 123-91-1 | 1,4-Dioxan GC Grade | 2932.99.00 | 250 ml. | 1126 | 18% |
| 30630LC100 | 123-91-1 | 1,4-Dioxan HPLC & Spectroscopy | 2932.99.00 | 1 ltr. | 2692 | 18% |
| 20633SG100 | 122-39-4 | Diphenylamine AR | 2921.11.10 | 100 gm. | 772 | 18% |
| 20637SG010 | 140-22-7 | 1,5-Diphenyl carbazide AR | 2928.00.10 | 10 gm. | 1567 | 18% |
| 20637SG025 | | | | 25 gm. | 2976 | 18% |
| 20638SG005 | 538-62-5 | 1,5-Diphenyl carbazone AR (reagent for mercury) | 2928.00.10 | 5 gm. | 1414 | 18% |
| 20638SG025 | | | | 25 gm. | 5039 | 18% |
| 11280SG500 | 7758-11-4 | di-Potassium hydrogen o-phosphate extrapure | 2835.24.00 | 500 gm. | 1096 | 18% |
| 11280SR500 | | | | 5 kg. | 8846 | 18% |
| 11280S5K50 | | | | 50 kg. | 118076 | 18% |
| 71280SG500 | 7758-11-4 | di-Potassium hydrogen ortho-phosphate AR/ACS | 2835.24.00 | 500 gm. | 1242 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|---------|--------|-----|
| 11610SG500 | 7558-79-4 | di-Sodium hydrogen o-phosphate anhydrous extrapure | 2835.22.00 | 500 gm. | 1114 | 18% |
| 11610SR500 | | | | 5 kg. | 6903 | 18% |
| 11610S5K50 | | | | 50 kg. | 88263 | 18% |
| 31610SG500 | 7558-79-4 | di-Sodium hydrogen o-phosphate anhydrous HPLC | 2835.22.00 | 500 gm. | 2039 | 18% |
| 21611SG500 | 10028-24-7 | di-Sodium hydrogen o-phosphate dihydrate AR | 2835.22.00 | 500 gm. | 1082 | 18% |
| 11611SG500 | 10028-24-7 | di-Sodium hydrogen o-phosphate dihydrate extrapure | 2835.22.00 | 500 gm. | 765 | 18% |
| 11611SR500 | | | | 5 kg. | 6212 | 18% |
| 11611S5K50 | | | | 50 kg. | 75460 | 18% |
| 71610SG500 | 7558-79-4 | di-Sodium hydrogen ortho-phosphate anhydrous AR/ACS | 2835.22.00 | 500 gm. | 1251 | 18% |
| 71610S5K50 | | | | 50 kg. | 150276 | 18% |
| 11697SG500 | 6106-24-7 | di-Sodium tartrate dihydrate EP | 2918.13.90 | 500 gm. | 1492 | 18% |
| 71700SG500 | 1303-96-4 | di-Sodium tetraborate decahydrate AR/ACS (Borax) | 2840.19.00 | 500 gm. | 812 | 18% |
| 11700SG500 | 1303-96-4 | di-Sodium tetraborate decahydrate extrapure | 2840.19.00 | 500 gm. | 494 | 18% |
| 11700SR500 | | | | 5 kg. | 3813 | 18% |
| 20650SG005 | 60-10-6 | Dithizone AR | 2930.90.99 | 5 gm. | 1771 | 18% |
| 20650SG025 | | | | 25 gm. | 7431 | 18% |
| 10864SG500 | 617-48-1 | DL-Malic acid extrapure | 2918.19.00 | 500 gm. | 1038 | 18% |

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|----------|--------|-----|
| 10655SG100 | 60-00-4 | EDTA acid extrapure | 2921.21.00 | 100 gm. | 433 | 18% |
| 10655SG500 | | | | 500 gm. | 1442 | 18% |
| 10657SG100 | 25102-12-9 | EDTA di-potassium salt extrapure | 2917.39.90 | 100 gm. | 492 | 18% |
| 10657SG500 | | | | 500 gm. | 1888 | 18% |
| 10660SG100 | 6381-92-6 | EDTA di-sodium salt dihydrate extrapure | 2921.21.00 | 100 gm. | 419 | 18% |
| 10660SG500 | | | | 500 gm. | 1465 | 18% |
| 10660SR500 | | | | 5 kg. | 11522 | 18% |
| 10660S5K50 | | | | 50 kg. | 69213 | 18% |
| 10659LM500 | NA | EDTA N/50 solution | 2921.21.00 | 500 ml. | 300 | 18% |
| 70660SG100 | 6381-92-6 | EDTA, di - sodium salt dihydrate AR/ACS | 2921.21.00 | 100 gm. | 537 | 18% |
| 70660SG500 | | | | 500 gm. | 1700 | 18% |
| 70660SR500 | | | | 5 kg. | 14045 | 18% |
| 10662SG025 | 17372-87-1 | Eosine yellow (w/s) for microscopy | 3204.14.81 | 25 gm. | 578 | 18% |
| 10662SG100 | | | | 100 gm. | 1712 | 18% |
| 20666SG250 | 8007-09-8 | Eschka's mixture AR | 2836.99.20 | 250 gm. | 26882 | 18% |
| 20665LM500 | 107-21-1 | Ethanediol AR | 2905.31.00 | 500 ml. | 845 | 18% |
| 20665LC250 | | | | 2.5 ltr. | 3002 | 18% |
| 20665LL025 | | | | 25 ltr. | 24850 | 18% |
| 10665LM500 | 107-21-1 | Ethanediol extrapure | 2905.31.00 | 500 ml. | 731 | 18% |
| 10665LC250 | | | | 2.5 ltr. | 2594 | 18% |
| 10665LL025 | | | | 25 ltr. | 24540 | 18% |
| 10670LM500 | 141-43-5 | Ethanolamine extrapure | 2921.11.10 | 500 ml. | 694 | 18% |
| 10670LC250 | | | | 2.5 ltr. | 2874 | 18% |
| 20673SG001 | 1239-45-8 | Ethidium bromide AR | 2827.59.90 | 1 gm. | 2066 | 18% |
| 20679LM500 | 141-78-6 | Ethyl Acetate AR DRY SOLVENT | 2915.31.00 | 500 ml. | 602 | 18% |
| 70680LM500 | 141-78-6 | Ethyl acetate AR/ACS | 2915.31.00 | 500 ml. | 544 | 18% |
| 70680LC250 | | | | 2.5 ltr. | 2153 | 18% |
| 70680LL025 | | | | 25 ltr. | 16891 | 18% |
| 40680LL025 | 141-78-6 | Ethyl acetate CR Grade | 2915.31.00 | 25 ltr. | 19215 | 18% |
| 40680LL200 | | | | 200 ltr. | 135226 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|---|------------|----------|--------|-----|
| 10680LM500 | 141-78-6 | Ethyl Acetate extrapure | 2915.31.00 | 500 ml. | 488 | 18% |
| 10680LC250 | | | | 2.5 ltr. | 1548 | 18% |
| 10680LL025 | | | | 25 ltr. | 17150 | 18% |
| 10680GL025 | 141-78-6 | Ethyl Acetate extrapure - GI Drum | 2915.31.00 | 25 ltr. | 21093 | 18% |
| 10680PL200 | 141-78-6 | Ethyl Acetate extrapure - Plastic drum | 2915.31.00 | 200 ltr. | 134400 | 18% |
| 66680LM250 | 141-78-6 | Ethyl acetate GC Grade | 2915.31.00 | 250 ml. | 1246 | 18% |
| 30680LC100 | 141-78-6 | Ethyl Acetate HPLC & Spectroscopy | 2915.31.00 | 1 ltr. | 1285 | 18% |
| 30680LC250 | | | | 2.5 ltr. | 2698 | 18% |
| 10677LM100 | 75-03-6 | Ethyl iodide for synthesis (Iodoethane) | 2903.22.00 | 100 ml. | 4082 | 18% |
| 10677LM250 | 75-03-6 | Ethyl iodide for synthesis | 2903.22.00 | 250 ml. | 7777 | 18% |
| 10681LM500 | 107-15-3 | Ethylenediamine extrapure | 2921.21.00 | 500 ml. | 1344 | 18% |
| 10681LC250 | | | | 2.5 ltr. | 4911 | 18% |

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 10682LM500 | NA | Fehlings solution A | 3822.00.90 | 500 ml. | 509 | 12% |
| 10683LM500 | NA | Fehlings solution B | 3822.00.90 | 500 ml. | 712 | 12% |
| 10685S5K50 | 7705-08-0 | Ferric chloride anhydrous extrapure (5X10) | 2827.10.00 | 50 kg. | 17763 | 18% |
| 10685SG500 | 7705-08-0 | Ferric chloride anhydrous extrapure | 2827.10.00 | 500 gm. | 450 | 18% |
| 10686SG500 | 10025-77-1 | Ferric chloride hexahydrate extrapure | 2827.10.00 | 500 gm. | 866 | 18% |
| 10689SG500 | 1309-37-1 | Ferric oxide red extrapure | 2821.10.10 | 500 gm. | 667 | 18% |
| 20690LM100 | 14634-91-4 | Ferroun solution AR | 3822.00.90 | 100 ml. | 1293 | 12% |
| 70700SG500 | 7782-63-0 | Ferrous sulphate heptahydrate AR/ACS | 2833.29.10 | 500 gm. | 664 | 18% |
| 10700S5K50 | 7782-63-0 | Ferrous sulphate heptahydrate extrapure (5X10) | 2833.29.10 | 50 kg. | 17255 | 18% |
| 10700SG500 | 7782-63-0 | Ferrous sulphate heptahydrate extrapure | 2833.29.10 | 500 gm. | 416 | 18% |
| 10700SR500 | | | | 5 kg. | 2591 | 18% |
| 10702SR100 | 1317-37-9 | Ferrous sulphide sticks | 2830.90.10 | 1 kg. | 1057 | 18% |
| 10705LC100 | 11043-90-6 | FINKLIN | 3402.19.00 | 1 ltr. | 848 | 18% |
| 10705LC500 | | | | 5 ltr. | 2080 | 18% |
| 10706LM125 | NA | Folin & Ciocalteu's phenol reagent | 3822.00.90 | 125 ml. | 767 | 12% |
| 10706LM250 | | | | 250 ml. | 1470 | 12% |
| 70710LC250 | 50-00-0 | Formaldehyde 37-41% soln AR/ACS | 2912.11.00 | 2.5 ltr. | 2070 | 18% |
| 70710LL030 | 50-00-0 | Formaldehyde 37-41% soln AR/ACS 30LT | 2912.11.00 | 30 ltr. | 8925 | 18% |
| 10710LL030 | 50-00-0 | Formaldehyde 37-41% soln extrapure 30LT | 2912.11.00 | 30 ltr. | 8820 | 18% |
| 10710LL100 | | | | 100 ltr. | 27370 | 18% |
| 10710LM500 | 50-00-0 | Formaldehyde 37-41% soln extrapure | 2912.11.00 | 500 ml. | 386 | 18% |
| 10710LC500 | | | | 5 ltr. | 1684 | 18% |
| 70710LM500 | 50-00-0 | Formaldehyde solution 37-41% solution AR/ACS | 2912.11.00 | 500 ml. | 488 | 18% |
| 10711LM500 | 75-12-7 | Formamide for synthesis | 2921.11.10 | 500 ml. | 976 | 18% |
| 10711LC250 | | | | 2.5 ltr. | 3514 | 18% |
| 10712LM500 | 64-18-6 | Formic acid 85% extrapure | 2915.11.00 | 500 ml. | 504 | 18% |
| 10712LC250 | | | | 2.5 ltr. | 2066 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|--|------------|---------|-------|-----|
| 10716SG025 | 632-99-5 | Fuchsin basic (CI 42510) for microscopy | 3204.13.43 | 25 gm. | 482 | 18% |
| 10716SG100 | | | | 100 gm. | 1261 | 18% |
| 10708SG500 | 110-17-8 | Fumaric acid for synthesis | 2917.19.90 | 500 gm. | 1074 | 18% |
| 10708S5K50 | | | | 50 kg. | 44562 | 18% |
| 10717SG500 | NA | Fusion mixture | 2836.40.00 | 500 gm. | 647 | 18% |



| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|--------|-----|
| 10748LM125 | NA | Gention violet aqueous stain solution for microscopy | 3204.13.00 | 125 ml. | 370 | 18% |
| 10747SG025 | 548-62-9 | Gention violet powder (M.S) (C.I.42535) | 3204.13.00 | 25 gm. | 546 | 18% |
| 10747SG100 | | | | 100 gm. | 1468 | 18% |
| 10718SG025 | 51811-82-6 | Giemsa stain for microscopy | 3204.19.90 | 25 gm. | 968 | 18% |
| 10718SG100 | | | | 100 gm. | 3355 | 18% |
| 10719LM125 | NA | Giemsa stain solution for microscopy | 3204.19.90 | 125 ml. | 438 | 18% |
| 10723LC100 | 111-30-8 | Glutaraldehyde soln. 25% for synthesis | 2912.19.10 | 1 ltr. | 1877 | 18% |
| 10723LM500 | 111-30-8 | Glutaraldehyde solution 25% for synthesis | 2912.19.10 | 500 ml. | 1074 | 18% |
| 70720LM500 | 56-81-5 | Glycerine AR/ACS | 2905.45.00 | 500 ml. | 1157 | 18% |
| 70720LC250 | | | | 2.5 ltr. | 4084 | 18% |
| 70720LL025 | | | | 25 ltr. | 19936 | 18% |
| 10720LM500 | 56-81-5 | Glycerine extrapure | 2905.45.00 | 500 ml. | 848 | 18% |
| 10720LC250 | | | | 2.5 ltr. | 3848 | 18% |
| 10720LC500 | | | | 500 ml. | 5771 | 18% |
| 10720LL025 | | | | 25 ltr. | 19740 | 18% |
| 70721SG250 | 56-40-6 | Glycine AR/ACS | 2922.49.10 | 250 gm. | 897 | 18% |
| 10721SG500 | 56-40-6 | Glycine extrapure | 2922.49.10 | 500 gm. | 1363 | 18% |
| 10721S5K50 | | | | 50 kg. | 139650 | 18% |
| 10724SG001 | 16903-35-8 | Gold chloride | 2843.30.00 | 1 gm. | 13511 | 18% |
| 10725LM125 | NA | Gram's iodine stain solution for microscopy | 2801.20.00 | 125 ml. | 340 | 12% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|----------|--------|-----|
| 10731SG025 | 517-28-2 | Haematoxylin stain certified | 3203.00.10 | 25 gm. | 21064 | 18% |
| 10731SG005 | 517-28-2 | Haematoxylin stain certified | 3203.00.10 | 5 gm. | 4639 | 18% |
| 10727LM500 | 142-82-5 | Heptane for synthesis | 2901.29.20 | 500 ml. | 735 | 18% |
| 10727LC250 | | | | 2.5 ltr. | 2869 | 18% |
| 10727LL025 | | | | 25 ltr. | 26308 | 18% |
| 37730SG025 | 22767-50-6 | 1-Heptane sulphonic acid sodium salt anhydrous HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 30730SG025 | 22767-50-6 | 1-Heptane sulphonic acid sodium salt monohydrate HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 10737SG500 | 100-97-0 | Hexamine extrapure | 2921.29.10 | 500 gm. | 814 | 18% |
| 70740LM500 | 110-54-3 | Hexane (fraction from petroleum) AR/ACS | 2710.19.90 | 500 ml. | 776 | 18% |
| 70740LC250 | | | | 2.5 ltr. | 2747 | 18% |
| 40740LL025 | 110-54-3 | Hexane CR Grade | 2710.19.90 | 25 ltr. | 14490 | 18% |
| 10740LM500 | 110-54-3 | Hexane extrapure | 2710.19.90 | 500 ml. | 650 | 18% |
| 10740LC250 | | | | 2.5 ltr. | 2082 | 18% |
| 30740LC250 | 110-54-3 | Hexane for HPLC & Spectroscopy | 2710.19.90 | 2.5 ltr. | 3856 | 18% |
| 30740LC100 | 110-54-3 | Hexane HPLC & Spectroscopy | 2710.19.90 | 1 ltr. | 1894 | 18% |
| 37760SG025 | 2832-45-3 | 1-Hexane sulphonic acid sodium salt anhydrous HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 30760SG025 | 2832-45-3 | 1-Hexane sulphonic acid sodium salt monohydrate HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 10762LM500 | 111-27-3 | 1-Hexanol for synthesis | 2906.13.90 | 500 ml.* | 1404 | 18% |
| 10763SG010 | 12055-62-8 | Holmium(III) oxide AR | 2822.00.10 | 10 gm. | 8260 | 18% |
| 10765LM500 | 7803-57-8 | Hydrazine hydrate 80% extrapure | 2825.10.20 | 500 ml. | 1619 | 18% |
| 10765LC250 | 7803-57-8 | | | 2.5 ltr. | 7037 | 18% |
| 10765LL025 | 7803-57-8 | | | 25 ltr. | 58065 | 18% |
| 10765LL200 | 7803-57-8 | | | 200 ltr. | 469000 | 18% |
| 10766LM500 | 7803-57-8 | Hydrazine hydrate 99 - 100% extrapure | 2825.10.20 | 500 ml. | 2146 | 18% |
| 10766LL025 | | | | 25 ltr. | 125965 | 18% |
| 20767SG100 | 10034-93-2 | Hydrazine Sulphate AR | 2825.10.00 | 100 gm. | 745 | 18% |
| 10767SG100 | 10034-93-2 | Hydrazine Sulphate extrapure | 2825.10.00 | 100 gm. | 691 | 18% |
| 10771LM500 | 7647-01-0 | Hydrochloric acid 0.1 N solution | 2806.10.00 | 500 ml. | 337 | 18% |
| 70770L2C50 | 7647-01-0 | Hydrochloric acid AR/ACS 2X5LT | 2806.10.00 | 2x5 ltr. | 3245 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|------------|--------|-----|
| 70770L8M50 | 7647-01-0 | Hydrochloric acid AR/ACS-Glass bottle | 2806.10.00 | 8x500 ml. | 3562 | 18% |
| 70770L4C25 | | | | 4x2.5 ltr. | 4968 | 18% |
| 10770L2C50 | 7647-01-0 | Hydrochloric acid extrapure | 2806.10.00 | 2x5 ltr. | 2361 | 18% |
| 10770L4C25 | | | | 4x2.5 ltr. | 3804 | 18% |
| 10770LL025 | | | | 25 ltr. | 6179 | 18% |
| 10770L8M50 | 7647-01-0 | Hydrochloric acid extrapure-Glass bottle | 2806.10.00 | 8x500 ml. | 3260 | 18% |
| 20773LM500 | 7664-39-3 | Hydrofluoric acid 40% AR | 2811.11.00 | 500 ml. | 1626 | 18% |
| 10773LM500 | 7664-39-3 | Hydrofluoric acid 40% extrapure | 2811.11.00 | 500 ml. | 947 | 18% |
| 10773LC500 | | | | 5 ltr. | 3749 | 18% |
| 70774LM500 | 7664-39-3 | Hydrofluoric acid 48% AR/ACS | 2811.11.00 | 500 ml. | 1646 | 18% |
| 10777L2M50 | 7722-84-1 | Hydrogen peroxide 50% w/v extrapure | 2847.00.00 | 20x500 ml. | 13350 | 12% |
| 10777LK050 | | | | 50 ltr | 84000 | 18% |
| 10776L2M50 | 7722-84-1 | Hydrogen peroxide solution 30% w/v extrapure | 2847.00.00 | 20x500 ml. | 11411 | 12% |
| 70776L2M50 | 7722-84-1 | Hydrogen peroxide solution 30% w/v AR/ACS | 2847.00.00 | 20x500 ml. | 12062 | 12% |
| 20778S5K50 | 123-31-9 | Hydroquinone AR (5X10) | 2907.22.00 | 50 kg. | 277550 | 18% |
| 20778SG500 | 123-31-9 | Hydroquinone AR | 2907.22.00 | 500 gm. | 3883 | 18% |
| 10778SG500 | 123-31-9 | Hydroquinone extrapure | 2907.22.00 | 500 gm. | 3232 | 18% |
| 20780SG100 | 5470-11-1 | Hydroxylamine hydrochloride AR | 2825.10.90 | 100 gm. | 784 | 18% |
| 20780SG500 | | | | 500 gm. | 3094 | 18% |
| 10780SG100 | 5470-11-1 | Hydroxylamine hydrochloride extrapure | 2825.10.90 | 100 gm. | 723 | 18% |
| 10780SG500 | | | | 500 gm. | 2907 | 18% |
| 10785SG100 | 10039-54-0 | Hydroxylamine sulphate extrapure | 2825.10.40 | 100 gm. | 424 | 18% |
| 10785SG500 | | | | 500 gm. | 1427 | 18% |
| 10787SG005 | 63451-35-4 | Hydroxynaphthol blue indicator | 2907.29.90 | 5 gm. | 624 | 18% |
| 10787SG025 | | | | 25 gm. | 2758 | 18% |
| 10759SR100 | 9004-63-5 | Hydroxypropyl Methyl Cellulose, LV Extrapure | 3912.39.19 | 1 kg. | 5716 | 18% |
| 20786SG100 | 148-24-3 | 8-Hydroxyquinoline AR | 2923.40.90 | 100 gm. | 2817 | 18% |
| 10786SG100 | 148-24-3 | 8-Hydroxyquinoline extrapure | 2923.40.90 | 100 gm. | 1374 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|--|------------|----------|-------|-----|
| 70790SG100 | 288-32-4 | Imidazole AR/ACS | 2933.29.90 | 100 gm. | 1173 | 18% |
| 10790SG100 | 288-32-4 | Imidazole extrapure | 2933.29.90 | 100 gm. | 956 | 18% |
| 10790SG500 | | | | 500 gm. | 2622 | 18% |
| 10795LM100 | NA | Immersion oil for microscopy | 3822.00.90 | 100 ml. | 2336 | 12% |
| 10795LM030 | NA | Immersion oil for microscopy 30ML | 3822.00.90 | 30 ml. | 780 | 12% |
| 20792SG025 | 860-22-0 | Indigo carmine AR | 3204.19.00 | 25 gm. | 552 | 18% |
| 20792SG100 | | | | 100 gm. | 1546 | 18% |
| 10796SG005 | 133-32-4 | Indole-3-butyric acid 99% for biochemistry | 2933.90.90 | 5 gm. | 1464 | 18% |
| 10796SG025 | | | | 25 gm. | 5568 | 18% |
| 70800SG100 | 7553-56-2 | Iodine (resublimed) AR/ACS | 2801.20.00 | 100 gm. | 4149 | 12% |
| 70800SG500 | | | | 500 gm. | 18158 | 12% |
| 10801LM500 | 7553-56-2 | Iodine 0.1N solution | 2801.20.00 | 500 ml. | 1525 | 12% |
| 10800SG100 | 7553-56-2 | Iodine extrapure | 2801.20.00 | 100 gm. | 3878 | 12% |
| 10800SG500 | | | | 500 gm. | 17600 | 12% |
| 10798SG050 | 7790-99-0 | Iodine monochloride extrapure | 2801.20.00 | 5x50 gm. | 8990 | 12% |
| 10799SG025 | 865-44-1 | Iodine trichloride extrapure | 2801.20.00 | 25 gm. | 4018 | 12% |
| 10803LC500 | NA | Iodophor | 3808.94.00 | 5 ltr. | 4706 | 18% |
| 10802SR100 | 7439-89-6 | Iron powder electrolytic | 7205.29.10 | 1 kg. | 1827 | 18% |
| 10802S5K50 | | | | 50 kg. | 97815 | 18% |
| 10202LM500 | 123-92-9 | iso-Amyl acetate purified | 2915.31.00 | 500 ml. | 856 | 18% |
| 10202LC250 | | | | 2.5 ltr. | 3257 | 18% |
| 10202LL025 | | | | 25 ltr. | 43430 | 18% |
| 10203LM500 | 123-51-3 | iso-Amyl alcohol (for synthesis and milk testing) purified | 2915.39.90 | 500 ml. | 689 | 18% |
| 10203LC250 | | | | 2.5 ltr. | 2757 | 18% |
| 10365LM500 | 78-83-1 | iso-Butyl alcohol extrapure | 2912.19.10 | 500 ml. | 579 | 18% |
| 10365LC250 | | | | 2.5 ltr. | 2198 | 18% |
| 21041LM500 | 540-84-1 | iso-Octane AR DRY SOLVENT | 2901.10.00 | 500 ml. | 1715 | 18% |
| 71040LM500 | 540-84-1 | iso-Octane AR/ACS | 2901.10.00 | 500 ml. | 1546 | 18% |
| 71040LC250 | | | | 2.5 ltr. | 6490 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|--------------------------------|------------|----------|------|-----|
| 11040LM500 | 540-84-1 | iso-Octane extrapure | 2901.10.00 | 500 ml. | 1301 | 18% |
| 11040LC250 | | | | 2.5 ltr. | 5351 | 18% |
| 66040LM250 | 540-84-1 | Iso-Octane GC Grade | 2901.10.00 | 250 ml. | 1368 | 18% |
| 31040LC100 | 540-84-1 | iso-Octane HPLC & Spectroscopy | 2901.10.00 | 1 ltr. | 3457 | 18% |

K

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|---------|----------------------|------------|-----------|------|-----|
| 20805LM500 | NA | Karl Fischer Reagent | 3822.00.90 | 500 ml. | 2838 | 12% |
| 20805LM250 | | | | 2x250 ml. | 2915 | 12% |

L

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|---------|------------------------------|------------|---------|-------|-----|
| 10220SG100 | 50-81-7 | L(+) Ascorbic acid extrapure | 2936.27.00 | 100 gm. | 1931 | 18% |
| 10220SG500 | | | | 500 gm. | 7021 | 18% |
| 10220SR500 | | | | 5 kg. | 20895 | 18% |
| 70220SG100 | 50-81-7 | L(+) Ascorbic acid AR/ACS | 2936.27.00 | 100 gm. | 2135 | 18% |
| 70220SG500 | | | | 500 gm. | 8034 | 18% |
| 11770SG500 | 87-69-4 | L(+) Tartaric acid extrapure | 2918.12.00 | 500 gm. | 2003 | 18% |
| 11770S5K50 | | | | 50 kg. | 75933 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|---------|-------|-----|
| 71770SG500 | 87-69-4 | L(+) Tartaric acid AR/ACS | 2918.12.00 | 500 gm. | 2277 | 18% |
| 20808LM500 | 50-21-5 | Lactic acid AR | 2918.11.10 | 500 ml. | 1405 | 18% |
| 10808LM500 | 50-21-5 | Lactic acid extrapure | 2918.11.10 | 500 ml. | 1143 | 18% |
| 10808LL025 | | | | 25 ltr. | 43141 | 18% |
| 10808LC500 | | | | 5 ltr. | 11320 | 18% |
| 10811SG500 | 10039-26-6 | Lactose monohydrate | 1702.11.10 | 500 gm. | 867 | 18% |
| 20812SG100 | 10025-84-0 | Lanthanum Chloride AR | 2827.39.00 | 100 gm. | 5558 | 18% |
| 20813SG100 | 10277-43-7 | Lanthanum nitrate AR | 2846.90.00 | 100 gm. | 6138 | 18% |
| 20816SG100 | 1312-81-8 | Lanthanum oxide AR | 2846.90.00 | 100 gm. | 5409 | 18% |
| 10846SR250 | 1335-32-6 | Lead acetate basic anhydrous (Lead sub acetate for Sugar Analysis) | 2915.29.90 | 2.5 kg. | 4522 | 18% |
| 10810SG500 | 6080-56-4 | Lead acetate trihydrate extrapure | 2915.29.90 | 500 gm. | 764 | 18% |
| 10815SG500 | 1344-36-1 | Lead carbonate Basic extrapure | 3109.09.01 | 500 gm. | 2144 | 18% |
| 10817SG500 | 1317-36-8 | Lead monoxide (Litharge) extrapure | 2824.10.10 | 500 gm. | 1058 | 18% |
| 10820SG500 | 10099-74-8 | Lead nitrate extrapure | 2834.29.90 | 500 gm. | 834 | 18% |
| 10824SG025 | 12627-53-1 | Leishman's stain for microscopy | 3204.19.00 | 25 gm. | 854 | 18% |
| 10824SG100 | | | | 100 gm. | 2584 | 18% |
| 10825LM250 | 12627-53-1 | Leishman's stain solution for microscopy | 3204.19.00 | 250 ml. | 401 | 18% |
| 10823SG100 | 16853-85-3 | Lithium aluminium hydride for synthesis | 2850.00.10 | 100 gm. | 34100 | 18% |
| 10826SG500 | 7550-35-8 | Lithium bromide anhydrous extrapure | 2827.59.90 | 500 gm. | 9426 | 18% |
| 10827SG250 | 554-13-2 | Lithium carbonate extrapure | 2836.91.00 | 250 gm. | 3389 | 18% |
| 10827SG500 | | | | 500 gm. | 5992 | 18% |
| 10828SG250 | 7447-41-8 | Lithium chloride anhydrous extrapure | 2827.39.90 | 250 gm. | 7040 | 18% |
| 10828SG500 | | | | 500 gm. | 11220 | 18% |
| 10829SG500 | 1310-66-3 | Lithium hydroxide monohydrate extrapure | 2825.20.00 | 500 gm. | 11440 | 18% |
| 10831SG500 | 7790-69-4 | Lithium nitrate anhydrous extrapure | 2834.29.90 | 500 gm. | 11330 | 18% |
| 10841BK010 | | Litmus paper blue | 3822.00.00 | 10 bkt. | 427 | 12% |
| 10842BK010 | | Litmus paper red | 3822.00.00 | 10 bkt. | 427 | 12% |
| 10839SG025 | 1393-92-6 | Litmus pH indicator | 2932.29.10 | 25 gm. | 642 | 18% |
| 10830LM125 | | Lugol's solution | 2827.60.10 | 125 ml. | 523 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|---------|-------|-----|
| 10834SG500 | 7439-95-4 | Magnesium (metal turning) | 2827.31.00 | 500 gm. | 1771 | 18% |
| 20836SG250 | 16674-78-5 | Magnesium acetate tetrahydrate AR | 2915.29.20 | 250 gm. | 4147 | 18% |
| 10836SG500 | 16674-78-5 | Magnesium acetate tetrahydrate extrapure | 2915.29.20 | 500 gm. | 1214 | 18% |
| 10838SG500 | 39409-82-0 | Magnesium carbonate light | 2836.99.20 | 500 gm. | 1777 | 18% |
| 70840SG500 | 7791-18-6 | Magnesium chloride hexahydrate AR/ACS | 2827.31.00 | 500 gm. | 533 | 18% |
| 70840S5K50 | | | | 50 kg. | 13482 | 18% |
| 10840SG500 | 7791-18-6 | Magnesium chloride hexahydrate extrapure | 2827.31.00 | 500 gm. | 400 | 18% |
| 10840SR500 | | | | 5 kg. | 2589 | 18% |
| 10840S5K50 | | | | 50 kg. | 13622 | 18% |
| 10847SG500 | 13446-18-9 | Magnesium nitrate extrapure | 2834.29.20 | 500 gm. | 451 | 18% |
| 10844SG500 | 1309-48-4 | Magnesium oxide light extrapure | 2836.99.20 | 500 gm. | 1600 | 18% |
| 10844SK020 | | | | 20 kg. | 63301 | 18% |
| 70850SG500 | 10034-99-8 | Magnesium sulphate heptahydrate AR/ACS | 2833.21.00 | 500 gm. | 503 | 18% |
| 10850S5K50 | 10034-99-8 | Magnesium sulphate heptahydrate extrapure (5X10) | 2833.21.00 | 50 kg. | 16226 | 18% |
| 10850SG500 | 10034-99-8 | Magnesium sulphate heptahydrate extrapure | 2833.21.00 | 500 gm. | 338 | 18% |
| 10850SR500 | | | | 5 kg. | 2177 | 18% |
| 10860SG025 | 2437-29-8 | Malachite green (M.S.) extrapure | 3204.13.61 | 25 gm. | 370 | 18% |
| 10860SG100 | | | | 100 gm. | 928 | 18% |
| 10860SG500 | | | | 500 gm. | 3495 | 18% |
| 10863SG500 | 110-16-7 | Maleic acid extrapure | 2917.19.10 | 500 gm. | 975 | 18% |
| 10866SG250 | 6363-53-7 | Maltose monohydrate for biochemistry | 1702.50.00 | 250 gm. | 979 | 18% |
| 70870SG500 | 10034-96-5 | Manganese sulfate monohydrate AR/ACS | 2833.29.40 | 500 gm. | 1197 | 18% |
| 10867SG500 | 6156-78-1 | Manganous acetate tetrahydrate extrapure | 2915.29.20 | 500 gm. | 623 | 18% |
| 70868SG500 | 13446-34-9 | Manganous chloride tetrahydrate AR/ACS | | 500 gm. | 1037 | 18% |
| 10868SG500 | 13446-34-9 | Manganous chloride tetrahydrate extrapure | 2827.39.90 | 500 gm. | 697 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 10870SG500 | 10034-96-5 | Manganous sulphate monohydrate extrapure | 2833.29.40 | 500 gm. | 797 | 18% |
| 10870SR500 | | | | 5 kg. | 6765 | 18% |
| 10870S5K50 | | | | 50 kg. | 85708 | 18% |
| 20505LM500 | 108-39-4 | m-Cresol AR | 2907.12.90 | 500 ml. | 2343 | 18% |
| 10505LM500 | 108-39-4 | m-Cresol extrapure | 2907.12.90 | 500 ml. | 1920 | 18% |
| 10497SG001 | 2303-01-7 | m-Cresol purple indicator | 2932.29.90 | 1 gm. | 480 | 18% |
| 10497SG005 | | | | 5 gm. | 1846 | 18% |
| 20890SG100 | 1600-27-7 | Mercuric acetate AR | 2915.29.90 | 100 gm. | 8435 | 18% |
| 10890SG100 | 1600-27-7 | Mercuric acetate extrapure | 2915.29.90 | 100 gm. | 7698 | 18% |
| 10890SG500 | | | | 500 gm. | 29107 | 18% |
| 20895SG100 | 7789-47-1 | Mercuric bromide AR | 2827.59.90 | 100 gm. | 23522 | 18% |
| 10901LMI25 | NA | Mercuric chloride 5% w/v solution | 2812.10.47 | 125 ml. | 505 | 18% |
| 20900SG100 | 7487-94-7 | Mercuric chloride AR | 2812.10.47 | 100 gm. | 6071 | 18% |
| 20900SG250 | | | | 250 gm. | 14460 | 18% |
| 10900SG100 | 7487-94-7 | Mercuric chloride extrapure | 2812.10.47 | 100 gm. | 5928 | 18% |
| 10900SG250 | | | | 250 gm. | 13073 | 18% |
| 10904SG100 | 7774-29-0 | Mercuric iodide red extrapure | 2812.10.47 | 100 gm. | 5934 | 18% |
| 10910SG100 | 7783-34-8 | Mercuric nitrate monohydrate extrapure | 2834.29.90 | 100 gm. | 6266 | 18% |
| 10910SG500 | | | | 500 gm. | 25849 | 18% |
| 10911SG100 | 21908-53-2 | Mercuric oxide red extrapure | 2817.00.10 | 100 gm. | 7001 | 18% |
| 10912SG100 | 21908-53-2 | Mercuric oxide yellow extrapure | 2817.00.10 | 100 gm. | 6835 | 18% |
| 20920SG250 | 7783-35-9 | Mercuric sulphate AR | 2833.21.00 | 250 gm. | 15355 | 18% |
| 10920SG250 | 7783-35-9 | Mercuric sulphate extrapure | 2833.21.00 | 250 gm. | 14805 | 18% |
| 10925SG100 | 592-85-8 | Mercuric Thiocyanate extrapure | 2833.21.00 | 100 gm. | 7190 | 18% |
| 20921SG100 | 10112-91-1 | Mercurous chloride AR | 2827.39.20 | 100 gm. | 6675 | 18% |
| 10922SG100 | 10415-75-5 | Mercurous nitrate extrapure | 2834.29.30 | 100 gm. | 7063 | 18% |
| 20888SG250 | 7439-97-6 | Mercury (metal) for polarography, 99.99% AR | 2805.40.00 | 250 gm. | 15644 | 18% |
| 10888SG250 | 7439-97-6 | Mercury (metal) pure | 2805.40.00 | 250 gm. | 12349 | 18% |
| 10926LM500 | 75-75-2 | Methane sulphonic acid for synthesis | 2904.90.20 | 500 ml. | 2042 | 18% |
| 10926LC250 | | | | 2.5 ltr. | 7040 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|--|------------|----------|-------|-----|
| 10928LM500 | 124-63-0 | Methanesulphonyl chloride for synthesis | 2904.90.20 | 500 ml. | 3283 | 18% |
| 10928LC250 | | | | 2.5 ltr. | 8280 | 18% |
| 10931LM500 | 67-56-1 | Methanol anhydrous (suitable for Karl Fischer) | 2905.11.00 | 500 ml. | 535 | 18% |
| 20940LM500 | 67-56-1 | Methanol anhydrous AR DRY SOLVENT | 2905.11.00 | 500 ml. | 551 | 18% |
| 20940LC250 | | | | 2.5 ltr. | 1626 | 18% |
| 70930LM500 | 67-56-1 | Methanol AR/ACS | 2905.11.00 | 500 ml. | 405 | 18% |
| 70930LC250 | | | | 2.5 ltr. | 925 | 18% |
| 70930LL025 | | | | 25 ltr. | 6811 | 18% |
| 30932LC250 | 67-56-1 | Methanol EL grade | 2905.11.00 | 2.5 ltr. | 1361 | 18% |
| 30932LL025 | | | | 25 ltr. | 8656 | 18% |
| 10930LM500 | 67-56-1 | Methanol extrapure | 2905.11.00 | 500 ml. | 331 | 18% |
| 10930LC250 | | | | 2.5 ltr. | 709 | 18% |
| 10930LL025 | | | | 25 ltr. | 6650 | 18% |
| 10930GL025 | 67-56-1 | Methanol extrapure - GI drum | 2905.11.00 | 25 ltr. | 9333 | 18% |
| 10930PL200 | 67-56-1 | Methanol extrapure - Plastic drum | 2905.11.00 | 200 ltr. | 49000 | 18% |
| 10929LM500 | 67-56-1 | Methanol for synthesis | 2905.11.00 | 500 ml. | 266 | 18% |
| 10929LC250 | | | | 2.5 ltr. | 697 | 18% |
| 66930LM250 | 67-56-1 | Methanol GC grade | 2905.11.00 | 250 ml. | 884 | 18% |
| 40930LC250 | 67-56-1 | Methanol Gradient Grade | 2905.11.00 | 2.5 ltr. | 963 | 18% |
| 30930LC250 | 67-56-1 | Methanol HPLC & Spectroscopy | 2905.11.00 | 2.5 ltr. | 940 | 18% |
| 10935LC250 | 107-98-2 | 1-Methoxy-2-Propanol extrapure | 2909.49.00 | 2.5 ltr. | 4111 | 18% |
| 10933LM500 | 79-20-9 | Methyl acetate for synthesis | 2905.11.00 | 500 ml. | 1408 | 18% |
| 10945LM500 | 109-86-4 | Methyl cellosolve extrapure | 2909.43.00 | 500 ml. | 818 | 18% |
| 10945LC250 | | | | 2.5 ltr. | 3544 | 18% |
| 10944SG500 | 9004-67-5 | Methyl cellulose 300-560 | 3912.39.12 | 500 gm. | 5962 | 18% |
| 10946SG500 | 9004-67-5 | Methyl cellulose 4000cps | 3912.39.12 | 500 gm. | 6212 | 18% |
| 10957LM250 | | | | 250 ml. | 17358 | 18% |
| 10960SG025 | 547-58-0 | Methyl orange (pH indicator) | 2927.00.90 | 25 gm. | 362 | 18% |
| 10960SG100 | | | | 100 gm. | 1139 | 18% |
| 10960SR100 | | | | 1 kg. | 9117 | 18% |
| 10961LM125 | 547-58-0 | Methyl orange solution | 2927.00.90 | 125 ml. | 231 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-------------|---|------------|----------|--------|-----|
| 70965LM500 | 108-10-1 | 4-Methylpentan-2-one AR/ACS (Methyl isobutyl ketone) | 2914.13.00 | 500 ml. | 966 | 18% |
| 70965LC250 | | | | 2.5 ltr. | 3355 | 18% |
| 10965LM500 | 108-10-1 | 4-Methylpentan-2-one extrapure | | 500 ml. | 687 | 18% |
| 10965LC250 | | | | 2.5 ltr. | 2972 | 18% |
| 10965LL025 | | | | 25 ltr. | 29715 | 18% |
| 10980SG025 | 493-52-7 | Methyl red (pH indicator) | 2927.00.90 | 25 gm. | 427 | 18% |
| 10980SG100 | | | | 100 gm. | 1542 | 18% |
| 10980SR100 | | | | 1 kg. | 13782 | 18% |
| 10981LM125 | 493-52-7 | Methyl red solution | 2927.00.90 | 125 ml. | 255 | 18% |
| 10962LM500 | 119-36-8 | Methyl salicylate extrapure | 2918.23.10 | 500 ml. | 1313 | 18% |
| 10982LC250 | 96-47-9 | 2-Methyl tetrahydrofuran 98%, for synthesis | 2932.11.00 | 2.5 ltr. | 11675 | 18% |
| 10982LM500 | 96-47-9 | 2-Methyl tetrahydrofuran, 98% for synthesis | | 500 ml. | 3198 | 18% |
| 10982LL025 | | | | 25 ltr. | 102277 | 18% |
| 10950SG025 | 61-73-4 | Methylene blue for microscopy | 3822.00.90 | 25 gm. | 787 | 12% |
| 10950SG100 | | | | 100 gm. | 2532 | 12% |
| 10951LM125 | 61-73-4 | Methylene blue solution | 3822.00.90 | 125 ml. | 278 | 12% |
| 10972SG250 | 308080-99-1 | Molecular Sieves 3A 1/16" pellets | 3823.90.98 | 250 gm. | 996 | 18% |
| 10971SG250 | 308080-99-1 | Molecular Sieves 3A 1/8" pellets | 3823.90.98 | 250 gm. | 996 | 18% |
| 10974SG250 | 70995-01-0 | Molecular Sieves 4A 1/16" pellets | 3823.90.98 | 250 gm. | 1105 | 18% |
| 10973SG250 | 70995-01-0 | Molecular sieves 4A 1/8" pellets | 3823.90.98 | 250 gm. | 1020 | 18% |
| 10976SG250 | 1344-01-0 | Molecular Sieves 5A 1/16" pellets | 3823.90.98 | 250 gm. | 1020 | 18% |
| 10975SG250 | 1344-01-0 | Molecular Sieves 5A 1/8" pellets | 3823.90.98 | 250 gm. | 1020 | 18% |
| 10985SG100 | 1313-27-5 | Molybdenum trioxide extrapure | 2825.70.10 | 100 gm. | 4259 | 18% |
| 70990SG100 | 7782-91-4 | Molybdic acid hydrate AR/ACS | 2825.70.20 | 100 gm. | 3262 | 18% |
| 10990SG100 | 7782-91-4 | Molybdic acid hydrate extrapure | 2825.70.20 | 100 gm. | 2522 | 18% |
| 10990SG500 | | | | 500 gm. | 10179 | 18% |
| 10995SG100 | 96-27-5 | Monothioglycerol | 2905.59.00 | 100 gm. | 6007 | 18% |
| 10995SR100 | | | | 1 kg. | 33789 | 18% |
| 10996LM500 | 110-91-8 | Morpholine for synthesis | 2933.39.17 | 500 ml. | 1138 | 18% |
| 10996LC250 | | | | 2.5 ltr. | 4262 | 18% |
| 10996LL025 | | | | 25 ltr. | 38290 | 18% |

(*till stock lasts)

M

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---|------------|---------|------|-----|
| 70999SG005 | 3051-09-0 | Murexide AR/ACS (Ammonium Purpurate) | 2923.90.00 | 5 gm. | 699 | 18% |
| 70999SG025 | | | | 25 gm. | 2015 | 18% |
| 10999SG005 | 3051-09-0 | Murexide extrapure | 2923.90.00 | 5 gm. | 605 | 18% |
| 10999SG025 | | | | 25 gm. | 1720 | 18% |

N

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---|------------|----------|-------|-----|
| 70970LM500 | 872-50-4 | N - Methyl - 2 - pyrrolidone AR/ACS | 2933.79.00 | 500 ml. | 3127 | 18% |
| 71014SG005 | 1465-25-4 | N-(1-Naphthyl)ethylenediamine dihydrochloride 99% AR/ACS | 2921.21.00 | 5 gm. | 1744 | 18% |
| 71014SG010 | 1465-25-4 | N-(1-Naphthyl)ethylenediamine dihydrochloride 99% AR/ACS | 2921.21.00 | 10 gm. | 3295 | 18% |
| 71014SG025 | | | | 25 gm. | 7273 | 18% |
| 20590LM500 | 127-19-5 | N,N-Dimethyl acetamide AR | 2924.29.90 | 500 ml. | 933 | 18% |
| 20590LC250 | | | | 2.5 ltr. | 3127 | 18% |
| 20591LM500 | 127-19-5 | N,N-Dimethyl acetamide AR DRY SOLVENT | 2924.29.90 | 500 ml. | 1005 | 18% |
| 10590LM500 | 127-19-5 | N,N-Dimethyl acetamide extrapure | 2924.29.90 | 500 ml. | 839 | 18% |
| 10590LC250 | | | | 2.5 ltr. | 2192 | 18% |
| 10590LL025 | | | | 25 ltr. | 37646 | 18% |
| 90590LM500 | 127-19-5 | N,N-Dimethyl acetamide GC-HS | 2924.29.90 | 500 ml. | 1359 | 18% |
| 90590LC100 | | | | 1 ltr. | 2317 | 18% |
| 30590LC100 | 127-19-5 | N,N-Dimethyl acetamide HPLC & Spectroscopy | 2924.29.90 | 1 ltr. | 1920 | 18% |
| 20599LM500 | 68-12-2 | N,N-Dimethyl formamide AR DRY SOLVENT | 2924.29.90 | 500 ml. | 699 | 18% |
| 20599LC250 | | | | 2.5 ltr. | 3202 | 18% |
| 70600LM500 | 68-12-2 | N,N-Dimethyl formamide AR/ACS | 2924.29.90 | 500 ml. | 683 | 18% |
| 70600LC250 | | | | 2.5 ltr. | 2785 | 18% |
| 70600LL025 | | | | 25 ltr. | 21478 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|--|------------|----------|--------|-----|
| 10600LM500 | 68-12-2 | N,N-Dimethyl formamide extrapure | 2924.29.90 | 500 ml. | 639 | 18% |
| 10600LC250 | | | | 2.5 ltr. | 2415 | 18% |
| 10600LL025 | | | | 25 ltr. | 18608 | 18% |
| 10600LL200 | | | | 200 ltr. | 133126 | 18% |
| 90600LM500 | 68-12-2 | N,N-Dimethyl formamide GC-HS | 2924.29.90 | 500 ml. | 1504 | 18% |
| 90600LC100 | | | | 1 ltr. | 2510 | 18% |
| 30600LC100 | 68-12-2 | N,N-Dimethyl formamide HPLC & Spectroscopy | 2924.29.90 | 1 ltr. | 1473 | 18% |
| 10597LM500 | 121-69-7 | N,N-Dimethylaniline extrapure | 2921.49.90 | 500 ml. | 846 | 18% |
| 11000SG500 | 91-20-3 | Naphthalene extrapure | 2902.90.40 | 500 gm. | 850 | 18% |
| 11012SG005 | 145-50-6 | 1-Naphtholbenzein | 2914.50.00 | 5 gm. | 1501 | 18% |
| 11012SG025 | | | | 25 gm. | 5711 | 18% |
| 11003SG025 | 86-87-3 | 1-Naphthylacetic acid | 2916.39.00 | 25 gm. | 467 | 18% |
| 11003SG100 | | | | 100 gm. | 1209 | 18% |
| 10360LM500 | 123-86-4 | n-Butyl acetate extrapure | 2915.33.00 | 500 ml. | 595 | 18% |
| 10360LC250 | | | | 2.5 ltr. | 2333 | 18% |
| 11007LM125 | 7774-29-0 | Nessler's reagent (for ammonia & ammonium) | 3822.00.90 | 125 ml. | 525 | 12% |
| 20728LM500 | 142-82-5 | n-Heptane 95% AR | 2901.29.90 | 500 ml. | 1238 | 18% |
| 20728LC250 | | | | 2.5 ltr. | 4554 | 18% |
| 20729LM500 | 142-82-5 | n-Heptane 99% AR | 2901.29.90 | 500 ml. | 1828 | 18% |
| 20729LC250 | | | | 2.5 ltr. | 7991 | 18% |
| 20729LL025 | | | | 25 ltr. | 29843 | 18% |
| 10729LM500 | 142-82-5 | n-Heptane 99% for synthesis | 2901.29.90 | 500 ml. | 1506 | 18% |
| 10729LC250 | | | | 2.5 ltr. | 6949 | 18% |
| 10729LL025 | | | | 25 ltr. | 28583 | 18% |
| 30729LC100 | 142-82-5 | n-Heptane 99% HPLC & Spectroscopy | 2901.29.90 | 1 ltr. | 3019 | 18% |
| 30729LC250 | | | | 2.5 ltr. | 5907 | 18% |
| 66729LM250 | 142-82-5 | n-Heptane GC Grade | 2901.29.90 | 250 ml. | 1677 | 18% |
| 20749LM500 | 110-54-3 | n-Hexane 95% AR | 2901.10.00 | 500 ml. | 1322 | 18% |
| 20749LC250 | | | | 2.5 ltr. | 4638 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|------------|-------|-----|
| 20750LM500 | 110-54-3 | n-Hexane 99% AR | 2901.10.00 | 500 ml. | 1527 | 18% |
| 20750LC250 | | | | 2.5 ltr. | 4916 | 18% |
| 20750LL025 | | | | 25 ltr. | 42840 | 18% |
| 20751LM500 | 110-54-3 | n-Hexane 99% AR DRY SOLVENT | 2901.10.00 | 500 ml. | 1944 | 18% |
| 10750LM500 | 110-54-3 | n-Hexane 99% for synthesis | 2901.10.00 | 500 ml. | 1482 | 18% |
| 10750LC250 | | | | 2.5 ltr. | 4796 | 18% |
| 10750LL025 | | | | 25 ltr. | 42840 | 18% |
| 66750LM250 | 110-54-3 | n-Hexane 99% GC grade | 2901.10.00 | 250 ml. | 1502 | 18% |
| 30750LC100 | 110-54-3 | n-Hexane 99% HPLC & Spectroscopy | 2901.10.00 | 1 ltr. | 2246 | 18% |
| 30750LC250 | | | | 2.5 ltr. | 5328 | 18% |
| 11013SG500 | 373-02-4 | Nickel acetate tetrahydrate extrapure | 2915.29.90 | 500 gm. | 6027 | 18% |
| 21017SG500 | 12635-29-9 | Nickel aluminium alloy powder 50% AR | 3815.11.00 | 500 gm. | 6084 | 18% |
| 11008SG500 | 39430-27-8 | Nickel carbonate extrapure | 2836.99.90 | 500 gm. | 5822 | 18% |
| 11009SG500 | 7791-20-0 | Nickel chloride hexahydrate extrapure | 2827.35.00 | 500 gm. | 3719 | 18% |
| 11010SG500 | 10101-97-0 | Nickel sulphate hexahydrate extrapure | 2833.24.00 | 500 gm. | 3191 | 18% |
| 71020SG010 | 485-47-2 | Ninhydrin AR/ACS | 2914.50.00 | 10 gm. | 1254 | 18% |
| 71020SG025 | | | | 25 gm. | 3071 | 18% |
| 71030L8M50 | 7697-37-2 | Nitric acid AR/ACS | 2808.00.10 | 8x500 ml. | 4301 | 18% |
| 71030L4C25 | | | | 4x2.5 ltr. | 6118 | 18% |
| 11030L8M50 | 7697-37-2 | Nitric acid extrapure | 2808.00.10 | 8x500 ml. | 3542 | 18% |
| 11030L4C25 | | | | 4x2.5 ltr. | 5376 | 18% |
| 71037LM500 | 75-52-5 | Nitromethane AR/ACS | 2929.90.00 | 500 ml. | 2565 | 18% |
| 21039SG005 | 4264-83-9 | 4-Nitrophenyl phosphate disodium salt AR | 2918.21.90 | 5 gm. | 1854 | 18% |
| 10970LM500 | 872-50-4 | N-Methyl-2-Pyrrolidone extrapure | 2933.79.00 | 500 ml. | 2279 | 18% |
| 10970LC250 | | | | 2.5 ltr. | 8802 | 18% |
| 10970LL025 | | | | 25 ltr. | 75318 | 18% |
| 90970LM500 | 872-50-4 | N-Methyl-2-pyrrolidone GC-HS | 2933.79.00 | 500 ml. | 5268 | 18% |
| 90970LC100 | | | | 1 ltr. | 8776 | 18% |
| 30970LC100 | 872-50-4 | N-Methyl-2-pyrrolidone HPLC & Spectroscopy | 2933.79.00 | 1 ltr. | 7286 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|-------------------------------|------------|----------|-------|-----|
| 21088LM500 | 109-66-0 | n-Pentane AR | 2902.90.90 | 500 ml. | 1135 | 18% |
| 21088LC250 | | | | 2.5 ltr. | 4252 | 18% |
| 21088LL025 | | | | 25 ltr. | 49726 | 18% |
| 21089LM500 | 109-66-0 | n-Pentane AR DRY SOLVENT | 2902.90.90 | 500 ml. | 1512 | 18% |
| 11088LM500 | 109-66-0 | n-Pentane for synthesis | 2902.90.90 | 500 ml. | 969 | 18% |
| 11088LC250 | | | | 2.5 ltr. | 3872 | 18% |
| 66088LM250 | 109-66-0 | n-Pentane GC Grade | 2902.90.90 | 250 ml. | 1709 | 18% |
| 31088LC100 | 109-66-0 | n-Pentane HPLC & Spectroscopy | 2902.90.90 | 1 ltr. | 4118 | 18% |



| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 37050SG025 | 5324-84-5 | 1-Octane sulphononic acid sodium salt anhydrous HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 31050SG025 | 5324-84-5 | 1-Octane sulphononic acid sodium salt monohydrate HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 11052LM500 | 111-87-5 | 1-Octanol for synthesis | 2905.16.90 | 500 ml. | 1627 | 18% |
| 11052LC250 | | | | 2.5 ltr. | 7961 | 18% |
| 11055LC250 | 112-80-1 | Oleic acid extrapure | 3823.12.00 | 2.5 ltr. | 3126 | 18% |
| 11057SG025 | 1936-15-8 | Orange G for microscopy | 2927.00.00 | 25 gm. | 684 | 18% |
| 11059SG010 | 6153-39-5 | Orcinol Monohydrate | 2907.29.00 | 10 gm. | 2862 | 18% |
| 11190LM500 | 7664-38-2 | ortho-Phosphoric acid | 2809.20.10 | 500 ml. | 991 | 18% |
| 11190LC250 | | | | 2.5 ltr. | 4243 | 18% |
| 11190LL025 | | | | 25 ltr. | 36108 | 18% |
| 71190LM500 | 7664-38-2 | ortho-Phosphoric acid 85% AR/ACS | 2809.20.10 | 500 ml. | 1200 | 18% |
| 71190LC250 | | | | 2.5 ltr. | 5092 | 18% |
| 71190LL025 | | | | 25 ltr. | 36108 | 18% |
| 21191LM500 | 7664-38-2 | ortho-Phosphoric acid 88% AR | 2809.20.10 | 500 ml. | 1245 | 18% |
| 21191LC250 | | | | 2.5 ltr. | 5407 | 18% |
| 21191LL025 | | | | 25 ltr. | 70000 | 18% |
| 11191LM500 | 7664-38-2 | ortho-Phosphoric acid 88% extrapure | 2809.20.10 | 500 ml. | 1129 | 18% |
| 11191LC250 | | | | 2.5 ltr. | 4918 | 18% |
| 11191LL025 | | | | 25 ltr. | 68810 | 18% |
| 31190LM500 | 7664-38-2 | ortho-Phosphoric acid for HPLC | 2809.20.10 | 500 ml. | 3483 | 18% |
| 31190LC100 | 7664-38-2 | ortho-Phosphoric acid HPLC | 2809.20.10 | 1 ltr. | 5687 | 18% |
| 71060SG001 | 20816-12-0 | Osmic acid 99.95% AR/ACS | 2843.90.19 | 1 gm. | 14513 | 18% |
| 11867LM500 | 119-93-7 | o-Tolidine Reagent | 2921.43.30 | 500 ml. | 381 | 18% |
| 11867LL025 | | | | 25 ltr. | 39375 | 18% |
| 11071LM500 | 6153-56-6 | Oxalic acid 0.1N solution | 2917.11.10 | 500 ml. | 366 | 18% |
| 71070SG500 | 6153-56-6 | Oxalic acid dihydrate AR/ACS | 2917.11.10 | 500 gm. | 609 | 18% |
| 11070SG500 | 6153-56-6 | Oxalic acid dihydrate extrapure | 2917.11.10 | 500 gm. | 538 | 18% |
| 11070SR500 | | | | 5 kg. | 4001 | 18% |
| 11070S5K50 | | | | 50 kg. | 46456 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|------------|-------|-----|
| 11074SG001 | 7647-10-1 | Palladium [II] chloride extrapure | 2843.90.12 | 1 gm. | 14951 | 18% |
| 11085SG010 | NA | Palladium-charcoal activated (5% pd) Hydrogenation catalyst | 3815.12.90 | 10 gm. | 9517 | 18% |
| 11078LM500 | 8012-95-1 | Paraffin liquid heavy | 2710.19.90 | 500 ml. | 736 | 18% |
| 11078LC250 | | | | 2.5 ltr. | 2502 | 18% |
| 11079LM500 | 8012-95-1 | Paraffin liquid light | 2710.19.90 | 500 ml. | 500 | 18% |
| 11079LC250 | | | | 2.5 ltr. | 1721 | 18% |
| 11079LL025 | | | | 25 ltr. | 17064 | 18% |
| 11080SG500 | 8002-74-2 | Paraffin wax 58-60°C | 2712.20.90 | 500 gm. | 961 | 18% |
| 11080SR200 | | | | 2 kg. | 3077 | 18% |
| 11082SG500 | 8002-74-2 | Paraffin wax 60-62°C | 2712.20.90 | 500 gm. | 947 | 18% |
| 11082SR200 | | | | 2 kg. | 3432 | 18% |
| 11087SG500 | 8002-74-2 | Paraffin wax with Cerasin congealing point 60°C | 2712.20.90 | 500 gm. | 1060 | 18% |
| 11083SG500 | 30525-89-4 | Paraformaldehyde 96% for synthesis | 2912.60.00 | 500 gm. | 565 | 18% |
| 21086SG005 | 3737-95-9 | Patton & reeder s reagent AR | 2927.00.90 | 5 gm. | 543 | 18% |
| 21086SG025 | | | | 25 gm. | 1455 | 18% |
| 37090SG025 | 22767-49-3 | 1-Pentane sulphonic acid sodium salt anhydrous HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 31090SG025 | 22767-49-3 | 1-Pentane sulphonic acid sodium salt monohydrate HPLC | 2904.10.90 | 25 gm. | 1563 | 18% |
| 11121LM500 | 7601-90-3 | Perchloric acid 0.1 N solution | 2829.90.10 | 500 ml. | 1648 | 18% |
| 71110L8M50 | 7601-90-3 | Perchloric acid 60% AR/ACS | 2829.90.10 | 8x500 ml. | 15263 | 18% |
| 71120L8M50 | 7601-90-3 | Perchloric acid 70% AR/ACS | 2829.90.10 | 8x500 ml. | 16841 | 18% |
| 71120L4C25 | | | | 4x2.5 ltr. | 35436 | 18% |
| 21124SG025 | 10450-60-9 | Periodic acid AR | 2829.90.30 | 25 gm. | 3081 | 18% |
| 21124SG100 | | | | 100 gm. | 9504 | 18% |
| 71130LC250 | 8032-32-4 | Petroleum ether 35°-60°C AR/ACS | 2710.19.80 | 2.5 ltr. | 3604 | 18% |
| 71130LM500 | 8032-32-4 | Petroleum ether 35°-60°C AR/ACS | 2710.19.80 | 500 ml. | 948 | 18% |
| 21130LM500 | 8032-32-4 | Petroleum ether 40°-60°C AR | 2710.19.80 | 500 ml. | 1127 | 18% |
| 21130LC250 | | | | 2.5 ltr. | 5010 | 18% |
| 21130LL025 | | | | 25 ltr. | 44835 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|--|------------|----------|--------|-----|
| 11130LM500 | 8032-32-4 | Petroleum ether 40°-60°C extrapure | 2710.19.80 | 500 ml. | 849 | 18% |
| 11130LC250 | | | | 2.5 ltr. | 2969 | 18% |
| 11130LL025 | | | | 25 ltr. | 44415 | 18% |
| 11130LL200 | | | | 200 ltr. | 329000 | 18% |
| 31130LC100 | 8032-32-4 | Petroleum ether 40°-60°C HPLC & Spectroscopy | 2710.19.80 | 1 ltr. | 2045 | 18% |
| 21140LM500 | 8032-32-4 | Petroleum ether 60-80 AR | 2710.19.80 | 500 ml. | 557 | 18% |
| 21140LC250 | | | | 2.5 ltr. | 1856 | 18% |
| 11140LM500 | 8032-32-4 | Petroleum ether 60-80 extrapure | 2710.19.80 | 500 ml. | 471 | 18% |
| 11140LC250 | | | | 2.5 ltr. | 1760 | 18% |
| 11140LL025 | | | | 25 ltr. | 15400 | 18% |
| 11140LL200 | | | | 200 ltr. | 100154 | 18% |
| 11141LM500 | 8032-32-4 | Petroleum ether 80-100°C for synthesis | 2710.19.80 | 500 ml. | 593 | 18% |
| 11141LC250 | | | | 2.5 ltr. | 2143 | 18% |
| 21148BK010 | NA | pH indicator paper full range pH 1.0-14.0 | 3822.00.00 | 10 bkt. | 396 | 12% |
| 11148BK060 | NA | pH Indicator paper multi range (each of the above 6 ranges 10 bkt) | 3822.00.00 | 60 bkt. | 2119 | 12% |
| 11142BK010 | NA | pH indicator paper pH 2.0-4.5 | 3822.00.00 | 10 bkt. | 395 | 12% |
| 11143BK010 | NA | pH Indicator paper pH 3.5-6.0 | 3822.00.00 | 10 bkt. | 395 | 12% |
| 11144BK010 | NA | pH indicator paper pH 5.0-7.5 | 3822.00.00 | 10 bkt. | 395 | 12% |
| 11145BK010 | NA | pH Indicator paper pH 6.5-9.0 | 3822.00.00 | 10 bkt. | 395 | 12% |
| 11146BK010 | NA | pH Indicator paper pH 8.0-10.5 | 3822.00.00 | 10 bkt. | 395 | 12% |
| 11147BK010 | NA | pH Indicator paper wide range pH 2.0-10.5 | 3822.00.00 | 10 bkt. | 396 | 12% |
| 21149SG005 | 3829-86-5 | 1,10 Phenanthroline hydrochloride AR | 2933.99.00 | 5 gm. | 1465 | 18% |
| 21149SG025 | | | | 25 gm. | 6129 | 18% |
| 21150SG005 | 5144-89-8 | 1,10 Phenanthroline monohydrate AR | 2933.99.00 | 5 gm. | 1462 | 18% |
| 21150SG025 | | | | 25 gm. | 5765 | 18% |
| 71160SG500 | 108-95-2 | Phenol AR/ACS | 2907.11.10 | 500 gm. | 887 | 18% |
| 11160SG500 | 108-95-2 | Phenol extrapure | 2907.11.10 | 500 gm. | 841 | 18% |
| 21180SG025 | 143-74-8 | Phenol red AR | 2934.99.00 | 25 gm. | 820 | 18% |
| 21180SG100 | | | | 100 gm. | 2496 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|---------|--------|-----|
| 11181LM125 | 143-74-8 | Phenol red indicator solution | 2934.99.00 | 125 ml. | 235 | 18% |
| 71170SG100 | 77-09-8 | Phenolphthalein AR/ACS | 2932.29.10 | 100 gm. | 1175 | 18% |
| 11170SG100 | 77-09-8 | Phenolphthalein indicator powder | 2932.29.10 | 100 gm. | 854 | 18% |
| 11170SR100 | | | | 50 kg. | 511203 | 18% |
| 11170S5K25 | 77-09-8 | Phenolphthalein indicator powder (5X5) | 2932.29.10 | 25 kg. | 239608 | 18% |
| 11171LM125 | 77-09-8 | Phenolphthalein indicator solution | 2932.29.10 | 125 ml. | 217 | 18% |
| 11173LM500 | 122-99-6 | Phenoxyethanol (2-Phenoxyethanol) | 2909.49.90 | 500 ml. | 1930 | 18% |
| 11179LM250 | 100-63-0 | Phenyl hydrazine for synthesis | 2928.00.00 | 250 ml. | 2155 | 18% |
| 21185SG100 | 8003-05-2 | Phenylmercuric nitrate (basic) AR | 2931.00.90 | 100 gm. | 77653 | 18% |
| 21187SG025 | 108-73-6 | Phloroglucinol AR | 2905.19.90 | 25 gm. | 2183 | 18% |
| 21187SG100 | | | | 100 gm. | 8136 | 18% |
| 21188SG025 | 51429-74-4 | Phosphomolybdic acid AR | 2825.90.90 | 25 gm. | 2344 | 18% |
| 11193SG500 | 1314-56-3 | Phosphorous pentoxide | 2809.10.00 | 500 gm. | 1588 | 18% |
| 21193SG500 | 1314-56-3 | Phosphorous pentoxide AR | 2809.10.00 | 500 gm. | 2001 | 18% |
| 11189SG100 | 12501-23-4 | Phosphotungstic acid extrapure | 2825.90.90 | 100 gm. | 3524 | 18% |
| 11201SG500 | 88-99-3 | Phthalic acid extrapure | 2917.39.90 | 500 gm. | 1022 | 18% |
| 11203SG500 | 88-89-1 | Picric acid (moistened with water) extrapure | 2808.00.10 | 500 gm. | 1425 | 18% |
| 11204LM500 | 88-89-8 | Picric acid solution (saturated aqueous solution) | 2808.00.10 | 500 ml. | 1410 | 18% |
| 20971LM250 | | Platinum Cobalt Color Reference Solution(Hazen 500) 250ML | 3822.00.19 | 250 ml. | 32500 | 18% |
| 11216LM500 | 25322-68-3 | Polyethylene glycol 400 for synthesis | 3907.20.90 | 500 ml. | 1320 | 18% |
| 11219SG500 | 25322-68-3 | Polyethylene glycol 6000 for synthesis | 3907.20.90 | 500 gm. | 1153 | 18% |
| 71205SG500 | 127-08-2 | Potassium acetate AR/ACS | 2915.29.90 | 500 gm. | 1075 | 18% |
| 71205S5K50 | 127-08-2 | Potassium acetate AR/ACS(5X10) | | 50 kg. | 118563 | 18% |
| 11205SG500 | 127-08-2 | Potassium acetate extrapure | 2915.29.90 | 500 gm. | 933 | 18% |
| 11206SG250 | 28300-74-5 | Potassium antimony tartrate | 2918.13.00 | 250 gm. | 1714 | 18% |
| 11209S5K50 | 298-14-6 | Potassium bicarbonate extrapure (5X10) | | 50 kg. | 59976 | 18% |
| 11209SG500 | 298-14-6 | Potassium bicarbonate extrapure | 2836.99.90 | 500 gm. | 806 | 18% |
| 11212SG500 | 7646-93-7 | Potassium bisulphate purified | 2833.29.90 | 500 gm. | 804 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|---------|--------|-----|
| 11213SG500 | 868-14-4 | Potassium bitartrate extrapure | 2918.13.90 | 500 gm. | 1662 | 18% |
| 11220SG500 | 7758-01-2 | Potassium bromate extrapure | 2829.90.20 | 500 gm. | 1446 | 18% |
| 71230SG500 | 7758-02-3 | Potassium bromide AR/ACS | 2827.51.20 | 500 gm. | 1390 | 18% |
| 11230SG500 | 7758-02-3 | Potassium bromide extrapure | 2827.51.20 | 500 gm. | 1135 | 18% |
| 11230SR500 | | | | 5 kg. | 10781 | 18% |
| 31230SG100 | 7758-02-3 | Potassium bromide IR Spectroscopy | 2827.51.20 | 100 gm. | 6431 | 18% |
| 21235SG500 | 584-08-7 | Potassium carbonate granular AR | 2836.40.00 | 500 gm. | 804 | 18% |
| 11235S5K50 | 584-08-7 | Potassium carbonate granular extrapure (5X10) | 2836.40.00 | 50 kg. | 54394 | 18% |
| 11235SG500 | 584-08-7 | Potassium carbonate granular extrapure | 2836.40.00 | 500 gm. | 556 | 18% |
| 11235SR500 | | | | 5 kg. | 4223 | 18% |
| 21240SG500 | 7447-40-7 | Potassium chloride AR | 3104.20.00 | 500 gm. | 411 | 18% |
| 21240SR500 | | | | 5 kg. | 3831 | 18% |
| 11240S5K50 | 7447-40-7 | Potassium chloride extrapure (5X10) | 3104.20.00 | 50 kg. | 29400 | 18% |
| 11240SG500 | 7447-40-7 | Potassium chloride extrapure | 3104.20.00 | 500 gm. | 332 | 18% |
| 11240SR500 | | | | 5 kg. | 2362 | 18% |
| 11241SG001 | 16921-30-5 | Potassium chloroplatinate | 2843.90.19 | 1 gm. | 11348 | 18% |
| 21250SG500 | 7789-00-6 | Potassium chromate AR | 2841.50.90 | 500 gm. | 1576 | 18% |
| 11250SG500 | 7789-00-6 | Potassium chromate extrapure | 2841.50.90 | 500 gm. | 1459 | 18% |
| 11261LM500 | NA | Potassium dichromate 0.25 N solution | 2841.50.90 | 500 ml. | 822 | 18% |
| 71260SG500 | 7778-50-9 | Potassium dichromate AR/ACS | 2841.50.90 | 500 gm. | 1527 | 18% |
| 11260SG500 | 7778-50-9 | Potassium dichromate extrapure | 2841.50.90 | 500 gm. | 1394 | 18% |
| 11260SR500 | | | | 5 kg. | 12620 | 18% |
| 11270S5K50 | 7778-77-0 | Potassium di-hydrogen o-phosphate extrapure (5X10) | 2835.24.00 | 50 kg. | 85260 | 18% |
| 11270SG500 | 7778-77-0 | Potassium di-hydrogen o-phosphate extrapure | 2835.24.00 | 500 gm. | 832 | 18% |
| 11270SR500 | | | | 5 kg. | 6029 | 18% |
| 71270SG500 | 7778-77-0 | Potassium di-hydrogen ortho-phosphate AR/ACS | 2835.24.00 | 500 gm. | 1093 | 18% |
| 71270SR500 | | | | 5 kg. | 7636 | 18% |
| 31270SG500 | 7778-77-0 | Potassium di-hydrogen ortho-phosphate for HPLC | 2835.24.00 | 500 gm. | 1470 | 18% |
| 11275S5K25 | 13746-66-2 | Potassium ferricyanide EP (5X5) | 2837.20.20 | 25 kg. | 131740 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|---------|-------|-----|
| 11275SG500 | 13746-66-2 | Potassium ferricyanide extrapure | 2837.20.20 | 500 gm. | 2252 | 18% |
| 11276SG500 | 14459-95-1 | Potassium ferrocyanide extrapure | 2837.20.30 | 500 gm. | 1152 | 18% |
| 21278SG500 | 7789-23-3 | Potassium fluoride AR | 2826.19.90 | 500 gm. | 2144 | 18% |
| 11278SG500 | 7789-23-3 | Potassium fluoride extrapure | 2826.19.90 | 500 gm. | 1454 | 18% |
| 71282SG500 | 877-24-7 | Potassium hydrogen phthalate AR/ACS | 2917.39.90 | 500 gm. | 1105 | 18% |
| 11282SG500 | 877-24-7 | Potassium hydrogen phthalate extrapure | 2917.39.90 | 500 gm. | 936 | 18% |
| 11291LM500 | NA | Potassium hydroxide 40% w/v solution | 2815.20.00 | 500 ml. | 346 | 18% |
| 11289SG500 | 1310-58-3 | Potassium hydroxide flakes extrapure | 2815.20.00 | 500 gm. | 542 | 18% |
| 11289SR500 | | | | 5 kg. | 3541 | 18% |
| 11289S5K50 | | | | 50 kg. | 51226 | 18% |
| 21290SG500 | 1310-58-3 | Potassium hydroxide pellets AR | 2815.20.00 | 500 gm. | 789 | 18% |
| 21290SR500 | | | | 5 kg. | 5488 | 18% |
| 21290S5K50 | | | | 50 kg. | 55076 | 18% |
| 11290SG500 | 1310-58-3 | Potassium hydroxide pellets extrapure | 2815.20.00 | 500 gm. | 639 | 18% |
| 11290SR500 | | | | 5 kg. | 4576 | 18% |
| 11290S5K50 | | | | 50 kg. | 52500 | 18% |
| 21300SG100 | 7758-05-6 | Potassium iodate AR | 2829.90.30 | 100 gm. | 2594 | 18% |
| 11300SG100 | 7758-05-6 | Potassium iodate extrapure | 2829.90.30 | 100 gm. | 2220 | 18% |
| 11300SG250 | | | | 250 gm. | 5033 | 18% |
| 21310SG100 | 7681-11-0 | Potassium iodide AR | 2827.60.10 | 100 gm. | 3237 | 18% |
| 21310SG250 | | | | 250 gm. | 7180 | 18% |
| 21310SG500 | | | | 500 gm. | 12476 | 18% |
| 11310SG100 | 7681-11-0 | Potassium iodide extrapure | 2827.60.10 | 100 gm. | 3025 | 18% |
| 11310SG250 | | | | 250 gm. | 6710 | 18% |
| 11310SG500 | | | | 500 gm. | 11660 | 18% |
| 11315SG500 | 4429-42-9 | Potassium metabisulphite extrapure | 2832.20.10 | 500 gm. | 727 | 18% |
| 21320SG500 | 7757-79-1 | Potassium nitrate AR | 2834.21.00 | 500 gm. | 737 | 18% |
| 11320SG500 | 7757-79-1 | Potassium nitrate extrapure | 2834.21.00 | 500 gm. | 666 | 18% |
| 11325SG500 | 6487-48-5 | Potassium oxalate extrapure | 2917.11.90 | 500 gm. | 822 | 18% |
| 11331LM500 | NA | Potassium permanganate 0.1N solution | 2841.61.00 | 500 ml. | 674 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 71330SG500 | 7722-64-7 | Potassium permanganate crystals AR/ACS | 2841.61.00 | 500 gm. | 986 | 18% |
| 11330SIK50 | 7722-64-7 | Potassium permanganate crystals extrapure (1X50) | | 50 kg. | 76976 | 18% |
| 11330S5K50 | 7722-64-7 | Potassium permanganate crystals extrapure (5X10) | 2841.61.00 | 50 kg. | 72450 | 18% |
| 11330SG500 | 7722-64-7 | Potassium permanganate crystals extrapure | 2841.61.00 | 500 gm. | 825 | 18% |
| 11330SR500 | | | | 5 kg. | 7849 | 18% |
| 11335S5K50 | 7727-21-1 | Potassium persulphate extrapure (5X10) | | 50 kg. | 68054 | 18% |
| 11335SG500 | 7727-21-1 | Potassium persulphate extrapure | 2833.40.00 | 500 gm. | 877 | 18% |
| 11340SG500 | 6381-59-5 | Potassium sodium(+)-tartrate extrapure | 2918.13.90 | 500 gm. | 1464 | 18% |
| 71340SG500 | 6381-59-5 | Potassium sodium tartarate AR/ACS | 2918.13.90 | 500 gm. | 1627 | 18% |
| 11350SR100 | 24634-61-5 | Potassium sorbate extrapure | 2916.19.30 | 1 kg. | 3211 | 18% |
| 11350S5K25 | | | | 25 kg. | 66064 | 18% |
| 71360SG500 | 7778-80-5 | Potassium sulphate AR/ACS | 3104.30.00 | 500 gm. | 795 | 18% |
| 11360SG500 | 7778-80-5 | Potassium sulphate extrapure | 3104.30.00 | 500 gm. | 570 | 18% |
| 11360SR500 | | | | 5 kg. | 4768 | 18% |
| 11360S5K50 | | | | 50 kg. | 66721 | 18% |
| 11233SG250 | 865-47-4 | Potassium tert-Butoxide for synthesis | 2905.10.00 | 250 gm. | 2503 | 18% |
| 71370SG500 | 333-20-0 | Potassium thiocyanate AR/ACS | 2833.30.10 | 500 gm. | 1752 | 18% |
| 11370SG500 | 333-20-0 | Potassium thiocyanate extrapure | 2833.30.10 | 500 gm. | 1529 | 18% |
| 71380LM500 | 71-23-8 | Propan-1-ol AR/ACS (n-Propanol) | 2905.12.10 | 500 ml. | 778 | 18% |
| 71380LC250 | 71-23-8 | Propan-1-ol AR/ACS | 2905.12.10 | 2.5 ltr. | 3115 | 18% |
| 11380LM500 | 71-23-8 | Propan-1-ol extrapure | 2905.12.10 | 500 ml. | 637 | 18% |
| 11380LC250 | | | | 2.5 ltr. | 2611 | 18% |
| 31380LC100 | 71-23-8 | Propan-1-ol HPLC & Spectroscopy | 2905.12.10 | 1 ltr. | 4346 | 18% |
| 12390LM500 | 67-63-0 | Propan-2-ol 70% v/v | 2905.12.20 | 500 ml. | 707 | 18% |
| 12390LC250 | | | | 2.5 ltr. | 2188 | 18% |
| 12390LL025 | | | | 25 ltr. | 21875 | 18% |
| 12390GL025 | NA | Propan-2-ol 70% v/v- GI drum | 2905.12.20 | 25 ltr. | 21875 | 18% |
| 21391LM500 | 67-63-0 | Propan-2-ol AR DRY SOLVENT | 2905.12.20 | 500 ml. | 975 | 18% |
| 71390LL025 | 67-63-0 | Propan-2-ol AR/ACS | 2905.12.20 | 25 ltr. | 16765 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|--------------------------------------|------------|----------|--------|-----|
| 71390LM500 | 67-63-0 | Propan-2-ol AR/ACS (IPA) | 2905.12.20 | 500 ml. | 882 | 18% |
| 71390LC250 | | | | 2.5 ltr. | 2470 | 18% |
| 71390LC500 | | | | 5 ltr. | 3743 | 18% |
| 41390LC250 | 67-63-0 | Propan-2-ol EL Grade | 2905.12.20 | 2.5 ltr. | 2799 | 18% |
| 11390LM500 | 67-63-0 | Propan-2-ol extrapure | 2905.12.20 | 500 ml. | 670 | 18% |
| 11390LC250 | | | | 2.5 ltr. | 2164 | 18% |
| 11390LC500 | | | | 5 ltr. | 3175 | 18% |
| 11390LL025 | | | | 25 ltr. | 17150 | 18% |
| 11390GL025 | 67-63-0 | Propan-2-ol extrapure-GI drum | 2905.12.20 | 25 ltr. | 21000 | 18% |
| 11390LL200 | 67-63-0 | Propan-2-ol extrapure - MS drum | 2905.12.20 | 200 ltr. | 134400 | 18% |
| 11390PL200 | 67-63-0 | Propan-2-ol extrapure - Plastic drum | 2905.12.20 | 200 ltr. | 134400 | 18% |
| 66390LM250 | 67-63-0 | Propan-2-ol GC grade | 2905.12.20 | 250 ml. | 1021 | 18% |
| 31390LC100 | 67-63-0 | Propan-2-ol HPLC & Spectroscopy | 2905.12.20 | 1 ltr. | 1392 | 18% |
| 31390LC250 | | | | 2.5 ltr. | 2665 | 18% |
| 11397LM500 | 108-32-7 | Propylene carbonate extrapure | 2920.90.10 | 500 ml. | 2077 | 18% |
| 11397LC250 | | | | 2.5 ltr. | 7782 | 18% |
| 71400LM500 | 57-55-6 | Propylene glycol AR/ACS | 2905.32.00 | 500 ml. | 883 | 18% |
| 71400LC250 | | | | 2.5 ltr. | 3609 | 18% |
| 71400LL025 | | | | 25 ltr. | 32410 | 18% |
| 11400LM500 | 57-55-6 | Propylene glycol extrapure | 2905.32.00 | 500 ml. | 847 | 18% |
| 11400LC250 | | | | 2.5 ltr. | 3191 | 18% |
| 11400LL025 | | | | 25 ltr. | 31990 | 18% |
| 21421LM500 | 110-86-1 | Pyridine AR DRY SOLVENT | 2933.31.00 | 500 ml. | 1673 | 18% |
| 21421LC250 | | | | 2.5 ltr. | 7360 | 18% |
| 71420LM500 | 110-86-1 | Pyridine AR/ACS | 2933.31.00 | 500 ml. | 1580 | 18% |
| 71420LC250 | | | | 2.5 ltr. | 7118 | 18% |
| 11420LM500 | 110-86-1 | Pyridine extrapure | 2933.31.00 | 500 ml. | 1479 | 18% |
| 11420LC250 | | | | 2.5 ltr. | 5550 | 18% |
| 31420LC100 | 110-86-1 | Pyridine HPLC & Spectroscopy | 2933.31.00 | 1 ltr. | 3008 | 18% |
| 21430SG100 | 87-66-1 | Pyrogallol AR | 2903.59.00 | 100 gm. | 4057 | 18% |
| 21430SG500 | | | | 500 gm. | 14387 | 18% |

(*till stock lasts)

Q

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|---------|---------------------|------------|---------|------|-----|
| 11435LM500 | 91-22-5 | Quinoline extrapure | 2933.39.19 | 500 ml. | 2973 | 18% |

R

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|----------|---|------------|---------|------|-----|
| 71440SG250 | 108-46-3 | Resorcinol (1,3-Benzenediol) AR/ ACS reagent | 2907.21.00 | 250 gm. | 4143 | 18% |
| 11440SG250 | 108-46-3 | Resorcinol extrapure | 2907.21.00 | 250 gm. | 3252 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-------------|--|------------|----------|-------|-----|
| 11441SG025 | 477-73-6 | Safranin stain M.S. (C.I.50240) for microscopy | 3204.13.39 | 25 gm. | 1106 | 18% |
| 11442LM125 | | Safranin stain solution M.S. for microscopy | 3204.13.39 | 125 ml. | 381 | 18% |
| 11436LM100 | 90-02-8 | Salicylaldehyde 99% for synthesis | 2912.49.90 | 100 ml. | 2591 | 18% |
| 11436LM500 | | | | 500 ml. | 8076 | 18% |
| 11443SG500 | 69-72-7 | Salicylic acid for synthesis | 2918.21.10 | 500 gm. | 1035 | 18% |
| 11444LM500 | | Schiff's reagent | 3822.00.90 | 500 ml. | 774 | 12% |
| 11445SG100 | 7446-08-4 | Selenium Dioxide extrapure | 2811.29.90 | 100 gm. | 9522 | 18% |
| 21452SG025 | 7782-49-2 | Selenium metal powder 99.5% AR | 2804.90.00 | 25 gm. | 2871 | 18% |
| 21452SG100 | | | | 100 gm. | 8834 | 18% |
| 11448SG500 | 63231-67-4 | Silica gel 100-200 mesh | 3824.90.25 | 500 gm. | 1178 | 18% |
| 11448S5K25 | | | | 25 kg. | 36225 | 18% |
| 11457SG500 | 63231-67-4 | Silica gel 230-400 mesh | 3824.90.25 | 500 gm. | 1484 | 18% |
| 11457S5K25 | | | | 25 kg. | 44706 | 18% |
| 11450SG500 | 63231-67-4 | Silica gel 5-8 mesh | 3824.90.25 | 500 gm. | 736 | 18% |
| 11446SG500 | 63231-67-4 | Silica gel 60-120 mesh | 3824.90.25 | 500 gm. | 1097 | 18% |
| 11446S5K25 | | | | 25 kg. | 33600 | 18% |
| 11453SG500 | 112926-00-8 | Silica gel G for TLC | 3824.90.25 | 500 gm. | 1106 | 18% |
| 11454SG500 | 63231-67-4 | Silica gel GF 254 for TLC | 3824.90.25 | 500 gm. | 6065 | 18% |
| 11455SG500 | 63231-67-4 | Silica gel H for TLC (without binder) | 3824.90.25 | 500 gm. | 1129 | 18% |
| 11456SG500 | 63231-67-4 | Silica gel HF 254 for TLC | 3824.90.25 | 500 gm. | 5551 | 18% |
| 11460SG050 | | Silicon grease for high vacuum | 3910.00.90 | 50 gm. | 534 | 18% |
| 11461LM250 | 63148-62-9 | Silicon oil | 3910.00.90 | 250 ml. | 1367 | 18% |
| 11461LC250 | | | | 2.5 ltr. | 10387 | 18% |
| 11462SG100 | 12027-43-9 | Silicotungstic acid | 2811.19.00 | 100 gm. | 5128 | 18% |
| 11463SG025 | 563-63-3 | Silver acetate extrapure | 2843.29.00 | 25 gm. | 7697 | 18% |
| 11464SG025 | 7785-23-1 | Silver bromide extrapure | 2843.29.00 | 25 gm. | 9836 | 18% |
| 11465SG025 | 534-16-7 | Silver carbonate extrapure | 2843.29.00 | 25 gm. | 8396 | 18% |
| 11466SG025 | 7783-90-6 | Silver chloride extrapure | 2843.29.00 | 25 gm. | 8228 | 18% |
| 21467SG025 | 1470-61-7 | Silver diethyl dithiocarbamate AR | 2843.29.00 | 25 gm. | 7868 | 18% |
| 11471LM500 | 7761-88-8 | Silver nitrate 0.02 N solution | 2843.21.00 | 500 ml. | 1134 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|-------------------------------------|------------|---------|--------|-----|
| 71470SG010 | 7761-88-8 | Silver nitrate AR/ACS | 2843.21.00 | 10 gm. | 3069 | 18% |
| 71470SG025 | | | | 25 gm. | 7085 | 18% |
| 71470SG100 | | | | 100 gm. | 26060 | 18% |
| 11470SG010 | 7761-88-8 | Silver nitrate extrapure | 2843.21.00 | 10 gm. | 3008 | 18% |
| 11470SG025 | | | | 25 gm. | 7006 | 18% |
| 11470SG100 | | | | 100 gm. | 25686 | 18% |
| 11472LM500 | 7761-88-8 | Silver nitrate N/10 solution | 2843.21.00 | 500 ml. | 3420 | 18% |
| 11475SG025 | 20667-12-3 | Silver oxide extrapure | 2843.29.00 | 25 gm. | 12951 | 18% |
| 71480SG025 | 10294-26-5 | Silver sulphate AR/ACS | 2843.29.00 | 25 gm. | 9163 | 18% |
| 71480SR100 | | | | 1 kg. | 380756 | 18% |
| 11480SG025 | 10294-26-5 | Silver sulphate extrapure | 2843.29.00 | 25 gm. | 7985 | 18% |
| 71490SG250 | 127-09-3 | Sodium acetate anhydrous AR/ACS | 2915.29.90 | 250 gm. | 627 | 18% |
| 11490SG500 | 127-09-3 | Sodium acetate anhydrous extrapure | 2915.29.90 | 500 gm. | 605 | 18% |
| 11490SR500 | | | | 5 kg. | 4415 | 18% |
| 11490S5K50 | | | | 50 kg. | 43313 | 18% |
| 31490SG500 | 127-09-3 | Sodium acetate anhydrous HPLC | 2915.29.90 | 500 gm. | 3025 | 18% |
| 71500SG500 | 6131-90-4 | Sodium acetate trihydrate AR/ACS | 2915.29.90 | 500 gm. | 462 | 18% |
| 11500SG500 | 6131-90-4 | Sodium acetate trihydrate extrapure | 2915.29.90 | 500 gm. | 420 | 18% |
| 11500SR500 | | | | 5 kg. | 2651 | 18% |
| 11500S5K50 | | | | 50 kg. | 23713 | 18% |
| 31500SG500 | 6131-90-4 | Sodium acetate trihydrate for HPLC | 2915.29.90 | 500 gm. | 2727 | 18% |
| 11510SG500 | 9005-38-3 | Sodium alginate extrapure | 3913.10.10 | 500 gm. | 2630 | 18% |
| 11510S5K25 | | | | 25 kg. | 144340 | 18% |
| 21513SG250 | 10048-95-0 | Sodium arsenate AR | 2842.90.10 | 250 gm. | 5933 | 18% |
| 21514SG250 | 7784-46-5 | Sodium arsenite AR | 2842.90.10 | 250 gm. | 10813 | 18% |
| 21514SG500 | | | | 500 gm. | 21229 | 18% |
| 71520SG100 | 26628-22-8 | Sodium azide AR/ACS | 2850.00.30 | 100 gm. | 1391 | 18% |
| 71520SG500 | | | | 500 gm. | 4665 | 18% |
| 11520SG100 | 26628-22-8 | Sodium azide extrapure | 2850.00.30 | 100 gm. | 950 | 18% |
| 11520SG500 | | | | 500 gm. | 3811 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|---------|-------|-----|
| 11540SG500 | 532-32-1 | Sodium benzoate extrapure | 2916.31.40 | 500 gm. | 799 | 18% |
| 11540SR500 | | | | 5 kg. | 5873 | 18% |
| 11540S5K50 | | | | 50 kg. | 62745 | 18% |
| 71550SG500 | 144-55-8 | Sodium bicarbonate anhydrous AR/ACS 2836.30.00 | | 500 gm. | 432 | 18% |
| 71550SR500 | | | | 5 kg. | 2664 | 18% |
| 11550SG500 | 144-55-8 | Sodium bicarbonate anhydrous extrapure | 2836.30.00 | 500 gm. | 329 | 18% |
| 11550SR500 | | | | 5 kg. | 2252 | 18% |
| 11550S5K50 | | | | 50 kg. | 17276 | 18% |
| 11551SG500 | 10034-88-5 | Sodium bisulphate monohydrate extrapure | 2833.19.90 | 500 gm. | 607 | 18% |
| 11551SR500 | | | | 5 kg. | 4859 | 18% |
| 11552SG500 | 7631-90-5 | Sodium bisulphite extrapure | 2832.10.20 | 500 gm. | 534 | 18% |
| 11552SR500 | | | | 5 kg. | 4833 | 18% |
| 71553SG100 | 16940-66-2 | Sodium borohydride AR/ACS | 2850.00.10 | 100 gm. | 2868 | 18% |
| 71553SG500 | | | | 500 gm. | 9014 | 18% |
| 11553SG100 | 16940-66-2 | Sodium borohydride extrapure | 2850.00.10 | 100 gm. | 2190 | 18% |
| 11553SG500 | | | | 500 gm. | 9276 | 18% |
| 11554SG500 | 7789-38-0 | Sodium bromate extrapure | 2829.90.90 | 500 gm. | 1502 | 18% |
| 11555SG500 | 7647-15-6 | Sodium bromide extrapure | 2827.51.10 | 500 gm. | 1070 | 18% |
| 71560SG500 | 497-19-8 | Sodium carbonate anhydrous AR/ACS | 2836.20.20 | 500 gm. | 529 | 18% |
| 71560SR500 | | | | 5 kg. | 2835 | 18% |
| 11560SG500 | 497-19-8 | Sodium carbonate anhydrous extrapure | 2836.20.20 | 500 gm. | 417 | 18% |
| 11560SR500 | | | | 5 kg. | 2640 | 18% |
| 11560S5K50 | | | | 50 kg. | 34650 | 18% |
| 71570SG500 | 7647-14-5 | Sodium chloride AR/ACS | 2827.39.90 | 500 gm. | 283 | 18% |
| 71570SR500 | | | | 5 kg. | 1209 | 18% |
| 71570S5K50 | | | | 50 kg. | 14221 | 18% |
| 11570SG500 | 7647-14-5 | Sodium chloride extrapure | 2827.39.90 | 500 gm. | 256 | 18% |
| 11570SR500 | | | | 5 kg. | 1136 | 18% |
| 11570SK050 | | | | 50 kg. | 14126 | 18% |
| 11570S5K50 | | | | 50 kg. | 14126 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|---------|-------|-----|
| 11573SG500 | 7775-11-3 | Sodium chromate extrapure | 2841.50.00 | 500 gm. | 1018 | 18% |
| 11579SG100 | 13600-98-1 | Sodium cobaltinitrite | 2842.90.90 | 100 gm. | 3510 | 18% |
| 11582SG500 | 7789-12-0 | Sodium dichromate dihydrate extrapure | 2841.30.00 | 500 gm. | 1143 | 18% |
| 21587SG100 | 20624-25-3 | Sodium diethyl dithiocarbamate AR | 2930.20.00 | 100 gm. | 793 | 18% |
| 21587SG500 | | | | 500 gm. | 3370 | 18% |
| 31590SG500 | 13472-35-0 | Sodium di-hydrogen o-phosp dihyd HPLC | 2835.22.00 | 500 gm. | 1882 | 18% |
| 21590SG500 | 13472-35-0 | Sodium di-hydrogen o-phosphate dihydrate AR | 2835.22.00 | 500 gm. | 999 | 18% |
| 21590SR500 | | | | 5 kg. | 6412 | 18% |
| 11590SG500 | 13472-35-0 | Sodium di-hydrogen o-phosphate dihydrate extrapure | 2835.22.00 | 500 gm. | 898 | 18% |
| 11590SR500 | | | | 5 kg. | 4930 | 18% |
| 11590S5K50 | | | | 50 kg. | 52063 | 18% |
| 71591SG500 | 10049-21-5 | Sodium di-hydrogen ortho-phosphate monohydrate AR/ACS (Sodium acid phosphate) | 2835.22.00 | 500 gm. | 1162 | 18% |
| 11591SG500 | 10049-21-5 | Sodium di-hydrogen ortho-phosphate monohydrate extrapure | 2835.22.00 | 500 gm. | 833 | 18% |
| 21594SG005 | 6152-67-6 | Sodium diphenylamine sulphonate AR | 2921.11.10 | 5 gm. | 557 | 18% |
| 21594SG025 | | | | 25 gm. | 1915 | 18% |
| 11597SG500 | 7775-14-6 | Sodium dithionite | 2832.10.20 | 500 gm. | 1098 | 18% |
| 21599SG500 | 7681-49-4 | Sodium fluoride AR | 2826.11.00 | 500 gm. | 1420 | 18% |
| 11599SG500 | 7681-49-4 | Sodium fluoride extrapure | 2826.11.00 | 500 gm. | 907 | 18% |
| 11600SG500 | 6035-47-8 | Sodium formaldehyde sulfoxylate dihydrate extrapure | 2831.10.20 | 500 gm. | 1110 | 18% |
| 21601SG500 | 141-53-7 | Sodium formate AR | 2915.12.10 | 500 gm. | 697 | 18% |
| 11601SG500 | 141-53-7 | Sodium formate extrapure | 2915.12.10 | 500 gm. | 609 | 18% |
| 11602SG500 | 10124-56-8 | Sodium hexametaphosphate extrapure | 2835.29.40 | 500 gm. | 664 | 18% |
| 11602SR500 | | | | 5 kg. | 5113 | 18% |
| 11602S5K50 | | | | 50 kg. | 51576 | 18% |
| 11605SG100 | 7646-69-7 | Sodium hydride suspension | 2850.00.10 | 100 gm. | 1491 | 18% |
| 11605SG500 | | | | 500 gm. | 4504 | 18% |
| 11605SR500 | | | | 5 kg. | 6777 | 18% |
| 11618LM500 | 1310-73-2 | Sodium hydroxide 0.1 N solution | 2815.12.00 | 500 ml. | 315 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---|------------|---------|-------|-----|
| 11619SG500 | 1310-73-2 | Sodium hydroxide flakes extrapure | 2815.11.10 | 500 gm. | 358 | 18% |
| 11619SR500 | | | | 5 kg. | 2225 | 18% |
| 11619S5K50 | | | | 50 kg. | 22750 | 18% |
| 21620SG500 | 1310-73-2 | Sodium hydroxide pellets AR | 2815.11.90 | 500 gm. | 491 | 18% |
| 21620SR500 | | | | 5 kg. | 3271 | 18% |
| 21620S5K50 | | | | 50 kg. | 32326 | 18% |
| 11620SG500 | 1310-73-2 | Sodium hydroxide pellets extrapure | 2815.11.90 | 500 gm. | 442 | 18% |
| 11620SR500 | | | | 5 kg. | 2957 | 18% |
| 11620S5K50 | | | | 50 kg. | 30450 | 18% |
| 11621LC100 | 7681-52-9 | Sodium hypochlorite 4% - 6% w/v | 2828.90.19 | 1 ltr. | 512 | 18% |
| 11621LC500 | | | | 5 ltr. | 1042 | 18% |
| 11621LL025 | | | | 25 ltr. | 3969 | 18% |
| 11622SG100 | 7681-82-5 | Sodium iodide extrapure | 2827.60.20 | 100 gm. | 3166 | 18% |
| 11622SG250 | | | | 250 gm. | 6777 | 18% |
| 71625SG100 | 151-21-3 | Sodium lauryl sulphate AR/ACS | 3402.11.90 | 100 gm. | 4012 | 18% |
| 71625SG250 | | | | 250 gm. | 7906 | 18% |
| 71625SG500 | | | | 500 gm. | 13440 | 18% |
| 71625SR100 | | | | 1 kg. | 2034 | 18% |
| 11625SG500 | 151-21-3 | Sodium lauryl sulphate extrapure | 3402.11.90 | 500 gm. | 1265 | 18% |
| 11625SR500 | | | | 5 kg. | 11385 | 18% |
| 21630SG500 | 7681-57-4 | Sodium metabisulphite AR | 2832.20.90 | 500 gm. | 482 | 18% |
| 21630SR500 | | | | 5 kg. | 3364 | 18% |
| 11630SG500 | 7681-57-4 | Sodium metabisulphite extrapure | 2832.20.90 | 500 gm. | 385 | 18% |
| 11630SR500 | | | | 5 kg. | 2540 | 18% |
| 11630S5K50 | | | | 50 kg. | 26250 | 18% |
| 11485SG250 | 7440-23-5 | Sodium metal | 2805.11.00 | 250 gm. | 881 | 5% |
| 11485SG500 | | | | 500 gm. | 1404 | 5% |
| 71635SG100 | 7790-28-5 | Sodium metaperiodate AR/ACS (Sodium periodate) | 2829.90.30 | 100 gm. | 2317 | 18% |
| 11635SG100 | 7790-28-5 | Sodium metaperiodate extrapure | 2829.90.30 | 100 gm. | 1937 | 18% |
| 11635SG500 | | | | 500 gm. | 8064 | 18% |
| 11636SG500 | 124-41-4 | Sodium methoxide extrapure | 2815.11.90 | 500 gm. | 1386 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--------------------------------------|------------|---------|-------|-----|
| 71638SG100 | 10102-40-6 | Sodium molybdate AR/ACS | 2841.70.20 | 100 gm. | 1935 | 18% |
| 71638SG500 | | | | 500 gm. | 7724 | 18% |
| 11638SG100 | 10102-40-6 | Sodium molybdate dihydrate extrapure | 2841.70.20 | 100 gm. | 1743 | 18% |
| 11638SG500 | | | | 500 gm. | 7206 | 18% |
| 71640SG500 | 7631-99-4 | Sodium nitrate AR/ACS | 2834.10.10 | 500 gm. | 509 | 18% |
| 11640SG500 | 7631-99-4 | Sodium nitrate extrapure | 2834.10.10 | 500 gm. | 454 | 18% |
| 11640SR500 | | | | 5 kg. | 3490 | 18% |
| 11640S5K50 | | | | 50 kg. | 51576 | 18% |
| 71650SG500 | 7632-00-0 | Sodium nitrite AR/ACS | 2834.10.10 | 500 gm. | 618 | 18% |
| 71650SK050 | | | | 50 kg. | 51576 | 18% |
| 11650SG500 | 7632-00-0 | Sodium nitrite extrapure | 2834.10.10 | 500 gm. | 571 | 18% |
| 11650SR500 | | | | 5 kg. | 4192 | 18% |
| 11650S5K50 | | | | 50 kg. | 51576 | 18% |
| 11655SG100 | 13755-38-9 | Sodium nitroprusside extrapure | 2837.20.50 | 100 gm. | 936 | 18% |
| 11655SG500 | | | | 500 gm. | 4405 | 18% |
| 21675SG500 | 62-76-0 | Sodium oxalate AR | 2917.11.90 | 500 gm. | 901 | 18% |
| 11675SG500 | 62-76-0 | Sodium oxalate extrapure | 2917.11.90 | 500 gm. | 683 | 18% |
| 12661SG500 | 10486-00-7 | Sodium perborate tetrahydrate pure | 2840.19.90 | 500 gm. | 452 | 18% |
| 12661SR100 | | | | 1 kg. | 794 | 18% |
| 12661S5K25 | | | | 25 kg. | 26509 | 18% |
| 21658SG500 | 1313-60-6 | Sodium peroxide granular,95% AR | 2825.90.00 | 500 gm. | 28389 | 18% |
| 21679SG001 | 523-21-7 | Sodium rhodizonate AR | 2942.00.00 | 1 gm. | 1672 | 18% |
| 21679SG005 | | | | 5 gm. | 6279 | 18% |
| 11677SG500 | 54-21-7 | Sodium salicylate purified | 2915.90.90 | 500 gm. | 1021 | 18% |
| 21682SG100 | 13410-01-0 | Sodium Selenate AR | 2841.90.00 | 100 gm. | 7453 | 18% |
| 21683SG025 | 10102-18-8 | Sodium Selenite AR | 2841.90.00 | 25 gm. | 1941 | 18% |
| 11678SG500 | 10213-79-3 | Sodium silicate (meta) pentahydrate | 2839.11.00 | 500 gm. | 644 | 18% |
| 11681LC500 | NA | Sodium silicate solution EP | 2839.11.00 | 5 ltr. | 1450 | 18% |
| 71680SG500 | 7757-82-6 | Sodium sulphate anhydrous AR/ACS | 2833.11.00 | 500 gm. | 386 | 18% |
| 11680SG500 | 7757-82-6 | Sodium sulphate anhydrous extrapure | 2833.11.00 | 500 gm. | 320 | 18% |
| 11680SR500 | | | | 5 kg. | 2079 | 18% |
| 11680S5K50 | | | | 50 kg. | 52500 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|----------|-------|-----|
| 11685SG500 | 1313-84-4 | Sodium sulphide flakes extrapure | 2830.10.00 | 500 gm. | 454 | 18% |
| 21690SG500 | 7757-83-7 | Sodium sulphite anhydrous AR | 2832.10.00 | 500 gm. | 380 | 18% |
| 11690SG500 | 7757-83-7 | Sodium sulphite anhydrous extrapure | 2832.10.00 | 500 gm. | 346 | 18% |
| 11690SR500 | | | | 5 kg. | 2345 | 18% |
| 11690S5K50 | | | | 50 kg. | 51254 | 18% |
| 11698SG100 | 6106-24-7 | Sodium tartrate dihydrate water standard for Karl Fischer titration | 2918.13.90 | 100 gm.* | 324 | 18% |
| 21708SG010 | 143-66-8 | Sodium tetraphenyl boron AR | 2931.00.90 | 10 gm. | 2487 | 18% |
| 21708SG025 | | | | 25 gm. | 5092 | 18% |
| 11709SG500 | 540-72-7 | Sodium thiocyanate extrapure | 2837.20.90 | 500 gm. | 1200 | 18% |
| 11711LM500 | 10102-17-7 | Sodium thiosulphate 0.1N solution | 2832.30.10 | 500 ml. | 270 | 18% |
| 71710SG500 | 10102-17-7 | Sodium thiosulphate pentahydrate AR/ACS | 2832.30.10 | 500 gm. | 481 | 18% |
| 11710SG500 | 10102-17-7 | Sodium thiosulphate pentahydrate extrapure | 2832.30.10 | 500 gm. | 321 | 18% |
| 11710SR500 | | | | 5 kg. | 2093 | 18% |
| 11710S5K50 | | | | 50 kg. | 51576 | 18% |
| 11715SG500 | 7758-29-4 | Sodium tripolyphosphate extrapure | 2835.31.00 | 500 gm. | 436 | 18% |
| 11717SG100 | 10213-10-2 | Sodium tungstate hydrated extrapure | 2841.80.10 | 100 gm. | 1591 | 18% |
| 11717SG250 | | | | 250 gm. | 3827 | 18% |
| 11717S5K25 | | | | 25 kg. | 26210 | 18% |
| 11720SG025 | 1787-61-7 | Solochrome black-T extrapure | 3204.12.67 | 25 gm. | 617 | 18% |
| 11720SG100 | | | | 100 gm. | 1506 | 18% |
| 11722SG500 | 110-44-1 | Sorbic acid pure | 2916.19.60 | 500 gm. | 2662 | 18% |
| 71730SG100 | 10025-69-1 | Stannous chloride dihydrate AR/ACS | 2827.39.90 | 100 gm. | 1645 | 18% |
| 71730SG250 | | | | 250 gm. | 3514 | 18% |
| 11730SG100 | 10025-69-1 | Stannous chloride dihydrate extrapure | 2827.39.90 | 100 gm. | 1387 | 18% |
| 11730SG250 | | | | 250 gm. | 2461 | 18% |
| 21740SG500 | 9005-84-9 | Starch AR | 3505.10.90 | 500 gm. | 2059 | 18% |
| 71740SG500 | 9005-84-9 | Starch AR/ACS | | 500 gm. | 757 | 18% |
| 11740SG500 | 9005-84-9 | Starch extrapure | 3505.10.90 | 500 gm. | 1366 | 18% |
| 11735BK005 | 9005-25-8 | Starch iodide papers 5BK | 3822.00.90 | 5 bkt. | 261 | 12% |
| 71740SG100 | 9005-84-9 | Starch, soluble AR/ACS | 3505.10.90 | 100 gm. | 637 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|----------------------------------|------------|------------|-------|-----|
| 11737SG500 | 57-11-4 | Stearic acid for synthesis | 2915.70.90 | 500 gm. | 640 | 18% |
| 11741SG500 | 1633-05-2 | Strontium carbonate extrapure | 2836.92.00 | 500 gm. | 1295 | 18% |
| 11742SG500 | 10025-70-4 | Strontium chloride extrapure | 2827.39.30 | 500 gm. | 734 | 18% |
| 11743SG500 | 10042-76-9 | Strontium nitrate extrapure | 2834.29.10 | 500 gm. | 764 | 18% |
| 11744SG500 | 110-15-6 | Succinic acid extrapure | 2917.19.30 | 500 gm. | 1045 | 18% |
| 11744S5K50 | | | | 50 kg. | 52063 | 18% |
| 21746SG500 | 57-50-1 | Sucrose AR | 2940.00.00 | 500 gm. | 871 | 18% |
| 11746SG500 | 57-50-1 | Sucrose extrapure | 2940.00.00 | 500 gm. | 477 | 18% |
| 11745SG500 | 5329-14-6 | Sulphamic acid extrapure | 2935.00.90 | 500 gm. | 545 | 18% |
| 11745SR500 | | | | 5 kg. | 3965 | 18% |
| 11745SK050 | | | | 50 kg. | 35826 | 18% |
| 11745S5K50 | | | | 50 kg. | 37660 | 18% |
| 11753SG500 | 63-74-1 | Sulphanilamide | 2935.00.00 | 500 gm. | 4979 | 18% |
| 21753SG100 | 63-74-1 | Sulphanilamide AR | 2935.00.00 | 100 gm. | 1525 | 18% |
| 21753SG500 | | | | 500 gm. | 6277 | 18% |
| 11757SG500 | 127-57-3 | Sulphanilic acid | 2921.42.34 | 500 gm. | 1596 | 18% |
| 21757SG100 | 127-57-3 | Sulphanilic acid AR | 2921.42.34 | 100 gm. | 1017 | 18% |
| 11751SG500 | 5965-83-3 | 5-Sulphosalicylic acid extrapure | 2918.21.90 | 500 gm. | 2000 | 18% |
| 11759SG500 | 7704-34-9 | Sulphur powder | 2802.00.00 | 500 gm. | 533 | 18% |
| 11761LM500 | 7664-93-9 | Sulphuric acid 0.02N solution | 2807.00.10 | 500 ml. | 309 | 18% |
| 71760L8M50 | 7664-93-9 | Sulphuric acid AR/ACS | 2807.00.10 | 8x500 ml. | 4364 | 18% |
| 71760L2C50 | | | | 2x5 ltr. | 3878 | 18% |
| 71760L4C25 | | | | 4x2.5 ltr. | 6662 | 18% |
| 11760L8M50 | 7664-93-9 | Sulphuric acid extrapure | 2807.00.10 | 8x500 ml. | 3790 | 18% |
| 11760L2C50 | | | | 2x5 ltr. | 3372 | 18% |
| 11760L4C25 | | | | 4x2.5 ltr. | 5586 | 18% |
| 11760L020 | | | | 20 ltr. | 19764 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|----------|--------|-----|
| 70334LM500 | 75-65-0 | tert-Butanol AR/ACS | 2905.19.10 | 500 ml. | 1277 | 18% |
| 30334LC100 | 75-65-0 | tert-Butanol HPLC & Spectroscopy | 2905.19.10 | 1 ltr. | 2809 | 18% |
| 70376LC250 | 1634-04-4 | tert-Butyl methyl ether AR/ACS | 2909.19.00 | 2.5 ltr. | 4303 | 18% |
| 70376LM500 | 1634-04-4 | tert-Butyl methyl ether AR/ACS (Methyl tert-butyl ether) | 2909.19.00 | 500 ml. | 1153 | 18% |
| 10376LM500 | 1634-04-4 | tert-Butyl methyl ether extrapure | 2909.19.00 | 500 ml. | 937 | 18% |
| 10376LC250 | | | | 2.5 ltr. | 3422 | 18% |
| 10376LL025 | | | | 25 ltr. | 28583 | 18% |
| 66376LM250 | 1634-04-4 | Tert-Butyl methyl ether GC Grade | 2909.19.00 | 250 ml. | 1765 | 18% |
| 30376LC100 | 1634-04-4 | tert-Butyl methyl ether HPLC | 2909.19.00 | 1 ltr. | 3842 | 18% |
| 30376LC250 | | | | 2.5 ltr. | 3404 | 18% |
| 21802LM500 | NA | Tetrabutyl ammonium hydroxide 0.1N in methanol AR | 2923.10.90 | 500 ml. | 3505 | 18% |
| 21780SG500 | 1643-19-2 | Tetrabutylammonium bromide AR | 2827.59.90 | 500 gm. | 2523 | 18% |
| 21790SG100 | 32503-27-8 | Tetrabutylammonium hydrogen sulphate AR | 3815.12.90 | 100 gm. | 3081 | 18% |
| 21812LM250 | 2052-49-5 | Tetrabutylammonium hydroxide 40% in Methanol AR | 2923.10.90 | 250 ml. | 9013 | 18% |
| 71800LM500 | 127-18-4 | Tetrachloroethylene AR/ACS | 2903.23.00 | 500 ml. | 2952 | 18% |
| 11800LM500 | 127-18-4 | Tetrachloroethylene extrapure | 2903.23.00 | 500 ml. | 2266 | 18% |
| 11800LC250 | | | | 2.5 ltr. | 9617 | 18% |
| 21813SG005 | 14937-42-9 | Tetradecylammonium bromide, 99% AR for HPLC (Tetrakis (decyl) ammonium bromide) | 2804.70.20 | 5 gm. | 8557 | 18% |
| 21813SG025 | | | | 25 gm. | 20997 | 18% |
| 31810SG025 | 4368-51-8 | Tetraheptylammonium bromide HPLC | 2827.59.90 | 25 gm. | 14042 | 18% |
| 31810SG100 | | | | 100 gm. | 45798 | 18% |
| 21821LM500 | 109-99-9 | Tetrahydrofuran AR DRY SOLVENT | 2932.11.00 | 500 ml. | 1681 | 18% |
| 21821LC250 | | | | 2.5 ltr. | 7089 | 18% |
| 71820LM500 | 109-99-9 | Tetrahydrofuran AR/ACS | 2932.11.00 | 500 ml. | 1531 | 18% |
| 71820LC250 | | | | 2.5 ltr. | 6783 | 18% |
| 71820LL200 | | | | 200 ltr. | 259826 | 18% |
| 41820LL025 | 109-99-9 | Tetrahydrofuran DRY-S | 2932.11.00 | 25 ltr. | 35408 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 11820LM500 | 109-99-9 | Tetrahydrofuran extrapure | 2932.11.00 | 500 ml. | 1464 | 18% |
| 11820LC250 | | | | 2.5 ltr. | 6248 | 18% |
| 11820LL025 | | | | 25 ltr. | 35408 | 18% |
| 66820LM250 | 109-99-9 | Tetrahydrofuran GC Grade | 2932.11.00 | 250 ml. | 1176 | 18% |
| 31820LM500 | 109-99-9 | Tetrahydrofuran HPLC & Spectroscopy | 2932.11.00 | 500 ml. | 1617 | 18% |
| 31820LC100 | | | | 1 ltr. | 3072 | 18% |
| 31820LC250 | | | | 2.5 ltr. | 7220 | 18% |
| 32820LM500 | 109-99-9 | Tetrahydrofuran HPLC & Spectroscopy (without stabilizer) | 2932.11.00 | 500 ml. | 2036 | 18% |
| 11840SG025 | 54-64-8 | Thimerosal extrapure | 2930.20.00 | 25 gm. | 18690 | 18% |
| 11840SG100 | | | | 100 gm. | 68997 | 18% |
| 11842SG100 | 62-55-5 | Thioacetamide purified | 2930.90.70 | 100 gm. | 2304 | 18% |
| 11841LM500 | 68-11-1 | Thioglycollic acid 80% extrapure | 2930.90.99 | 500 ml. | 2688 | 18% |
| 11843LM500 | 7719-09-7 | Thionyl chloride for synthesis | 2812.10.47 | 500 ml. | 850 | 18% |
| 71845SG500 | 62-56-6 | Thiourea AR/ACS | 2930.90.10 | 500 gm. | 1766 | 18% |
| 11845S5K25 | 62-56-6 | Thiourea extrapure (5X5) | 2930.90.10 | 25 kg. | 49786 | 18% |
| 11845SG500 | 62-56-6 | Thiourea extrapure | 2930.90.10 | 500 gm. | 1604 | 18% |
| 11846SG005 | 76-61-9 | Thymol blue indicator | 3822.00.90 | 5 gm. | 450 | 12% |
| 11846SG025 | | | | 25 gm. | 1519 | 12% |
| 11851SG100 | 89-83-8 | Thymol crystals 99% | 2907.19.30 | 100 gm. | 1180 | 18% |
| 11851SG250 | | | | 250 gm. | 2409 | 18% |
| 21850SG005 | 125-20-2 | Thymolphthalein AR | 2932.29.90 | 5 gm. | 576 | 18% |
| 21850SG025 | | | | 25 gm. | 1638 | 18% |
| 41850LM125 | 125-20-2 | Thymolphthalein solution | 2932.29.90 | 125 ml. | 304 | 18% |
| 21855SG100 | 7440-31-5 | Tin (metal) granular AR | 8007.00.90 | 100 gm. | 2279 | 18% |
| 21855SG500 | | | | 500 gm. | 10051 | 18% |
| 11855SG100 | 7440-31-5 | Tin (metal) granular extrapure | 8007.00.90 | 100 gm. | 1729 | 18% |
| 11855SG500 | | | | 500 gm. | 7556 | 18% |
| 10148LM500 | | TISAB III solution for fluoride determination 500ML | 3822.90.90 | 500 ml. | 12000 | 18% |
| 11856S5K50 | 13463-67-7 | Titanium dioxide (5X10) | 2823.00.10 | 50 kg. | 74550 | 18% |
| 11856SG500 | 13463-67-7 | Titanium dioxide | 2823.00.10 | 500 gm. | 751 | 18% |
| 21856SG500 | 13463-67-7 | Titanium dioxide AR | 2823.00.10 | 500 gm. | 1796 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|--|------------|----------|--------|-----|
| 11861LM500 | 108-88-3 | Toluene (Rectified) extrapure | 2902.30.00 | 500 ml. | 542 | 18% |
| 11861LC250 | | | | 2.5 ltr. | 1838 | 18% |
| 11861LL025 | | | | 25 ltr. | 31628 | 18% |
| 11860LM500 | 108-88-3 | Toluene (Sulphur / Thiophene free) extrapure | 2902.30.00 | 500 ml. | 527 | 18% |
| 11860LC250 | | | | 2.5 ltr. | 1652 | 18% |
| 11860LL025 | | | | 25 ltr. | 21583 | 18% |
| 11860LL200 | | | | 200 ltr. | 149800 | 18% |
| 21861LM500 | 108-88-3 | Toluene AR DRY SOLVENT | 2902.30.00 | 500 ml. | 626 | 18% |
| 71860LM500 | 108-88-3 | Toluene AR/ACS | 2902.30.00 | 500 ml. | 594 | 18% |
| 71860LC250 | | | | 2.5 ltr. | 1946 | 18% |
| 71860LL025 | | | | 25 ltr. | 21583 | 18% |
| 41860LL025 | 108-88-3 | Toluene DRY-S | | 25 ltr. | 21583 | 18% |
| 66860LM250 | 108-88-3 | Toluene GC Grade | 2902.30.00 | 250 ml. | 975 | 18% |
| 31860LC100 | 108-88-3 | Toluene HPLC & Spectroscopy | 2902.30.00 | 1 ltr. | 1389 | 18% |
| 21863LM500 | 102-76-1 | Triacetin AR | 2917.19.90 | 500 ml. | 1288 | 18% |
| 21863LL025 | | | | 25 ltr. | 101456 | 18% |
| 10098SG500 | 3458-72-8 | tri-Ammonium citrate extrapure | 2918.15.00 | 500 gm. | 1936 | 18% |
| 10444SG500 | 7758-87-4 | Tricalcium phosphate extrapure | 2835.26.00 | 500 gm. | 660 | 18% |
| 11864SG100 | 76-03-9 | Trichloroacetic acid extrapure | 2915.40.30 | 100 gm. | 307 | 18% |
| 11864SG500 | | | | 500 gm. | 905 | 18% |
| 71865LC250 | 79-01-6 | Trichloroethylene AR/ACS | 2903.22.00 | 2.5 ltr. | 4461 | 18% |
| 71865LM500 | 79-01-6 | Trichloroethylene AR/ACS | 2903.22.00 | 500 ml. | 1234 | 18% |
| 11865LM500 | 79-01-6 | Trichloroethylene extrapure | 2903.22.00 | 500 ml. | 1149 | 18% |
| 11865LC250 | | | | 2.5 ltr. | 4251 | 18% |
| 11865LL025 | | | | 25 ltr. | 12250 | 18% |
| 21870LM500 | 102-71-6 | Triethanolamine AR | 2922.13.00 | 500 ml. | 1164 | 18% |
| 11870LM500 | 102-71-6 | Triethanolamine extrapure | 2922.13.00 | 500 ml. | 928 | 18% |
| 11870LC250 | | | | 2.5 ltr. | 3909 | 18% |
| 21873LM500 | 77-93-0 | Triethyl citrate AR | 2918.15.90 | 500 ml. | 4013 | 18% |
| 21869LM500 | 121-44-8 | Triethylamine DRY SOLVENT | 2921.11.10 | 500 ml. | 1978 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|--|------------|----------|-------|-----|
| 11871LM500 | 121-44-8 | Triethylamine for synthesis | 2921.11.10 | 500 ml. | 1611 | 18% |
| 11871LC250 | | | | 2.5 ltr. | 6093 | 18% |
| 31871LC100 | 121-44-8 | Triethylamine HPLC & Spectroscopy | 2921.11.10 | 1 ltr. | 6526 | 18% |
| 21874LM100 | 76-05-1 | Trifluoroacetic acid AR | 2826.90.00 | 100 ml. | 2683 | 18% |
| 11874LM100 | 76-05-1 | Trifluoroacetic acid extrapure | 2826.90.00 | 100 ml. | 2266 | 18% |
| 11874LM500 | | | | 500 ml. | 6168 | 18% |
| 31874LM100 | 76-05-1 | Trifluoroacetic acid HPLC & Spectroscopy | 2826.90.00 | 100 ml. | 2981 | 18% |
| 31874LM500 | | | | 500 ml. | 8134 | 18% |
| 31874LC100 | | | | 1 ltr. | 13414 | 18% |
| 21255SG500 | 6100-05-6 | tri-Potassium citrate AR | 2918.15.10 | 500 gm. | 758 | 18% |
| 11255SG500 | 6100-05-6 | tri-Potassium citrate extrapure | 2918.15.10 | 500 gm. | 717 | 18% |
| 21890SG100 | 77-86-1 | Tris-Buffer AR | 2914.70.90 | 100 gm. | 1673 | 18% |
| 21890SG500 | | | | 500 gm. | 5955 | 18% |
| 11890SG100 | 77-86-1 | Tris-Buffer extrapure | 2914.70.90 | 100 gm. | 1301 | 18% |
| 11890SG500 | | | | 500 gm. | 4434 | 18% |
| 71580SG500 | 6132-04-3 | tri-Sodium citrate dihydrate AR/ACS | 2918.15.20 | 500 gm. | 794 | 18% |
| 11580SG500 | 6132-04-3 | tri-Sodium citrate dihydrate extrapure | 2918.15.20 | 500 gm. | 566 | 18% |
| 11580S5K50 | | | | 50 kg. | 46326 | 18% |
| 21670SG500 | 10101-89-0 | tri-Sodium ortho-phosphate dodecahydrate AR | 2835.22.00 | 500 gm. | 547 | 18% |
| 21670SR500 | | | | 5 kg. | 4467 | 18% |
| 21670S5K50 | | | | 50 kg. | 50313 | 18% |
| 11670SG500 | 10101-89-0 | tri-Sodium ortho-phosphate dodecahydrate extrapure | 2835.22.00 | 500 gm. | 466 | 18% |
| 11670SR500 | | | | 5 kg. | 3476 | 18% |
| 11670S5K50 | | | | 50 kg. | 40138 | 18% |
| 11891LM500 | 9002-93-1 | Triton X-100 for biochemistry | 3402.13.00 | 500 ml. | 1860 | 18% |
| 11898LM500 | 9005-65-6 | Tween 80 for synthesis | 3402.90.99 | 500 ml. | 1545 | 18% |
| 11898LC250 | | | | 2.5 ltr. | 7106 | 18% |
| 1189811025 | | | | 25 ltr. | 11451 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|---------|---|------------|---------|-----|-----|
| 11900LM125 | NA | Universal indicator solution (pH 4-11) incl.colourscale | 3822.00.90 | 125 ml. | 194 | 12% |
| 11900LM500 | | | | 500 ml. | 392 | 12% |
| 71920SG500 | 57-13-6 | Urea AR/ACS | 3102.10.00 | 500 gm. | 685 | 18% |
| 11920SG500 | 57-13-6 | Urea extrapure | 3102.10.00 | 500 gm. | 500 | 18% |

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---------------------------|------------|---------|-----|-----|
| 31930LC100 | 7732-18-5 | Water HPLC & Spectroscopy | 2201.10.10 | 1 ltr. | 579 | 18% |

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---------------------|------------|----------|-------|-----|
| 21940LM500 | 1330-20-7 | Xylene AR | 2902.44.00 | 500 ml. | 632 | 18% |
| 21940LC250 | | | | 2.5 ltr. | 2199 | 18% |
| 21940LL025 | | | | 25 ltr. | 10785 | 18% |
| 11940LM500 | 1330-20-7 | Xylene extrapure | 2902.44.00 | 500 ml. | 570 | 18% |
| 11940LC250 | | | | 2.5 ltr. | 1940 | 18% |
| 11940LL025 | | | | 25 ltr. | 10675 | 18% |
| 11940LL200 | | | | 200 ltr. | 75354 | 18% |
| 21942SG005 | 3618-43-7 | Xylenol orange AR | 2927.00.90 | 5 gm. | 1089 | 18% |
| 21942SG010 | | | | 10 gm. | 1863 | 18% |

(*till stock lasts)

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|------------|---|------------|---------|-------|-----|
| 21970SG500 | 7440-66-6 | Zinc (metal) dust AR | 7903.10.00 | 500 gm. | 1148 | 18% |
| 11970SG500 | 7440-66-6 | Zinc (metal) dust extrapure | 7903.10.00 | 500 gm. | 1080 | 18% |
| 11970SR500 | | | | 5 kg. | 9501 | 18% |
| 11970S5K50 | | | | 50 kg. | 49329 | 18% |
| 11971SG500 | 7440-66-6 | Zinc (metal) granular (As free) extrapure | 7903.10.00 | 500 gm. | 1584 | 18% |
| 11972SG500 | 5970-45-6 | Zinc acetate extrapure | 2915.29.90 | 500 gm. | 801 | 18% |
| 11972S5K50 | | | | 50 kg. | 52063 | 18% |
| 11975SG500 | 12539-71-8 | Zinc carbonate basic extrapure | 2836.99.90 | 500 gm. | 1912 | 18% |
| 11980SG500 | 7646-85-7 | Zinc chloride anhydrous extrapure | 2827.20.00 | 500 gm. | 681 | 18% |
| 11980S5K50 | | | | 50 kg. | 52063 | 18% |
| 71980SG500 | 7646-85-7 | Zinc chloride AR/ACS | 2827.20.00 | 500 gm. | 1321 | 18% |
| 21990SG500 | 1314-13-2 | Zinc oxide AR | 2817.00.10 | 500 gm. | 1017 | 18% |
| 11990SG500 | 1314-13-2 | Zinc oxide extrapure | 2817.00.10 | 500 gm. | 934 | 18% |
| 11990SR500 | | | | 5 kg. | 7756 | 18% |
| 11990S5K50 | | | | 50 kg. | 92078 | 18% |
| 31990S5K25 | 1314-13-2 | Zinc oxide Special for electroplating | | 25 kg. | 27305 | 18% |
| 72000SG500 | 7446-20-0 | Zinc sulphate heptahydrate AR/ACS | 2833.21.00 | 500 gm. | 758 | 18% |
| 12000SG500 | 7446-20-0 | Zinc sulphate heptahydrate extrapure | 2833.21.00 | 500 gm. | 486 | 18% |
| 12000S5K50 | | | | 50 kg. | 39410 | 18% |
| 12010SG500 | 7446-19-7 | Zinc sulphate monohydrate extrapure | 2833.21.00 | 500 gm. | 608 | 18% |

Special Reagents

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|---|------------|---------|-------|-----|
| 61030LM500 | 7697-37-2 | Nitric acid TracepurA (For trace metal analysis) | 28.08.0010 | 500 ml. | 17545 | 18% |
| 61030LC100 | 7697-37-2 | Nitric acid TracepurA (For trace metal analysis) | 28.08.0010 | 1 Ltr. | 24079 | 18% |
| 60770LM500 | 7647-01-0 | Hydrochloric acid TracepurA (For trace metal analysis) | 28.06.1000 | 500 ml. | 16940 | 18% |
| 60770LC100 | 7647-01-0 | Hydrochloric acid TracepurA (For trace metal analysis) | 28.06.1000 | 1 Ltr. | 22990 | 18% |
| 60776LM500 | 7722-84-1 | Hydrogen peroxide solution 30% w/v TracepurA | 2847.00.00 | 500 ml. | 30782 | 18% |
| 60760LM500 | 7664-93-9 | Sulphuric acid TracepurA (For trace metal analysis) | 28.07.0010 | 500 ml. | 17727 | 18% |
| 60760LC100 | 7664-93-9 | Sulphuric acid TracepurA (For trace metal analysis) | 28.07.0010 | 1 Ltr. | 23111 | 18% |

Part-II

High Purity Solvents

High Purity Solvents

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|----------|-----------|-----------------------------------|------------|---------|-------|-----|
| 7518 | 75-05-8 | Acetonitrile Supra Gradient Grade | 29.26.9001 | 4 Ltr. | 11955 | 18% |
| 3040 | 75-05-8 | Acetonitrile LCMS Grade | 29.26.9001 | 4 Ltr. | 13964 | 18% |
| 7517 | 67-56-1 | Methanol Supra Gradient Grade | 29.05.1100 | 4 Ltr. | 4508 | 18% |
| 3041 | 67-56-1 | Methanol LCMS Grade | 29.05.1100 | 4 Ltr. | 10835 | 18% |
| 3042 | 7732-18-5 | Water LCMS Grade | 28.53.0010 | 4 Ltr. | 7777 | 18% |

(*till stock lasts)

Part-III

Inorganic Ventures' Product Range

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---|----------------------------|------------|---------|-------|-----|
| ICP & ICP-MS Single-element Standards - 1 ppm | | | | | |
| MSHG-1PPM-125ML | 1 ppm Mercury for ICP | 38.22.0090 | 125 ml. | 11879 | 12% |
| ICP & ICP-MS Single-element Standards - 10 ppm | | | | | |
| MSAG-10PPM-125ML | 10ug/mL SILVER 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSAL-10PPM-125ML | 10ug/mL ALUMINUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSAS-10PPM-125ML | 10ug/mL ARSENIC 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSAU-10PPM-125ML | 10ug/mL GOLD HCl 125mL | 38.22.0090 | 125 ml. | 36336 | 12% |
| MSB-10PPM-125ML | 10ug/mL BORON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSBA-10PPM-125ML | 10ug/mL BARIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSBE-10PPM-125ML | 10ug/mL BERYLLIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSBI-10PPM-125ML | 10ug/mL BISMUTH 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCA-10PPM-125ML | 10ug/mL CALCIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCD-10PPM-125ML | 10ug/mL CADMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCE-10PPM-125ML | 10ug/mL CERIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCO-10PPM-125ML | 10ug/mL COBALT 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCS-10PPM-125ML | 10ug/mL CESIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCU-10PPM-125ML | 10ug/mL COPPER 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSFE-10PPM-125ML | 10ug/mL IRON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSGE-10PPM-125ML | 10ug/mL GERMANIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHF-10PPM-125ML | 10ug/mL HAFNIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHG-10PPM-125ML | 10ug/mL MERCURY HCl 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHGN-10PPM-125ML | 10ug/mL MERCURY HNO3 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHO-10PPM-125ML | 10ug/mL HOLMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSIN-10PPM-125ML | 10ug/mL INDIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSK-10PPM-125ML | 10ug/mL POTASSIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSLI-10PPM-125ML | 10ug/mL LITHIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSMG-10PPM-125ML | 10ug/mL MAGNESIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSMN-10PPM-125ML | 10ug/mL MANGANESE 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSMO-10PPM-125ML | 10ug/mL MOLYBDENUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSNA-10PPM-125ML | 10ug/mL SODIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSNI-10PPM-125ML | 10ug/mL NICKEL 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSOS-10PPM-125ML | 10ug/mL OSMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSP-10PPM-125ML | 10ug/mL PHOSPHORUS 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|--|----------------------------|------------|---------|-------|-----|
| MSPB-10PPM-125ML | 10ug/mL LEAD 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSPT-10PPM-125ML | 10ug/mL PLATINUM 125mL | 38.22.0090 | 125 ml. | 30280 | 12% |
| MSRH-10PPM-125ML | 10ug/mL RHODIUM HCl 125mL | 38.22.0090 | 125 ml. | 39364 | 12% |
| MSRH-10PPM-125ML | 10ug/mL RHODIUM HNO3 125mL | 38.22.0090 | 125 ml. | 39364 | 12% |
| MSS-10PPM-125ML | 10ug/mL SULFUR 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSB-10PPM-125ML | 10ug/mL ANTIMONY 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSC-10PPM-125ML | 10ug/mL SCANDIUM 125mL | 38.22.0090 | 125 ml. | 27951 | 12% |
| MSSE-10PPM-125ML | 10ug/mL SELENIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSI-10PPM-125ML | 10ug/mL SILICON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSN-10PPM-125ML | 10ug/mL TIN 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSR-10PPM-125ML | 10ug/mL STRONTIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTB-10PPM-125ML | 10ug/mL TERBIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTEN-10PPM-125ML | 10ug/mL TELLURIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTI-10PPM-125ML | 10ug/mL TITANIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTL-10PPM-125ML | 10ug/mL THALLIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSV-10PPM-125ML | 10ug/mL VANADIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSW-10PPM-125ML | 10ug/mL TUNGSTEN 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSY-10PPM-125ML | 10ug/mL YTTRIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSZN-10PPM-125ML | 10ug/mL ZINC 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| ICP & ICP-MS Single-element Standards - 100 ppm | | | | | |
| MSAG-100PPM-125ML | 100ug/mL SILVER 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSAL-100PPM-125ML | 100ug/mL ALUMINUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSAS-100PPM-125ML | 100ug/mL ARSENIC 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSAU-100PPM-125ML | 100ug/mL GOLD HCl 125mL | 38.22.0090 | 125 ml. | 30979 | 12% |
| MSB-100PPM-125ML | 100ug/mL BORON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSBA-100PPM-125ML | 100ug/mL BARIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSBE-100PPM-125ML | 100ug/mL BERYLLIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSBI-100PPM-125ML | 100ug/mL BISMUTH 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCA-100PPM-125ML | 100ug/mL CALCIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCD-100PPM-125ML | 100ug/mL CADMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCE-100PPM-125ML | 100ug/mL CERIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCO-100PPM-125ML | 100ug/mL COBALT 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCS-100PPM-125ML | 100ug/mL CESIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|--------------------|-----------------------------|------------|---------|-------|-----|
| MSCU-100PPM-125ML | 100ug/mL COPPER 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSFE-100PPM-125ML | 100ug/mL IRON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSGE-100PPM-125ML | 100ug/mL GERMANIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHF-100PPM-125ML | 100ug/mL HAFNIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHG-100PPM-125ML | 100ug/mL MERCURY HCl 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHGN-100PPM-125ML | 100ug/mL MERCURY HNO3 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSHO-100PPM-125ML | 100ug/mL HOLMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSIN-100PPM-125ML | 100ug/mL INDIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSK-100PPM-125ML | 100ug/mL POTASSIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSLI-100PPM-125ML | 100ug/mL LITHIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSMG-100PPM-125ML | 100ug/mL MAGNESIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSMN-100PPM-125ML | 100ug/mL MANGANESE 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSMO-100PPM-125ML | 100ug/mL MOLYBDENUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSNA-100PPM-125ML | 100ug/mL SODIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSNI-100PPM-125ML | 100ug/mL NICKEL 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSOS-100PPM-125ML | 100ug/mL OSMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSP-100PPM-125ML | 100ug/mL PHOSPHORUS 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSPB-100PPM-125ML | 100ug/mL LEAD 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSPT-100PPM-125ML | 100ug/mL PLATINUM 125mL | 38.22.0090 | 125 ml. | 27019 | 12% |
| MSRH-100PPM-125ML | 100ug/mL RHODIUM HCl 125mL | 38.22.0090 | 125 ml. | 34473 | 12% |
| MSRHN-100PPM-125ML | 100ug/mL RHODIUM HNO3 125mL | 38.22.0090 | 125 ml. | 34473 | 12% |
| MSS-100PPM-125ML | 100ug/mL SULFUR 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSB-100PPM-125ML | 100ug/mL ANTIMONY 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSC-100PPM-125ML | 100ug/mL SCANDIUM 125mL | 38.22.0090 | 125 ml. | 25622 | 12% |
| MSSE-100PPM-125ML | 100ug/mL SELENIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSI-100PPM-125ML | 100ug/mL SILICON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSN-100PPM-125ML | 100ug/mL TIN 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSSR-100PPM-125ML | 100ug/mL STRONTIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTB-100PPM-125ML | 100ug/mL TERBIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTEN-100PPM-125ML | 100ug/mL TELLURIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTI-100PPM-125ML | 100ug/mL TITANIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSTL-100PPM-125ML | 100ug/mL THALLIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSV-100PPM-125ML | 100ug/mL VANADIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---|-----------------------------|------------|---------|-------|-----|
| MSW-100PPM-125ML | 100ug/mL TUNGSTEN 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSY-100PPM-125ML | 100ug/mL YTTRIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSZN-100PPM-125ML | 100ug/mL ZINC 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| ICP & ICP-MS Single-element Standards - 1000 ppm | | | | | |
| CGAGI-30ML | 1000ug/mL SILVER 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGALI-30ML | 1000ug/mL ALUMINUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGALCLI-30ML | 1000ug/mL ALUMINUM/HCL 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGASI-30ML | 1000ug/mL ARSENIC 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGAUI-30ML | 1000ug/mL GOLD 30mL | 38.22.0090 | 30 ml. | 10399 | 12% |
| CGAUNI-30ML | 1000ug/mL GOLD/NIT 30mL | 38.22.0090 | 30 ml. | 11242 | 12% |
| CGBI-30ML | 1000ug/mL BORON 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGBAI-30ML | 1000ug/mL BARIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGBEI-30ML | 1000ug/mL BERYLLIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGBII-30ML | 1000ug/mL BISMUTH 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGCAI-30ML | 1000ug/mL CALCIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGCDI-30ML | 1000ug/mL CADMIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGCEI-30ML | 1000ug/mL CERIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGCOI-30ML | 1000ug/mL COBALT 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGCSI-30ML | 1000ug/mL CESIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGCUI-30ML | 1000ug/mL COPPER 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGDYI-30ML | 1000ug/mL DYSPROSIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGERI-30ML | 1000ug/mL ERBIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGEUI-30ML | 1000ug/mL EUROPIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGFEI-30ML | 1000ug/mL IRON 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGGAI-30ML | 1000ug/mL GALLIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGGDI-30ML | 1000ug/mL GADOLINIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGGEI-30ML | 1000ug/mL GERMANIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGHFI-30ML | 1000ug/mL HAFNIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGHGI-30ML | 1000ug/mL MERCURY 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGHOI-30ML | 1000ug/mL HOLMIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGICBRI-30ML | 1000ug/mL NH4 BROMIDE 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGICII-30ML | 1000ug/mL IODIDE 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGINI-30ML | 1000ug/mL INDIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---------------|--------------------------------|------------|---------|-------|-----|
| CGIRI-30ML | 1000ug/mL IRIIDIUM 30mL | 38.22.0090 | 30 ml. | 10961 | 12% |
| CGKI-30ML | 1000ug/mL POTASSIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGLAI-30ML | 1000ug/mL LANTHANUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGLII-30ML | 1000ug/mL LITHIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGLUI-30ML | 1000ug/mL LUTETIUM 30mL | 38.22.0090 | 30 ml. | 17144 | 12% |
| CGMGI-30ML | 1000ug/mL MAGNESIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGMNI-30ML | 1000ug/mL MANGANESE 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGMOI-30ML | 1000ug/mL MOLYBDENUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGNAI-30ML | 1000ug/mL SODIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGNBI-30ML | 1000ug/mL NIOBIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGNDI-30ML | 1000ug/mL NEODYMIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGNII-30ML | 1000ug/mL NICKEL 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGOSI-30ML | 1000ug/mL OSMIUM 30mL | 38.22.0090 | 30 ml. | 17144 | 12% |
| CGPI-30ML | 1000ug/mL PHOSPHORUS 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGPBI-30ML | 1000ug/mL LEAD 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGPDI-30ML | 1000ug/mL PALLADIUM 30mL | 38.22.0090 | 30 ml. | 10118 | 12% |
| CGPDNI-30ML | 1000ug/mL PALLADIUM HNO3 30mL | 38.22.0090 | 30 ml. | 19674 | 12% |
| CGPRI-30ML | 1000ug/mL PRASEODYMIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGPTI-30ML | 1000ug/mL PLATINUM 30mL | 38.22.0090 | 30 ml. | 16582 | 12% |
| CGPTNI-30ML | 1000ug/mL PLATINUM HNO3 30mL | 38.22.0090 | 30 ml. | 26981 | 12% |
| CGPTNO3I-30ML | 1000ug/mL CI FREE PLATINUM HNO | 38.22.0090 | 30 ml. | 26981 | 12% |
| CGRBI-30ML | 1000ug/mL RUBIDIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGREI-30ML | 1000ug/mL RHENIUM 30mL | 38.22.0090 | 30 ml. | 12366 | 12% |
| CGRHI-30ML | 1000ug/mL RHODIUM 30mL | 38.22.0090 | 30 ml. | 32321 | 12% |
| CGRHNI-30ML | 1000ug/mL RHODIUM HNO3 30mL | 38.22.0090 | 30 ml. | 55367 | 12% |
| CGRUI-30ML | 1000ug/mL RUTHENIUM 30mL | 38.22.0090 | 30 ml. | 9556 | 12% |
| CGSI-30ML | 1000ug/mL SULFUR 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGSBI-30ML | 1000ug/mL ANTIMONY 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGSCI-30ML | 1000ug/mL SCANDIUM 30mL | 38.22.0090 | 30 ml. | 12366 | 12% |
| CGSII-30ML | 1000ug/mL SILICON 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGSIOI-30ML | 1000ug/mL SILICA 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGSMI-30ML | 1000ug/mL SAMARIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGSNI-30ML | 1000ug/mL TIN 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---------------|-------------------------------|------------|---------|-------|-----|
| CGSNCLI-30ML | 1000ug/mL TIN/HCL 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGSRI-30ML | 1000ug/mL STRONTIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGTAI-30ML | 1000ug/mL TANTALUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGTBI-30ML | 1000ug/mL TERBIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGTEI-30ML | 1000ug/mL TELLURIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGTENI-30ML | 1000ug/mL TELLURIUM HNO3 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGTII-30ML | 1000ug/mL TITANIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGTLI-30ML | 1000ug/mL THALLIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGTMI-30ML | 1000ug/mL THULIUM 30mL | 38.22.0090 | 30 ml. | 8432 | 12% |
| CGVI-30ML | 1000ug/mL VANADIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGWI-30ML | 1000ug/mL TUNGSTEN 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGYI-30ML | 1000ug/mL YTTRIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGYBI-30ML | 1000ug/mL YTTERBIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGZNI-30ML | 1000ug/mL ZINC 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGZRI-30ML | 1000ug/mL ZIRCONIUM 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGAGI-125ML | 1000ug/mL SILVER 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGALI-125ML | 1000ug/mL ALUMINIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGALCLI-125ML | 1000ug/mL ALUMINUM/HCL 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGASI-125ML | 1000ug/mL ARSENIC 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGAUI-125ML | 1000ug/mL GOLD 125mL | 38.22.0090 | 125 ml. | 25049 | 12% |
| CGAUNI-125ML | 1000ug/mL GOLD/NIT 125mL | 38.22.0090 | 125 ml. | 25389 | 12% |
| CGBI-125ML | 1000ug/mL BORON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGBAI-125ML | 1000ug/mL BARIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGBEI-125ML | 1000ug/mL BERYLLIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGBII-125ML | 1000ug/mL BISMUTH 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCI-125ML | 1000ug/mL CARBON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCAI-125ML | 1000ug/mL CALCIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCDI-125ML | 1000ug/mL CADMIUM, 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCEI-125ML | 1000ug/mL CERIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCOI-125ML | 1000ug/mL COBALT 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCSI-125ML | 1000ug/mL CESIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGCUI-125ML | 1000ug/mL COPPER 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGDYI-125ML | 1000ug/mL DYSPROSIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|----------------|--------------------------------|------------|---------|-------|-----|
| CGER1-125ML | 1000ug/mL ERBIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGEU1-125ML | 1000ug/mL EUROPIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGFE1-125ML | 1000ug/mL IRON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGGA1-125ML | 1000ug/mL GALLIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGGD1-125ML | 1000ug/mL GADOLINIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGGE1-125ML | 1000ug/mL GERMANIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGHF1-125ML | 1000ug/mL HAFNIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGHG1-125ML | 1000ug/mL MERCURY 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGHO1-125ML | 1000ug/mL HOLMIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGICBR1-125ML | 1000ug/mL NH4 BROMIDE 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGICCL1-125ML | 1000ug/mL CHLORIDE/NH4CL 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGICII-125ML | 1000ug/mL IODIDE 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGINI-125ML | 1000ug/mL INDIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGIR1-125ML | 1000ug/mL IRIIDIUM, 125mL | 38.22.0090 | 125 ml. | 29371 | 12% |
| CGK1-125ML | 1000ug/mL POTASSIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGLA1-125ML | 1000ug/mL LANTHANUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGLI1-125ML | 1000ug/mL LITHIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGLU1-125ML | 1000ug/mL LUTETIUM 125mL | 38.22.0090 | 125 ml. | 27485 | 12% |
| CGMG1-125ML | 1000ug/mL MAGNESIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGMNI-125ML | 1000ug/mL MANGANESE 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGMO1-125ML | 1000ug/mL MOLYBDENUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGMSA1-125ML | 1000ug/mL SULFUR MSA 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGNA1-125ML | 1000ug/mL SODIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGNBI-125ML | 1000ug/mL NIOBIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGNB2O51-125ML | 1000ug/mL NIOBIUM, Ta FREE 125 | 38.22.0090 | 125 ml. | 14208 | 12% |
| CGND1-125ML | 1000ug/mL NEODYMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGNII-125ML | 1000ug/mL NICKEL 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGOS1-125ML | 1000ug/mL OSMIUM 125mL | 38.22.0090 | 125 ml. | 27485 | 12% |
| CGPI-125ML | 1000ug/mL PHOSPHORUS 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGPBI-125ML | 1000ug/mL LEAD 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGPD1-125ML | 1000ug/mL PALLADIUM 125mL | 38.22.0090 | 125 ml. | 30047 | 12% |
| CGPDNI-125ML | 1000ug/mL PALLADIUM HNO3 125mL | 38.22.0090 | 125 ml. | 62191 | 12% |
| CGPRI-125ML | 1000ug/mL PRASEODYMIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|----------------|--------------------------------|------------|---------|--------|-----|
| CGPTI-125ML | 1000ug/mL PLATINUM 125mL | 38.22.0090 | 125 ml. | 42625 | 12% |
| CGPTNI-125ML | 1000ug/mL PLATINUM HNO3 125mL | 38.22.0090 | 125 ml. | 77564 | 12% |
| CGPTNO3I-125ML | 1000ug/mL CI FREE PLATINUM HNO | 38.22.0090 | 125 ml. | 77564 | 12% |
| CGRBI-125ML | 1000ug/mL RUBIDIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGREI-125ML | 1000ug/mL RHENIUM 125mL | 38.22.0090 | 125 ml. | 20032 | 12% |
| CGRHI-125ML | 1000ug/mL RHODIUM 125mL | 38.22.0090 | 125 ml. | 96431 | 12% |
| CGRHNI-125ML | 1000ug/mL RHODIUM HNO3 125mL | 38.22.0090 | 125 ml. | 172830 | 12% |
| CGRUI-125ML | 1000ug/mL RUTHENIUM 125mL | 38.22.0090 | 125 ml. | 26107 | 12% |
| CGSI-125ML | 1000ug/mL SULFUR 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSBI-125ML | 1000ug/mL ANTIMONY 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSBFI-125ML | 1000ug/mL ANTIMONY HF 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSCI-125ML | 1000ug/mL SCANDIUM, 125mL | 38.22.0090 | 125 ml. | 20032 | 12% |
| CGSII-125ML | 1000ug/mL SILICON 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSINA1-125ML | 1000ug/mL SILICON/NaOH 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSIO1-125ML | 1000ug/mL SILICA 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSIONA1-125ML | 1000ug/mL SILICA/NaOH 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSMI-125ML | 1000ug/mL SAMARIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGSNI-125ML | 1000ug/mL TIN 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSNCL1-125ML | 1000ug/mL TIN/HCL 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSRI-125ML | 1000ug/mL STRONTIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGTA1-125ML | 1000ug/mL TANTALUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGTBI-125ML | 1000ug/mL TERBIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGTEI-125ML | 1000ug/mL TELLURIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGTENI-125ML | 1000ug/mL TELLURIUM HNO3 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGTII-125ML | 1000ug/mL TITANIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGTLI-125ML | 1000ug/mL THALLIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGTMI-125ML | 1000ug/mL THULIUM 125mL | 38.22.0090 | 125 ml. | 13277 | 12% |
| CGVI-125ML | 1000ug/mL VANADIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGWI-125ML | 1000ug/mL TUNGSTEN 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGWH2O1-125ML | 1000ug/mL TUNGSTEN IN H2O 125m | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGYI-125ML | 1000ug/mL YTTRIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGYBI-125ML | 1000ug/mL YTTERBIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGZNI-125ML | 1000ug/mL ZINC 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|--|-------------------------------|------------|---------|-------|-----|
| CGZRI-125ML | 1000ug/mL ZIRCONIUM 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| ICP & ICP-MS Single-element Standards - 10000 ppm | | | | | |
| CGAL10-30ML | 10,000ug/mL ALUMINUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGAS10-30ML | 10,000ug/mL ARSENIC 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGAU10-30ML | 10,000ug/mL GOLD 30mL | 38.22.0090 | 30 ml. | POR | 12% |
| CGBI10-30ML | 10,000ug/mL BISMUTH 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGCA10-30ML | 10,000ug/mL CALCIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGCE10-30ML | 10,000ug/mL CERIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGCO10-30ML | 10,000ug/mL COBALT 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGCR(3)10-30ML | 10,000ug/mL CHROMIUM(+3) 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGCU10-30ML | 10,000ug/mL COPPER 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGDY10-30ML | 10,000ug/mL DYSPROSIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGER10-30ML | 10,000ug/mL ERBIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGEU10-30ML | 10,000ug/mL EUROPIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGFE10-30ML | 10,000ug/mL IRON 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGGD10-30ML | 10,000ug/mL GADOLINIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGHO10-30ML | 10,000ug/mL HOLMIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGIR10-30ML | 10,000ug/mL IRIDIUM 30mL | 38.22.0090 | 30 ml. | POR | 12% |
| CGK10-30ML | 10,000ug/mL POTASSIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGLA10-30ML | 10,000ug/mL LANTHANUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGLI10-30ML | 10,000ug/mL LITHIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGLU10-30ML | 10,000ug/mL LUTETIUM 30mL | 38.22.0090 | 30 ml. | 35131 | 12% |
| CGMG10-30ML | 10,000ug/mL MAGNESIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGMN10-30ML | 10,000ug/mL MANGANESE 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGMO10-30ML | 10000ug/mL MOLYBDENUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGNA10-30ML | 10,000ug/mL SODIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGND10-30ML | 10,000ug/mL NEODYMIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGNI10-30ML | 10,000ug/mL NICKEL 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGP10-30ML | 10,000ug/mL PHOSPHORUS 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGPB10-30ML | 10,000ug/mL LEAD 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGPD10-30ML | 10,000ug/mL PALLADIUM 30mL | 38.22.0090 | 30 ml. | POR | 12% |
| CGPR10-30ML | 10,000ug/mL PRASEODYMIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGPT10-30ML | 10,000ug/mL PLATINUM 30mL | 38.22.0090 | 30 ml. | POR | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|-----------------|--------------------------------|------------|---------|-------|-----|
| CGRH10-30ML | 10,000ug/mL RHODIUM 30mL | 38.22.0090 | 30 ml. | POR | 12% |
| CGRU10-30ML | 10,000ug/mL RUTHENIUM 30mL | 38.22.0090 | 30 ml. | POR | 12% |
| CGS10-30ML | 10,000ug/mL SULFUR 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGSCI0-30ML | 10,000ug/mL SCANDIUM 30mL | 38.22.0090 | 30 ml. | 35131 | 12% |
| CGSE10-30ML | 10,000ug/mL SELENIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGS110-30ML | 10,000ug/mL SILICON 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGSM10-30ML | 10,000ug/mL SAMARIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGSN10-30ML | 10,000ug/mL TIN 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGTB10-30ML | 10,000ug/mL TERBIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGT110-30ML | 10,000ug/mL TITANIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGTM10-30ML | 10,000ug/mL THULIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGU10-30ML | 10,000ug/mL URANIUM 30mL | 38.22.0090 | 30 ml. | 15739 | 12% |
| CGV10-30ML | 10,000ug/mL VANADIUM 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGY10-30ML | 10,000ug/mL YTTRIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGYB10-30ML | 10,000ug/mL YTTERBIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGZN10-30ML | 10,000ug/mL ZINC 30mL | 38.22.0090 | 30 ml. | 13490 | 12% |
| CGZR10-30ML | 10,000ug/mL ZIRCONIUM 30mL | 38.22.0090 | 30 ml. | 14053 | 12% |
| CGAG10-125ML | 10,000ug/mL SILVER 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGAL10-125ML | 10,000ug/mL ALUMINUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGAS10-125ML | 10,000ug/mL ARSENIC 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGAU10-125ML | 10,000ug/mL GOLD 125mL | 38.22.0090 | 125 ml. | POR | 12% |
| CGB10-125ML | 10,000ug/mL BORON 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGBA10-125ML | 10,000ug/mL BARIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGBE10-125ML | 10,000ug/mL BERYLLIUM 125mL | 38.22.0090 | 125 ml. | 76166 | 12% |
| CGB110-125ML | 10,000ug/mL BISMUTH 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGC10-125ML | 10,000ug/mL CARBON 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGCA10-125ML | 10,000ug/mL CALCIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGCD10-125ML | 10,000ug/mL CADMIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGCE10-125ML | 10,000ug/mL CERIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGCO10-125ML | 10,000ug/mL COBALT 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGCR(3)10-125ML | 10,000ug/mL CHROMIUM(+3) 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGCS10-125ML | 10,000ug/mL CESIUM 125mL | 38.22.0090 | 125 ml. | 23525 | 12% |
| CGCU10-125ML | 10,000ug/mL COPPER 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|-----------------|--------------------------------|------------|---------|-------|-----|
| CGDY10-125ML | 10,000ug/mL DYSPROSIUM, 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGER10-125ML | 10,000ug/mL ERBIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGEU10-125ML | 10,000ug/mL EUROPIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGFE10-125ML | 10,000ug/mL IRON 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGGA10-125ML | 10,000ug/mL GALLIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGGD10-125ML | 10,000ug/mL GADOLINIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGGE10-125ML | 10,000ug/mL GERMANIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGHF10-125ML | 10,000ug/mL HAFNIUM 125mL | 38.22.0090 | 125 ml. | 52408 | 12% |
| CGHG10-125ML | 10,000ug/mL MERCURY 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGHO10-125ML | 10,000ug/mL HOLMIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGIN10-125ML | 10,000ug/mL INDIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGIR10-125ML | 10,000ug/mL IRIDIUM 125mL | 38.22.0090 | 125 ml. | POR | 12% |
| CGK10-125ML | 10,000ug/mL POTASSIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGLA10-125ML | 10,000ug/mL LANTHANUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGLI10-125ML | 10,000ug/ml LITHIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGLU10-125ML | 10,000ug/mL LUTETIUM 125mL | 38.22.0090 | 125 ml. | 52408 | 12% |
| CGMG10-125ML | 10,000ug/mL MAGNESIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGMN10-125ML | 10,000ug/mL MANGANESE 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGMO10-125ML | 10,000ug/mL MOLYBDENUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGMSA10-125ML | 10,000ug/mL SULFUR MSA 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGNA10-125ML | 10,000ug/mL SODIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGNB10-125ML | 10,000ug/mL NIOBIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGNB2O510-125ML | 10,000ug/mL NIOBIUM, Ta FREE 1 | 38.22.0090 | 125 ml. | 26321 | 12% |
| CGND10-125ML | 10,000ug/mL NEODYMIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGNI10-125ML | 10,000ug/mL NICKEL 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGPI10-125ML | 10,000ug/mL PHOSPHORUS 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGPB10-125ML | 10,000ug/mL LEAD 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGPD10-125ML | 10,000ug/mL PALLADIUM 125mL | 38.22.0090 | 125 ml. | POR | 12% |
| CGPR10-125ML | 10,000ug/mL PRASEODYMIUM125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGPT10-125ML | 10,000ug/mL PLATINUM 125mL | 38.22.0090 | 125 ml. | POR | 12% |
| CGRB10-125ML | 10,000ug/mL RUBIDIUM, 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGRE10-125ML | 10,000ug/mL RHENIUM 125mL | 38.22.0090 | 125 ml. | 52408 | 12% |
| CGRH10-125ML | 10,000ug/mL RHODIUM 125mL | 38.22.0090 | 125 ml. | POR | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|--|--------------------------------|------------|---------|-------|-----|
| CGRU10-125ML | 10,000ug/mL RUTHENIUM 125mL | 38.22.0090 | 125 ml. | POR | 12% |
| CGS10-125ML | 10,000ug/mL SULFUR 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGSB10-125ML | 10,000ug/mL ANTIMONY 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGSC10-125ML | 10,000ug/mL SCANDIUM 125mL | 38.22.0090 | 125 ml. | 52408 | 12% |
| CGSE10-125ML | 10,000ug/mL SELENIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGSI10-125ML | 10,000ug/mL SILICON 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGSM10-125ML | 10,000ug/mL SAMARIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGSN10-125ML | 10,000ug/mL TIN 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGSR10-125ML | 10,000ug/mL STRONTIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGTA10-125ML | 10,000ug/mL TANTALUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGTB10-125ML | 10,000ug/mL TERBIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGTE10-125ML | 10,000ug/mL TELLURIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGTI10-125ML | 10,000ug/mL TITANIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGTL10-125ML | 10,000ug/mL THALLIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGTM10-125ML | 10,000ug/mL THULIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGV10-125ML | 10,000ug/mL VANADIUM 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGW10-125ML | 10,000ug/mL TUNGSTEN 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGY10-125ML | 10,000ug/mL YTTRIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGYB10-125ML | 10,000ug/mL YTTERBIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGZN10-125ML | 10,000ug/mL ZINC 125mL | 38.22.0090 | 125 ml. | 21662 | 12% |
| CGZR10-125ML | 10,000ug/mL ZIRCONIUM 125mL | 38.22.0090 | 125 ml. | 23060 | 12% |
| CGZRCL10-125ML | 10,000ug/mL ZIRCONIUM/HCL 125m | 38.22.0090 | 125 ml. | 23060 | 12% |
| ICP & ICP-MS Isotopic Standards | | | | | |
| MS10B-10PPM-100ML | 10ug/mL 10 BORON 100mL | 38.22.0090 | 100 ml. | 28883 | 12% |
| MS11B-10PPM-100ML | 10ug/mL 11 BORON 100mL | 38.22.0090 | 100 ml. | 28883 | 12% |
| MS6LI-10PPM-125ML | 10ug/mL 6LITHIUM 125mL | 38.22.0090 | 125 ml. | 33541 | 12% |
| MS6LI-100PPM-125ML | 100ug/mL 6LITHIUM 125mL | 38.22.0090 | 125 ml. | 30280 | 12% |
| CG6LII-30ML | 1000ug/mL 6LITHIUM 30mL | 38.22.0090 | 30 ml. | 19111 | 12% |
| CG6LII-125ML | 1000ug/mL 6LITHIUM 125mL | 38.22.0090 | 125 ml. | 28184 | 12% |
| ICP & ICP-MS Speciation Standards | | | | | |
| MSCR(3)-10PPM-125ML | 10ug/mL CHROMIUM +3 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCR(6)-10PPM-125ML | 10ug/mL CHROMIUM +6 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| MSCR(3)-100PPM-125ML | 100ug/mL CHROMIUM +3 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---|---|------------|---------|-------|-----|
| MSCR(6)-100PPM-125ML | 100ug/mL CHROMIUM +6 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGAS(3)1-30ML | 1000ug/mL ARSENIC (+3) 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGAS(5)1-30ML | 1000ug/mL ARSENIC (+5) 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGCR(3)1-30ML | 1000ug/mL CHROMIUM(+3) 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGCR(6)1-30ML | 1000ug/mL CHROMIUM(+6) 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGSE(4)1-30ML | 1000ug/mL SELENIUM (+4) 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGSE(6)1-30ML | 1000ug/mL SELENIUM (+6) 30mL | 38.22.0090 | 30 ml. | 7588 | 12% |
| CGAS(3)1-125ML | 1000ug/mL ARSENIC (+3) 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGAS(5)1-125ML | 1000ug/mL ARSENIC (+5) 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCR(3)1-125ML | 1000ug/mL CHROMIUM(+3) 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGCR(6)1-125ML | 1000ug/mL CHROMIUM(+6) 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSE(4)1-125ML | 1000ug/mL SELENIUM (+4) 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| CGSE(6)1-125ML | 1000ug/mL SELENIUM (+6) 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| ICP & ICP-MS Multi-element Standards | | | | | |
| CCS-1-125ML | Rare Earths ICP-MS Standard | 38.22.0090 | 125 ml. | 45653 | 12% |
| CCS-2-125ML | Precious Metals ICP-MS Standard | 38.22.0090 | 125 ml. | 52641 | 12% |
| CCS-4-125ML | Alkali, Alkaline, Non-Transition ICP-MS | 38.22.0090 | 125 ml. | 42625 | 12% |
| CCS-5-125ML | Fluoride Soluble ICP-MS Standard | 38.22.0090 | 125 ml. | 42625 | 12% |
| CCS-6-125ML | Transition Elements ICP-MS Standard | 38.22.0090 | 125 ml. | 41461 | 12% |
| CIROS-OES-TS-125ML | CIROS ICP-OES Tuning Solution | 38.22.0090 | 125 ml. | 45886 | 12% |
| CMS-1-125ML | Rare Earth ICP-MS Standard | 38.22.0090 | 125 ml. | 34473 | 12% |
| CMS-2-125ML | Precious Metals ICP-MS Standard | 38.22.0090 | 125 ml. | 29349 | 12% |
| CMS-3-125ML | Fluoride Soluble ICP-MS Standard | 38.22.0090 | 125 ml. | 30280 | 12% |
| CMS-4-125ML | Hot-Plasma ICP-MS Complete Standard | 38.22.0090 | 125 ml. | 28650 | 12% |
| CMS-5-125ML | Cool-Plasma ICP-MS Complete Standard | 38.22.0090 | 125 ml. | 28650 | 12% |
| CSN-ISB-500ML | 1% Cesium Ionization Buffer for ICP | 38.22.0090 | 500 ml. | 28883 | 12% |
| CSN-ISB5-500ML | 5% Cesium Ionization Buffer for ICP | 38.22.0090 | 500 ml. | 92471 | 12% |
| GENESIS-ICAL-125ML | Genesis ICAL Calibration Standard | 38.22.0090 | 125 ml. | 30746 | 12% |
| IV-7-125ML | Quality Control Standard 7 | 38.22.0090 | 125 ml. | 30513 | 12% |
| IV-19-125ML | Quality Control Standard 19 | 38.22.0090 | 125 ml. | 61259 | 12% |
| IV-21-125ML | Quality Control Standard 21 | 38.22.0090 | 125 ml. | 67548 | 12% |
| IV-26-125ML | Quality Control Standard 26 | 38.22.0090 | 125 ml. | 74769 | 12% |
| IV-28-125ML | Quality Control Standard 28 | 38.22.0090 | 125 ml. | 79893 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---------------------|--|------------|---------|--------|-----|
| IV-ICPMS-71A-125ML | ICP-MS Complete Standard | 38.22.0090 | 125 ml. | 62890 | 12% |
| IV-ICPMS-71B-125ML | ICP-MS Refractory Elements Standard | 38.22.0090 | 125 ml. | 27019 | 12% |
| IV-ICPMS-71C-125ML | ICP-MS Precious Metals Standard | 38.22.0090 | 125 ml. | 54504 | 12% |
| IV-ICPMS-71D-125ML | ICP-MS Internal Standard | 38.22.0090 | 125 ml. | 39830 | 12% |
| IV-STOCK-1643-125ML | Trace Metals in Water | 38.22.0090 | 125 ml. | 89909 | 12% |
| IV-STOCK-2-125ML | ICP Calibration Standard 2 | 38.22.0090 | 125 ml. | 31911 | 12% |
| IV-STOCK-3-125ML | ICP Calibration Standard 3 | 38.22.0090 | 125 ml. | 21895 | 12% |
| IV-STOCK-4-125ML | ICP Calibration Standard 4 | 38.22.0090 | 125 ml. | 77098 | 12% |
| IV-STOCK-6-125ML | ICP-MS Calibration Standard 6 | 38.22.0090 | 125 ml. | 80592 | 12% |
| IV-STOCK-8-125ML | ICP Calibration Standard 8 | 38.22.0090 | 125 ml. | 63821 | 12% |
| IV-STOCK-9-125ML | ICP Calibration Standard - Toxic Element | 38.22.0090 | 125 ml. | 25622 | 12% |
| IV-STOCK-10-125ML | ICP Calibration Standard - Surface Water | 38.22.0090 | 125 ml. | 67082 | 12% |
| IV-STOCK-12-125ML | Stock Tuning Solution 12 | 38.22.0090 | 125 ml. | 31212 | 12% |
| IV-STOCK-13-125ML | ICP Calibration Standard - Trace Metals | 38.22.0090 | 125 ml. | 41228 | 12% |
| IV-STOCK-14-500ML | Wavelength Calibration standard | 38.22.0090 | 500 ml. | 72673 | 12% |
| IV-STOCK-17-125ML | ICP Calibration Standard - HCl Soluble Elements | 38.22.0090 | 125 ml. | 31445 | 12% |
| IV-STOCK-21-125ML | ICP-MS Calibration Standard 21 | 38.22.0090 | 125 ml. | 80126 | 12% |
| IV-STOCK-22-125ML | ICP-MS Calibration Standard 22 | 38.22.0090 | 125 ml. | 33774 | 12% |
| IV-STOCK-23-500ML | ICP-MS Calibration Standard 23 | 38.22.0090 | 500 ml. | 122286 | 12% |
| IV-STOCK-24-125ML | Stock Tuning Solution 24 | 38.22.0090 | 125 ml. | 42858 | 12% |
| IV-STOCK-26-125ML | ICP-MS Calibration Standard 26 | 38.22.0090 | 125 ml. | 34473 | 12% |
| IV-STOCK-27-125ML | ICP-MS Calibration Standard 27 | 38.22.0090 | 125 ml. | 74769 | 12% |
| IV-STOCK-28-125ML | ICP-MS Calibration Standard 28 | 38.22.0090 | 125 ml. | 65452 | 12% |
| IV-STOCK-29-125ML | ICP-MS Calibration Standard 29 | 38.22.0090 | 125 ml. | 38666 | 12% |
| IV-STOCK-30-125ML | ICP-MS Calibration Standard 30 | 38.22.0090 | 125 ml. | 29581 | 12% |
| IV-STOCK-31-125ML | ICP-MS Calibration Standard 31 | 38.22.0090 | 125 ml. | 39364 | 12% |
| IV-STOCK-33-125ML | ICP-MS Calibration Standard 33 | 38.22.0090 | 125 ml. | 25622 | 12% |
| IV-STOCK-34-125ML | ICP Calibration Standard 34 | 38.22.0090 | 125 ml. | 26553 | 12% |
| IV-STOCK-35-125ML | ICP Calibration Standard 35 | 38.22.0090 | 125 ml. | 23991 | 12% |
| IV-STOCK-36-125ML | ICP Calibration Standard 36 | 38.22.0090 | 125 ml. | 30513 | 12% |
| IV-STOCK-38-125ML | USP 232 Precious Metals Elemental Impurities | 38.22.0090 | 125 ml. | 52874 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|-----------------------|---|------------|---------|--------|-----|
| IV-STOCK-41-125ML | USP 232 Parenteral Elemental Impurities | 38.22.0090 | 125 ml. | 37035 | 12% |
| IV-STOCK-50-125ML | Environmental Calibration Standard | 38.22.0090 | 125 ml. | 71508 | 12% |
| IV-STOCK-51-125ML | 7500 Series PA Tuning Solution 1, 125 ml | 38.22.0090 | 125 ml. | 85251 | 12% |
| IV-STOCK-52-125ML | 7500 Series PA Tuning Solution 2, 125 mL | 38.22.0090 | 125 ml. | 40063 | 12% |
| IV-STOCK-53-125ML | Internal Standard | 38.22.0090 | 125 ml. | 35172 | 12% |
| IV-STOCK-54-125ML | Interference Check Standard | 38.22.0090 | 125 ml. | 28650 | 12% |
| IV-STOCK-55-125ML | Wavecal Standard, 125 mL | 38.22.0090 | 125 ml. | 26786 | 12% |
| IV-STOCK-56-125ML | ICP Calibration Standard, 125 mL | 38.22.0090 | 125 ml. | 20032 | 12% |
| IV-STOCK-57-125ML | ICP Calibration Standard, 125 mL | 38.22.0090 | 125 ml. | 17702 | 12% |
| IV-STOCK-58-125ML | ICP Calibration Standard | 38.22.0090 | 125 ml. | 65685 | 12% |
| IV-STOCK-60-125ML | USP 232 Drug Substance and Excipients | 38.22.0090 | 125 ml. | 40762 | 12% |
| IV-STOCK-72-125ML | 100 ppb Mercury for ICP/HNO ₃ | 38.22.0090 | 125 ml. | 29349 | 12% |
| IV-STOCK-73-125ML | 100 ppb Mercury for ICP/HCl | 38.22.0090 | 125 ml. | 29349 | 12% |
| IV-STOCK-74-500ML | ICP-MS Tuning solution | 38.22.0090 | 500 ml. | 76865 | 12% |
| IV-STOCK-75-125ML | ICP-MS Internal standard | 38.22.0090 | 125 ml. | 25622 | 12% |
| IV-STOCK-77-500ML | ICP-MS Tuning solution | 38.22.0090 | 500 ml. | 76865 | 12% |
| IV-STOCK-78-125ML | USP <232>/ICH Q3D(R2) Oral Elemental Impurities | 38.22.0090 | 125 ml. | 52874 | 12% |
| IV-STOCK-79-125ML | USP <232>/ICH Q3D(R2) Class 2B Oral Elemental Impurities | 38.22.0090 | 125 ml. | 52874 | 12% |
| PE-CHK-1-125ML | PE ICP-MS Check Standard 1 | 38.22.0090 | 125 ml. | 44023 | 12% |
| PE-TS-1-125ML | PE ICP-MS Tuning Solution | 38.22.0090 | 125 ml. | 45653 | 12% |
| PT-CO-500ML | Pt-Co (Platinum Cobalt) 500CU APHA/Hazen Color Standard | 38.22.0090 | 500 ml. | 64054 | 12% |
| THERMO-4AREV-500ML | THERMO TUNE B ICAP, 500 mL | 38.22.0090 | 500 ml. | 70110 | 12% |
| THERMO-4AREV-1L | THERMO TUNE B ICAP, 1L | 38.22.0090 | 1 Lit. | 98061 | 12% |
| THERMO-5A-125ML | THERMO ICAP Q | 38.22.0090 | 125 ml. | 52874 | 12% |
| THERMO-5A-250ML | THERMO ICAP Q | 38.22.0090 | 250 ml. | 83620 | 12% |
| THM-TS-1-125ML | THM ICP-MS Tuning Solution | 38.22.0090 | 125 ml. | 55203 | 12% |
| TUNE F-X-SERIES-125ML | TUNE F-X-SERIES-125ML | 38.22.0090 | 125 ml. | 110639 | 12% |
| VAR-CAL-1-125ML | VAR Calibration Standard 1 | 38.22.0090 | 125 ml. | 28650 | 12% |
| VAR-CAL-2-125ML | VAR Calibration Standard 2 | 38.22.0090 | 125 ml. | 72673 | 12% |
| VAR-CAL-7-125ML | VAR Calibration Standard 7 | 38.22.0090 | 125 ml. | 41927 | 12% |
| VAR-IS-1-125ML | VAR Internal Standard Solution | 38.22.0090 | 125 ml. | 42392 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|----------------------|--|------------|---------|-------|-----|
| VAR-TS-MS-125ML | VAR ICP-MS Tuning Solution | 38.22.0090 | 125 ml. | 42392 | 12% |
| WW-CAL-1A-125ML | Calibration Standard 1A | 38.22.0090 | 125 ml. | 41461 | 12% |
| WW-CAL-2-125ML | Calibration Standard 2 | 38.22.0090 | 125 ml. | 33075 | 12% |
| WW-CAL-3-125ML | Calibration Standard 3 | 38.22.0090 | 125 ml. | 27019 | 12% |
| WW-CAL-4A-125ML | Calibration Standard 4A | 38.22.0090 | 125 ml. | 25622 | 12% |
| WW-CAL-4B-125ML | Calibration Standard 4B | 38.22.0090 | 125 ml. | 13510 | 12% |
| WW-CAL-5-125ML | Calibration Standard 5 | 38.22.0090 | 125 ml. | 33774 | 12% |
| AGI-TS-1-125ML | AGI ICP-MS Tuning Solution 1 | 38.22.0090 | 125 ml. | 28650 | 12% |
| CLP-AES-CRQL-2-125ML | CRQL Standard 2 | 38.22.0090 | 125 ml. | 71741 | 12% |
| CLP-MS-SPK-125ML | Spike Standard | 38.22.0090 | 125 ml. | 58697 | 12% |
| CLP-MS-RINSE-125ML | 2% v/v HNO3 RINSE, 125mL | 38.22.0090 | 125 ml. | 16305 | 12% |
| CLP-MS-RINSE-500ML | 2% v/v HNO3 RINSE, 500mL | 38.22.0090 | 500 ml. | 36103 | 12% |
| CLPP-CAL-1-125ML | CLP CAL SOLUTION #1, 125mL | 38.22.0090 | 125 ml. | 89909 | 12% |
| CLPP-CAL-3-125ML | CLP CAL SOLUTION #3, 125mL | 38.22.0090 | 125 ml. | 30979 | 12% |
| CLPP-ICS-A-125ML | CLP INTERFER CHK SOLN A 125mL | 38.22.0090 | 125 ml. | 37501 | 12% |
| CLPP-ICS-B-125ML | CLP INTERFER CHK SOLN B, 125mL | 38.22.0090 | 125 ml. | 33774 | 12% |
| CLPP-ICS-B4-125ML | CLP INTERFER CHK SOLN B4,125mL | 38.22.0090 | 125 ml. | 29581 | 12% |
| CLPP-SPK-1-125ML | SOIL/WATER SPIKE SOLN 1, 125mL | 38.22.0090 | 125 ml. | 51244 | 12% |
| CLPP-SPK-2-125ML | SOIL/WATER SPIKE SOLN 2, 125mL | 38.22.0090 | 125 ml. | 11879 | 12% |
| 2007ICS-1-125ML | 200.7 Interference Check Standard 1 | 38.22.0090 | 125 ml. | 17935 | 12% |
| 2007ICS-3-125ML | 200.7 Interference Check Standard 3 | 38.22.0090 | 125 ml. | 42392 | 12% |
| 2007ICS-4-125ML | 200.7 Interference Check Standard 4 | 38.22.0090 | 125 ml. | 51476 | 12% |
| 2008CAL-1-125ML | 200.8 Calibration Standard 1 | 38.22.0090 | 125 ml. | 20032 | 12% |
| 2008CAL-2-125ML | 200.8 Calibration Standard 2 | 38.22.0090 | 125 ml. | 55203 | 12% |
| 2008ISS-125ML | 200.8 Internal Standard | 38.22.0090 | 125 ml. | 28184 | 12% |
| 2008TS-125ML | 200.8 Tuning Solution | 38.22.0090 | 125 ml. | 24690 | 12% |
| 6020CAL-1-125ML | 6020 Calibration Standard | 38.22.0090 | 125 ml. | 69645 | 12% |
| 6020ICS-0A-125ML | 6020 Rev-0 Interference Check Standard A | 38.22.0090 | 125 ml. | 61026 | 12% |
| 6020ICS-0B-125ML | 6020 Rev-0 Interference Check Standard B | 38.22.0090 | 125 ml. | 25156 | 12% |
| 6020ICS-8A-125ML | 6020 Ver-8 Interference Check Standard A | 38.22.0090 | 125 ml. | 65918 | 12% |
| 6020ICS-9A-125ML | 6020 Ver-9 Interference Check Standard A | 38.22.0090 | 125 ml. | 64054 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|-------------------|---|------------|---------|--------|-----|
| 6020ICS-9B-125ML | 6020 Ver-9 Interference Check Standard B | 38.22.0090 | 125 ml. | 25156 | 12% |
| 6020ISS-125ML | 6020 Internal Standard | 38.22.0090 | 125 ml. | 38200 | 12% |
| 6020SPK-W-125ML | 6020 Spike Standard - Water | 38.22.0090 | 125 ml. | 41228 | 12% |
| CLP-MS-SPK-125ML | Spike Standard | 38.22.0090 | 125 ml. | 58697 | 12% |
| CLPP-CAL-1-125ML | Calibration Standard 1 | 38.22.0090 | 125 ml. | 89909 | 12% |
| CLPP-CAL-3-125ML | Calibration Standard 3 | 38.22.0090 | 125 ml. | 30979 | 12% |
| CLPP-ICS-A-125ML | Interference Check Standard A | 38.22.0090 | 125 ml. | 37501 | 12% |
| CLPP-ICS-B-125ML | Interference Check Standard B | 38.22.0090 | 125 ml. | 33774 | 12% |
| CLPP-ICS-B4-125ML | Interference Check Standard B4 | 38.22.0090 | 125 ml. | 29581 | 12% |
| CLPP-SPK-1-125ML | Spike Standard 1 | 38.22.0090 | 125 ml. | 51244 | 12% |
| CLPP-SPK-2-125ML | Spike Standard 2 | 38.22.0090 | 125 ml. | 11879 | 12% |
| QCP-CICV-1-125ML | CICV Standard 1 | 38.22.0090 | 125 ml. | 75468 | 12% |
| QCP-CICV-3-125ML | CICV Standard 3 | 38.22.0090 | 125 ml. | 29814 | 12% |
| QCP-QCS-1-125ML | Quality Control Standard 1 | 38.22.0090 | 125 ml. | 85949 | 12% |
| QCP-QCS-2-125ML | Quality Control Standard 2 | 38.22.0090 | 125 ml. | 27252 | 12% |
| QCP-QCS-3-125ML | Quality Control Standard 3 | 38.22.0090 | 125 ml. | 61725 | 12% |
| QCP-QCS-4-125ML | Quality Control Standard 4 | 38.22.0090 | 125 ml. | 11879 | 12% |
| QCP-QCS-5-125ML | QC Standard / Instrument Performance Check | 38.22.0090 | 125 ml. | 31911 | 12% |
| WW-IPC-1-125ML | Instrument Performance Check 1 | 38.22.0090 | 125 ml. | 87347 | 12% |
| WW-IPC-2-125ML | Instrument Performance Check 2 | 38.22.0090 | 125 ml. | 30979 | 12% |
| WW-IPC-3-125ML | Instrument Performance Check 3 | 38.22.0090 | 125 ml. | 87347 | 12% |
| WW-LFS-1-125ML | Laboratory Fortified Stock Solution 1 | 38.22.0090 | 125 ml. | 87347 | 12% |
| WW-LFS-2-125ML | Laboratory Fortified Stock Solution 2 | 38.22.0090 | 125 ml. | 26321 | 12% |
| WW-MSCAL-1-125ML | Calibration Standard 3 | 38.22.0090 | 125 ml. | 10482 | 12% |
| WW-MSCAL-2-125ML | Calibration Standard 2A | 38.22.0090 | 125 ml. | 55203 | 12% |
| QCP-NUT-1-20ML | | 38.22.0090 | 20 ml. | 9240 | 12% |
| QCP-MIN-500ML | QCP MINERAL STANDARD, 20ML | 38.22.0090 | 500 ml. | 102795 | 12% |
| QCP-CN-20ML | QCP TOTAL CYANID STANDARD, 20ML | 38.22.0090 | 20 ml. | 20097 | 12% |
| QCP-HG-20ML | QC MERCURY STANDARD | 38.22.0090 | 20 ml. | 9702 | 12% |
| QCP-TRC-10ML | QCP Water standard Target range: 0.5-3.0mg/L; 10mL | 38.22.0090 | 10 ml. | 19404 | 12% |
| QCP-RAIN-125ML | QC RAINWATER STANDARD, 125mL | 38.22.0090 | 125 ml. | 56826 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|--|---|------------|---------|-------|-----|
| QCP-TMS-20ML | QC TRACE METALS STD, 20ml | 38.22.0090 | 20 ml. | 23331 | 12% |
| QCP-TURB-20ML | QC TURBIDITY STANDARD, 20ml | 38.22.0090 | 20 ml. | 10626 | 12% |
| QCP-OG-W-250ML | QCP Water standard Target weight range is 0.020-0.200g, 250ml | 38.22.0090 | 250 ml. | 39501 | 12% |
| QCP-SLD-450ML | QCP SOLID STANDARD, 450ML | 38.22.0090 | 450 ml. | 16632 | 12% |
| QCP-CAT-20ML | QC CATIONS STANDARD, 20ML | 38.22.0090 | 20 ml. | 16632 | 12% |
| QCP-CR6-20ML | QC CHROMIUM(+6) STANDARD, 20ML | 38.22.0090 | 20 ml. | 12012 | 12% |
| QCP-MTL-20ML | QCP METAL STANDARD, 20ML | 38.22.0090 | 20 ml. | 17556 | 12% |
| QCP-NT-20ML | QC NITRITE STANDARD, 20ML | 38.22.0090 | 20 ml. | 12936 | 12% |
| QCP-OG-A-20ML | QCP ACETONE STANDARD, 20ML | 38.22.0090 | 20 ml. | 27720 | 12% |
| QCP-WH-500ML | QCP WATER HARDNESS STANDARD, 500ML | 38.22.0090 | 500 ml. | 23562 | 12% |
| Ion Chromatography Anion & Cation Standards | | | | | |
| Ion Chromatography Anion Standards | | | | | |
| ICADPI-125ML | 1000µg/mL ADIPATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICBENI-125ML | 1000ug/mL BENZOATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICBRI-125ML | 1000ug/mL BROMIDE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICBRO3I-125ML | 1000ug/mL BROMATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICBTRI-125ML | 1000ug/mL BUTYRATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCITI-125ML | 1000ug/mL CITRATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCLI-125ML | 1000ug/mL CHLORIDE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCLO2I-125ML | 1000ug/mL CHLORITE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCLO3I-125ML | 1000ug/mL CHLORATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCLO4I-125ML | 1000ug/mL PERCHLORATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCO3I-125ML | 1000ug/mL CARBONATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCRO4I-125ML | 1000ug/mL CHROMATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICDEAI-125ML | 1000ug/mL DEA 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICDMAI-125ML | 1000ug/mL DIMETHYLAMINE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICFI-125ML | 1000ug/mL FLUORIDE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICGLYI-125ML | 1000ug/mL GLYCOLATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICGTRI-125ML | 1000µg/mL GLUTARATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICHCOI-125ML | 1000ug/mL FORMATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICII-125ML | 1000ug/mL IODIDE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICKHPI-125ML | 1000ug/mL PHTHALATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|--|--------------------------------|------------|---------|-------|-----|
| ICLCTI-125ML | 1000ug/mL LACTATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMEAI-125ML | 1000ug/mL MEA 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMLAI-125ML | 1000ug/mL MALATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMLEI-125ML | 1000ug/mL MALEATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMLOI-125ML | 1000ug/mL MALONATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMMAI-125ML | 1000ug/mL MMA 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMPAI-125ML | 1000ug/mL 3-METHOXYPROPYLAMINE | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMSAI-125ML | 1000ug/mL METHANE SULFONATE | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNNO2I-125ML | 1000ug/mL NITRITE AS N 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNNO3I-125ML | 1000ug/mL NITRATE AS N 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNO2I-125ML | 1000ug/mL NITRITE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNO3I-125ML | 1000ug/mL NITRATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNTAI-125ML | 1000ug/mL NITRILOTRIACETATE | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICOACI-125ML | 1000ug/mL ACETATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICOPRI-125ML | 1000ug/mL PROPIONATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICOXAI-125ML | 1000ug/mL OXALATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICPO4I-125ML | 1000ug/mL PHOSPHATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICPPO4I-125ML | 1000ug/mL PHOSPHATE AS P 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICS2O3I-125ML | 1000ug/mL THIOSULFATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICSCCI-125ML | 1000ug/mL SUCCINATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICSCNI-125ML | 1000ug/mL THIOCYANATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| IC-SCSI-125ML | 6 CATION CAL STD #1 125mL | 38.22.0090 | 125 ml. | 27951 | 12% |
| ICSO4I-125ML | 1000ug/mL SULFATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICTAI-125ML | 1000ug/mL TRIETHYLAMINE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICTEAI-125ML | 1000ug/mL TEA 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICTMAI-125ML | 1000ug/mL TMA 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICTMAHI-125ML | 1000ug/mL TMAH 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICTRTRI-125ML | 1000ug/mL TARTRATE 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| Ion Chromatography Cation Standards | | | | | |
| ICBAI-125ML | 1000ug/mL BARIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCAI-125ML | 1000ug/mL CALCIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICCSI-125ML | 1000ug/mL CESIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICKI-125ML | 1000ug/mL POTASSIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---|---|------------|---------|-------|-----|
| ICLI1-125ML | 1000ug/mL LITHIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICMGI-125ML | 1000ug/mL MAGNESIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNA1-125ML | 1000ug/mL SODIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNH41-125ML | 1000ug/mL AMMONIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICNNH41-125ML | 1000ug/mL AMMONIUM AS N 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICRB1-125ML | 1000ug/mL RUBIDIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| ICSRI-125ML | 1000ug/mL STRONTIUM 125mL | 38.22.0090 | 125 ml. | 13044 | 12% |
| Ion Chromatography Multi-ion Standards | | | | | |
| 300-CAL-A-125ML | 300.0 Calibration Standard | 38.22.0090 | 125 ml. | 25855 | 12% |
| 300-LFS-A-125ML | 300.0 Laboratory Fortification Solution | 38.22.0090 | 125 ml. | 41461 | 12% |
| IC-FAS-1A-125ML | Anion Calibration Standard | 38.22.0090 | 125 ml. | 27951 | 12% |
| IC-SCSI-125ML | 1000 ppm Cation Calibration Standard for IC | 38.22.0090 | 125 ml. | 27951 | 12% |
| IV-STOCK-59-125ML | Anion Calibration Standard, 125 mL | 38.22.0090 | 125 ml. | 27718 | 12% |
| IV-STOCK-7-125ML | Cation Calibration Standard 2 | 38.22.0090 | 125 ml. | 35638 | 12% |
| IV-STOCK-61-125ML | ANION MIX A | 38.22.0090 | 125 ml. | 19799 | 12% |
| IV-STOCK-62-125ML | CATION MIX B | 38.22.0090 | 125 ml. | 19799 | 12% |
| IV-STOCK-63-125ML | ANION MIX 4 | 38.22.0090 | 125 ml. | 19799 | 12% |
| IV-STOCK-64-125ML | ANION MIX 5 | 38.22.0090 | 125 ml. | 19799 | 12% |
| Ion Chromatography Eluent Concentrates | | | | | |
| BICARB-100ML | 0.5M BICARB ELUENT 100mL | 38.22.0090 | 500 ml. | 8851 | 12% |
| CARB-100ML | 0.5M NA CARB ELUENT 100mL | 38.22.0090 | 500 ml. | 8851 | 12% |
| ELUENT1817-500ML | .18M CARB/.17M BICARB 500mL | 38.22.0090 | 500 ml. | 31647 | 12% |
| ELUENT3510-500ML | .35M CARB/.10M BICARB 500mL | 38.22.0090 | 500 ml. | 31647 | 12% |
| ELUENT4514-500ML | 0.45M CARB / 0.14M BICARB | 38.22.0090 | 500 ml. | 41228 | 12% |
| MSAELUENT-100ML | 1.8M MSA 100mL | 38.22.0090 | 100 ml. | 8258 | 12% |
| MSAELUENT-500ML | 1.8M MSA 500mL | 38.22.0090 | 500 ml. | 31647 | 12% |
| ELUENT4514-500ML | 0.45M CARB / 0.14M BICARB | 38.22.0090 | 500 ml. | 41228 | 12% |
| Atomic Absorption Single-element Standards 1000ppm | | | | | |
| AAAG1-125ML | 1000ug/mL SILVER 125mL | 38.22.0090 | 125 ml. | 5982 | 12% |
| AAAL1-125ML | 1000ug/mL ALUMINUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAAS1-125ML | 1000ug/mL ARSENIC 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAAUI-125ML | 1000ug/mL GOLD 125mL | 38.22.0090 | 125ml. | 27053 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---------------|----------------------------|------------|---------|-------|-----|
| AABI-125ML | 1000ug/mL BORON 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AABA1-125ML | 1000ug/mL BARIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AABEI-125ML | 1000ug/mL BERYLLIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AABII-125ML | 1000ug/mL BISMUTH 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AACAI-125ML | 1000ug/mL CALCIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AACDI-125ML | 1000ug/mL CADMIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AACEI-125ML | 1000ug/mL CERIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AACOI-125ML | 1000ug/mL COBALT 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AACRI-125ML | 1000ug/mL CHROMIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AACSI-125ML | 1000ug/mL CESIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AACUI-125ML | 1000ug/mL COPPER 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AADY1-125ML | 1000ug/mL DYSPROSIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAERI-125ML | 1000ug/mL ERBIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAEUI-125ML | 1000ug/mL EUROPIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAFEI-125ML | 1000ug/mL IRON 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAGAI-125ML | 1000ug/mL GALLIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAGDI-125ML | 1000ug/mL GADOLINIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAGEI-125ML | 1000ug/mL GERMANIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAHFI-125ML | 1000ug/mL HAFNIUM 125mL | 38.22.0090 | 125ml. | 14353 | 12% |
| AAHG1-125ML | 1000ug/mL MERCURY 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAHO1-125ML | 1000ug/mL HOLMIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAINI-125ML | 1000ug/mL INDIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAIR1-125ML | 1000ug/mL IRIDIUM 125mL | 38.22.0090 | 125ml. | 31721 | 12% |
| AAK1-125ML | 1000ug/mL POTASSIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AALAI-125ML | 1000ug/mL LANTHANUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AALII-125ML | 1000ug/mL LITHIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AALUI-125ML | 1000ug/mL LUTETIUM 125mL | 38.22.0090 | 125ml. | 18155 | 12% |
| AAMG1-125ML | 1000ug/mL MAGNESIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAMNI-125ML | 1000ug/mL MANGANESE 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAMO1-125ML | 1000ug/mL MOLYBDENUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AANA1-125ML | 1000ug/mL SODIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AANBI-125ML | 1000ug/mL NIOBIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AANDI-125ML | 1000ug/mL NEODYMIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|----------------|-------------------------------|------------|---------|-------|-----|
| AANI1-125ML | 1000ug/mL NICKEL 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAP1-125ML | 1000ug/mL PHOSPHORUS 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAPBI-125ML | 1000ug/mL LEAD 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAPDI-125ML | 1000ug/mL PALLADIUM 125mL | 38.22.0090 | 125ml. | 27945 | 12% |
| AAPRI-125ML | 1000ug/mL PRASEODYMIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAPTI-125ML | 1000ug/mL PLATINUM 125mL | 38.22.0090 | 125ml. | 28196 | 12% |
| AARBI-125ML | 1000ug/mL RUBIDIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAREI-125ML | 1000ug/mL RHENIUM 125mL | 38.22.0090 | 125ml. | 15658 | 12% |
| AARHI-125ML | 1000ug/mL RHODIUM 125mL | 38.22.0090 | 125ml. | 79242 | 12% |
| AARUI-125ML | 1000ug/mL RUTHENIUM 125mL | 38.22.0090 | 125ml. | 28196 | 12% |
| AASI-125ML | 1000ug/mL SULFUR 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AASBI-125ML | 1000ug/mL ANTIMONY 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AASCI-125ML | 1000ug/mL SCANDIUM 125mL | 38.22.0090 | 125ml. | 13048 | 12% |
| AASEI-125ML | 1000ug/mL SELENIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AASII-125ML | 1000ug/mL SILICON 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AASMI-125ML | 1000ug/mL SAMARIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AASNI-125ML | 1000ug/mL TIN 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AASRI-125ML | 1000ug/mL STRONTIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AATAI-125ML | 1000ug/mL TANTALUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AATBI-125ML | 1000ug/mL TERBIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AATEI-125ML | 1000ug/mL TELLURIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AATII-125ML | 1000ug/mL TITANIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AATLI-125ML | 1000ug/mL THALLIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AATMI-125ML | 1000ug/mL THULIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAVI-125ML | 1000ug/mL VANADIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAWI-125ML | 1000ug/mL TUNGSTEN 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAYI-125ML | 1000ug/mL YTTRIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAYBI-125ML | 1000ug/mL YTTERBIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAZNI-125ML | 1000ug/mL ZINC 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| AAZRI-125ML | 1000ug/mL ZIRCONIUM 125mL | 38.22.0090 | 125ml. | 5982 | 12% |
| LACBI-125ML | 1% Lanthanum Releasing Agent | 38.22.0090 | 125ml. | POR | 12% |
| MM-MG-10-125ML | 1% Magnesium Nitrate Modifier | 38.22.0090 | 125ml. | 17935 | 12% |
| MM-P-40-125ML | 4% phosphate Modifier | 38.22.0090 | 125ml. | 11180 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|--|--|------------|---------|-------|-----|
| MM-PD-10-125ML | 1% Palladium Modifier | 38.22.0090 | 125ml. | POR | 12% |
| MM-PD-5-125ML | 0.5% Palladium Modifier | 38.22.0090 | 125ml. | POR | 12% |
| MM-PDMG-32-125ML | 0.3% Palladium/0.2% Magnesium Modifier | 38.22.0090 | 125ml. | POR | 12% |
| TCLP-IREV-125ML | TCLP Standard | 38.22.0090 | 125ml. | 30979 | 12% |
| TCLP-AA-HG-125ML | TCLP Mercury Standard | 38.22.0090 | 125ml. | 8385 | 12% |
| IV-STOCK-18-125ML | GFAA Calibration Standard | 38.22.0090 | 125ml. | 43790 | 12% |
| Potable & Waste Water Standards | | | | | |
| TOCKHPI-125ML | 1000 ppm Carbon for TOC made from KHP | 38.22.0090 | 125 ml. | 17004 | 12% |
| Conductivity standards | | | | | |
| CONI00000-25-125ML | 100,000 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI0000-25-125ML | 10,000 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI000-25-125ML | 1000 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI00-25-125ML | 100 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI0-25-125ML | 10 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI200-25-125ML | 1200 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI400-25-125ML | 1400 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI413-25-125ML | 1413 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI430-25-125ML | 1430 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI47-25-125ML | 147 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CON2-25-125ML | 2 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CON500-25-125ML | 500 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CON5-25-125ML | 5 CONDUCTIVITY, 125mL | 38.22.0090 | 125 ml. | 8851 | 12% |
| CON84-25-125ML | 84 CONDUCTIVITY, 125ML | 38.22.0090 | 125 ml. | 8851 | 12% |
| CONI00000-25-500ML | 100,000 CONDUCTIVITY, 500ML | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI0000-25-500ML | 10,000 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI000-25-500ML | 1000 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI00-25-500ML | 100 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI0-25-500ML | 10 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI200-25-500ML | 1200 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI400-25-500ML | 1400 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI413-25-500ML | 1413 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI430-25-500ML | 1430 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CONI47-25-500ML | 147 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |

| CATALOGUE NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|---------------------|---------------------------------|------------|---------|-------|-----|
| CON2-25-500ML | 2 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CON500-25-500ML | 500 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CON5-25-500ML | 5 CONDUCTIVITY, 500mL | 38.22.0090 | 500 ml. | 11879 | 12% |
| CON84-25-500ML | 84 CONDUCTIVITY, 500ML | 38.22.0090 | 500 ml. | 11879 | 12% |
| pH standards | | | | | |
| PH-1.68-500ML | pH 1.68 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 5948 | 12% |
| PH-10-500ML | pH 10 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 5948 | 12% |
| PH-11-500ML | pH 11 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 8250 | 12% |
| PH-12.47-500ML | pH 12.47 CALIBRATIONS STD, 500 | 38.22.0090 | 500 ml. | 8250 | 12% |
| PH-12-500ML | pH 12 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 7150 | 12% |
| PH-2-500ML | pH 2 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 7150 | 12% |
| PH-3-500ML | pH 3 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 7150 | 12% |
| PH-4-500ML | pH 4 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 5948 | 12% |
| PH-5-500ML | pH 5 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 7150 | 12% |
| PH-6.86-500ML | pH 6.86 CALIBRATIONS STD 500mL | 38.22.0090 | 500 ml. | 5948 | 12% |
| PH-6-500ML | pH 6 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 7150 | 12% |
| PH-7-500ML | pH 7 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 5948 | 12% |
| PH-8-500ML | pH 8 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 7150 | 12% |
| PH-9.18-500ML | pH 9.18 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 5948 | 12% |
| PH-9-500ML | pH 9 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 7150 | 12% |
| PH-10.01-500ML | pH 10.01 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 8250 | 12% |
| PH-12.45-500ML | pH 12.45 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 8250 | 12% |
| PH-4.01-500ML | pH 4.01 CALIBRATION STD, 500mL | 38.22.0090 | 500 ml. | 8250 | 12% |

New Products in 2024

| CATALOG NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|-------------------------|---|------------|----------|--------|-----|
| PCRM-IR-1000 | 1000 ug/g Ir solution standard in 10% v/v HCl 5x10ML AMPOULES | 3822.00.90 | 5x10 ml. | 300300 | 12% |
| PCRM-OS-1000 | 1000 ug/g Os solution standard in 15% v/v HCl 5x10ML AMPOULES | 3822.00.90 | 5x10 ml. | 300300 | 12% |
| DA-CL-COLOR-500ML | COMBINED COLOR REAGENT 500ML | 3822.00.90 | 500 ml. | 48510 | 12% |
| DA-NH3-DCI-500ML | DICHLOROISOCYANURATE (DCI) REAGENT 500ML | 3822.00.90 | 500 ml. | 48510 | 12% |
| QCP-CAT-20ML | QC CATIONS STANDARD, 20ML | 3822.00.90 | 20 ml. | 16632 | 12% |
| QCP-CR6-20ML | QC CHROMIUM(+6) STANDARD, 20ML | 3822.00.90 | 20 ml. | 12012 | 12% |
| QCP-HG-20ML | QC MERCURY STANDARD | 3822.00.90 | 20 ml. | 9702 | 12% |
| QCP-NT-20ML | QC NITRITE STANDARD, 20ML | 3822.00.90 | 20 ml. | 12936 | 12% |
| QCP-RAIN-125ML | QC RAINWATER STANDARD, 125ml | 3822.00.90 | 125 ml. | 56826 | 12% |
| QCP-TMS-20ML | QC TRACE METALS STD, 20ml | 3822.00.90 | 20 ml. | 23331 | 12% |
| QCP-TURB-20ML | QC TURBIDITY STANDARD, 20ml | 3822.00.90 | 20 ml. | 10626 | 12% |
| QCP-OG-A-20ML | QCP ACETONE STANDARD, 20ML | 3822.00.90 | 20 ml. | 27720 | 12% |
| QCP-MTL-20ML | QCP METAL STANDARD, 20ML | 3822.00.90 | 20 ml. | 17556 | 12% |
| QCP-MIN-500ML | QCP MINERAL STANDARD, 500ML | 3822.00.90 | 500 ml. | 102795 | 12% |
| QCP-SLD-450ML | QCP SOLID STANDARD, 450ML | 3822.00.90 | 450 ml. | 16632 | 12% |
| QCP-CN-20ML | QCP TOTAL CYANID STANDARD, 20ML | 3822.00.90 | 20 ml. | 20097 | 12% |
| QCP-WH-500ML | QCP WATER HARDNESS STANDARD, 500ML | 3822.00.90 | 500 ml. | 23562 | 12% |
| QCP-TRC-10ML | QCP Water standard Target range: 0.5-3.0mg/L; 10ml | 3822.00.90 | 10 ml. | 19404 | 12% |
| QCP-OG-W-250ML | QCP Water standard Target weight range is 0.020-0.200g, 250ml | 3822.00.90 | 250 ml. | 39501 | 12% |
| DA-NH3-SALICYLATE-500ML | SALICYLATE REAGENT 500ML | 3822.00.90 | 500 ml. | 48510 | 12% |
| DA-SO4-TURB-500ML | TURBIDIMETRIC REAGENT 500ML | 3822.00.90 | 500 ml. | 48510 | 12% |
| DA-PO4_ASCORBIC-500ML | WORKING ASCORBIC ACID 500ML | 3822.00.90 | 500 ml. | 48510 | 12% |
| DA-PO4-COLOR-500ML | WORKING PHOSPHATE COLOR REAGENT 500ML | 3822.00.90 | 500 ml. | 48510 | 12% |

Part-IV

Aquamicon's Product Range

From  Mitsubishi Chemical
Corporation

| CODE | DESCRIPTION | SPECIFICATION | HSN CODE | PACKING | ₹ | GST |
|---|---------------------------------|---|------------|--|-------|-----|
| AQUAMICRON Karl Fischer Reagents (Coulometric & Volumetric, Water standards) | | | | | | |
| COULOMETRIC REAGENTS | | | | | | |
| Anolyte Solution | | | | | | |
| APXAMIV | AQUAMICRON AXI | Moisture Max 0.15 mg H ₂ O/ml. | 3822.00.90 | 500 ml. | 13156 | 12% |
| Catholyte Solution | | | | | | |
| APCXUV | AQUAMICRON CXU | Moisture Max 0.6 mg H ₂ O/ml. | 3822.00.90 | 10x5 ml. ampule/case | 12298 | 12% |
| VOLUMETRIC REAGENTS | | | | | | |
| Aquamicon SS-Z Series | | | | | | |
| APSSZ10MV | Titrant SS-Z10M | Titre 0.7-1.2mg H ₂ O/ml. | 3822.00.90 | 500 ml. | 9743 | 12% |
| APSSZ30MV | Titrant SS-Z30M | Titre 2.5-3.5mg H ₂ O/ml. | 3822.00.90 | 500 ml. | POR | 12% |
| APSSZ50MV | Titrant SS-Z50M | Titre 4.5-5.5mg H ₂ O/ml. | 3822.00.90 | 500 ml. | 9743 | 12% |
| Solvent for Aldehydes & Ketones | | | | | | |
| APKTXV | Solvent KTX | Moisture Max. 0.5mg H ₂ O/ml. | 3822.00.90 | 500 ml. | 9743 | 12% |
| AQUAMICRON WATER STANDARDS | | | | | | |
| For Coulometric Titration | | | | | | |
| APAWS01V | AQUAMICRON Water Standard 0.1 | Titre 0.1± 0.01mg H ₂ O/gm | 2905.16.90 | 10x5 ml. ampule/case | 12355 | 18% |
| APAWS10V | AQUAMICRON Water Standard 1 | Titre 1.0± 0.05mg H ₂ O/gm | 2905.16.90 | 10x5 ml. ampule/case | 12412 | 18% |
| For Volumetric Titration | | | | | | |
| APAWS100V | AQUAMICRON Water Standard 10 | Titre 10.0 ± 0.5mg H ₂ O/gm | 2905.16.90 | 10x8 ml. ampule/case | 12729 | 18% |
| For Vaporization Method | | | | | | |
| APSWSV | AQUAMICRON Solid water Standard | Moisture 3.83 ± 0.1% | 2918.13** | 10 gm. | POR | 18% |
| Check Solution | | | | | | |
| APGCHPV | AQUAMICRON check solution P | Moisture 4.0± 2 mg H ₂ O/ml. | 2920.90.00 | 100 ml. glass bottle with Septum Cap | 11383 | 18% |

** HS Code

We declare that the items remain the same; it's a change in the product codes only. Please note the changes in the HSN Code and GST as well. An official declaration from MCC stating reasons is available.

| NEW CODE | OLD CODE | NEW HSN CODE | GST | OLD HSN CODE |
|-----------|----------|--------------|-----|--------------|
| APXAMIV | AXI | 3822.90.90 | 12% | 3822.00.90 |
| APCXUV | CXU | 3822.90.90 | 12% | |
| APSSZ10MV | SSZ10 | 3822.90.90 | 12% | |
| APSSZ30MV | SSZ30 | 3822.90.90 | 12% | |
| APSSZ50MV | SSZ50 | 3822.90.90 | 12% | |
| APKTXV | KTX | 3822.90.90 | 12% | |
| APAWS01V | AWS01 | 2905.16.90 | 18% | |
| APAWS10V | AWS10 | 2905.16.90 | 18% | |
| APAWS100V | AWS100 | 2905.16.90 | 18% | |
| APSWSV | SWS | 2918.13** | 18% | |
| APGCHPV | GCHP | 2920.90.00 | 18% | |

** HS Code

INCOMPATIBLE HAZARDOUS CHEMICALS – VIOLENT REACTIONS

Chemicals in left hand column should be stored under controlled conditions and these chemicals should not come in contact with the right hand column otherwise, a reaction will be violent and lead to an accident.

| Substance | Incompatible with |
|--|--|
| Acetic Acid | Chromic acid, Nitric acid, Peroxides, Permanganates |
| Acetic anhydride | Hydroxyl - containing compound, Ethylene glycol, Perchloric acid |
| Acetone | Concentrated nitric and sulphuric acid mixtures, Hydrogen peroxide |
| Acetylene | Chlorine, Bromine, Copper, Fluorine, Silver, Mercury |
| Alkali & alkaline earth metals such as Sodium Potassium, Lithium, Magnesium, Calcium | Carbon dioxide, Carbon tetrachloride and other chlorinated, hydrocarbons, halogens |
| Aluminium powder | Halogenated or oxygenated solvents |
| Ammonia (anhydrous) | Mercury, chlorine, calcium hypochlorite, iodine, bromine, hydrofluoric acid (anhydrous) |
| Ammonium nitrate | Acids, powdered metals, flammable liquids, chlorates, nitrates, sulfur, finely divided organic or combustible materials |
| Aniline | Nitric acid, hydrogen peroxide |
| Bromine | Ammonia, acetylene, butadiene, butane, methane, propane (or other petroleum gases), sodium carbide, benzene, finely divided metals, turpentine |
| Calcium oxide | Water |
| Carbon (activated) | Calcium hypochlorite, all oxidizing agents |
| Chlorates | Ammonium salts, Acids, metal powder, Phosphorous, Sulfur, finely divided organics or combustibles |
| Chromic acid | Acetic acid, naphthalene, camphor, glycerol, alcohol, chromium trioxide, flammable liquids in general |
| Chlorine | Ammonia, acetylene, butadiene, butane, methane, propane (or other petroleum gases), hydrogen, sodium carbide, benzene, finely divided metals, turpentine |
| Chlorine Dioxide | Ammonia, methane, phosphine, hydrogen sulphide |
| Copper | Acetylene, hydrogen peroxide |
| Fluorine | isolate from everything |
| Hydrazine | Hydrogen peroxide, Nitric acid, any other oxidant and heavy metal salts |
| Hydrocarbon (such as butane, benzene) | Fluorine, chlorine, bromine chromic acid, sodium peroxide |
| Hydrocyanic acid | Nitric acid, alkali |
| Hydrofluoric acid (anhydrous) | Ammonia (aqueous or anhydrous) |
| Hydrogen peroxides | Copper, chromium, iron, most metals or their salts, alcohols, acetone, organic materials, aniline, nitromethane, combustible materials |
| Hydrogen sulfide | Fuming nitric acid, oxidising gases |
| Iodine | Acetylene, ammonia (aqueous or anhydrous), hydrogen |
| Mercury | Acetylene, fulminic acids, ammonia, oxalic acid |
| Nitric acids | Acetic acid, aniline, chromic acid, hydrocyanic acid (concentrated) hydrogen sulfide, flammable liquids, flammable gases, copper, brass, any heavy metals |
| Nitromethane, lower nitroalkanes | Inorganic basis, amines, halogens |
| Oxalic acid | Silver, mercury |
| Oxygen | Oils, grease, hydrogen, flammable liquids, solids, gases |
| Perchloric acid | Acetic anhydride, bismuth & its alloys, alcohol, paper, wood, grease, oils |
| Peroxides, organic | Acids, (organic or mineral), avoid friction, store cold |
| Phosphinates | Any antioxidant |
| Phosphorus (white) | Air, oxygen alkalies, reducing agents |
| Potassium Chlorate | Sulfuric & other acids |
| Potassium perchlorate | (see sulfuric acids, other acids also chlorates) |
| Potassium permanganate | Glycerol, ethylene glycol, benzaldehyde, sulphuric acid |
| Silver | Acetylene, oxalic acids, tartaric acid, ammonium compounds, fulminic acid |
| Sodium | Carbon tetrachloride, carbon dioxide, water |
| Sodium Nitrate | Ammonium nitrate and other ammonium salts |
| Sodium Peroxide | Ethyl or methyl alcohol, glacial acetic acid, acetic anhydride, benzaldehyde, carbon disulphide, glycerine, ethylene glycol, ethyl acetate, methyl acetate, furfural |
| Sulphuric acid | Potassium chlorate, potassium perchlorate, potassium permanganate (similar compounds of light metals, such as sodium, lithium) |
| Thiocyanates | Metal nitrates, Nitrites, Oxidants |
| Trifluoromethane sulfonic acid | Perchlorate salts |

Solvent Properties & Miscibility Table

| SOLVENT POLARITY CHART | | Relative Polarity | | Compound Formula | | Group | | Representative Solvent Compounds | |
|------------------------|--|-----------------------|--|------------------|--|-----------------------|--|---|--|
| | | Nonpolar | | R - H | | Alkanes | | Petroleum ethers, Igitron, hexanes | |
| | | Ar - H | | R - O - R | | Aromatics | | Toluene, benzene | |
| | | R - X | | | | Ethers | | Diethyl ether | |
| | | R - COOR | | | | Alkyl halides | | Tetrachloroethane, chloroform | |
| | | R - CO - R | | | | Esters | | Ethylacetate | |
| | | R - NH | | | | Aldehydes and ketones | | Acetone, methyl ethyl ketone | |
| | | R - OH ¹ | | | | Amines | | Pyridine, triethylamine | |
| | | R - COOH ² | | | | Alcohols | | Methanol, ethanol, isopropanol, butanol | |
| | | H - OH | | | | Amides | | Dimethylformamide | |
| | | | | | | Carboxylic acids | | Ethanoic acid | |
| | | Water | | | | Water | | Water | |

| SOLVENT | POLARITY INDEX | REFRACTIVE INDEX @20°C | UV (nm) CUTOFF @1 AU | BOILING POINT (°C) | VISCOSITY (cPois) | SOLUBILITY IN WATER (% w/w) |
|-----------------------------------|----------------|------------------------|----------------------|--------------------|-------------------|-----------------------------|
| Acetic Acid | 6.2 | 1.372 | 230 | 118 | 1.26 | 100 |
| Acetone | 5.1 | 1.359 | 330 | 56 | 0.32 | 100 |
| Acetonitrile | 5.8 | 1.344 | 190 | 82 | 0.37 | 100 |
| Benzene | 2.7 | 1.501 | 280 | 80 | 0.65 | 0.18 |
| Butyl Acetate | 4.0 | 1.394 | 254 | 125 | 0.73 | 0.43 |
| n-Butanol | 3.9 | 1.399 | 215 | 118 | 2.98 | 7.81 |
| Carbon tetrachloride | 1.6 | 1.466 | 263 | 77 | 0.97 | 0.08 |
| Chloroform | 4.1 | 1.446 | 245 | 61 | 0.57 | 0.815 |
| Cyclohexane | 0.2 | 1.426 | 200 | 81 | 1.00 | 0.01 |
| 1,2-Dichloroethane ¹ | 3.5 | 1.424 | 225 | 84 | 0.79 | 0.81 |
| Dichloromethane ² | 3.1 | 1.424 | 235 | 41 | 0.44 | 1.6 |
| Dimethylformamide | 6.4 | 1.431 | 268 | 155 | 0.92 | 100 |
| Dimethyl Sulfoxide ³ | 7.2 | 1.478 | 268 | 189 | 2.00 | 100 |
| Dioxane | 4.8 | 1.422 | 215 | 101 | 1.54 | 100 |
| Ethyl Acetate | 4.4 | 1.372 | 260 | 77 | 0.45 | 8.7 |
| Ethanol | 5.2 | 1.360 | 210 | 78 | 1.20 | 100 |
| Diethyl Ether | 2.8 | 1.353 | 220 | 35 | 0.32 | 6.89 |
| Heptane | 0.0 | 1.387 | 200 | 98 | 0.39 | 0.0003 |
| Hexane | 0.0 | 1.375 | 200 | 69 | 0.33 | 0.001 |
| Methanol | 5.1 | 1.329 | 205 | 65 | 0.60 | 100 |
| Methyl-t-Butyl Ether ⁴ | 2.5 | 1.369 | 210 | 55 | 0.27 | 4.8 |
| Methyl Ethyl Ketone ⁵ | 4.7 | 1.379 | 329 | 80 | 0.45 | 24 |
| Pentane | 0.0 | 1.358 | 200 | 36 | 0.23 | 0.004 |
| n-Propanol | 4.0 | 1.384 | 210 | 97 | 2.27 | 100 |
| iso-Propanol ⁶ | 3.9 | 1.377 | 210 | 82 | 2.30 | 100 |
| di-iso-Propyl Ether | 2.2 | 1.368 | 220 | 68 | 0.37 | 100 |
| Tetrahydrofuran | 4.0 | 1.407 | 215 | 65 | 0.55 | 100 |
| Toluene | 2.4 | 1.496 | 285 | 111 | 0.59 | 0.051 |
| Trichloroethylene | 1.0 | 1.477 | 273 | 87 | 0.57 | 0.11 |
| Water | 9.0 | 1.333 | 200 | 100 | 1.00 | 100 |
| Xylene | 2.5 | 1.500 | 290 | 139 | 0.61 | 0.018 |

| SYNONYM TABLE | |
|-----------------------------|--|
| 1. Ethylene Chloride | |
| 2. Methylene Chloride | |
| 3. Methyl Sulfoxide | |
| 4. Ter-t-Butyl Methyl Ether | |
| 5. 2-Butanone | |
| 6. 2-Propanol | |

Immiscible
 Miscible

Immiscible means that in some proportions two phases will be produced.

Products – Till Stock Lasts

| SAP CODE | CAS NO. | PRODUCT DESCRIPTION | HSN CODE | PACKING | ₹ | GST |
|------------|-----------|------------------------------|----------|----------|------|-----|
| 10301SG100 | 592-10-4 | Bismuth carbonate extrapure* | 29161920 | 100 gm.* | 1398 | 18% |
| 10301SG500 | | | | 500 gm.* | 6219 | 18% |
| 10306SG250 | 7787-68-0 | Bismuth sulphate* | 28332990 | 250 gm.* | 5054 | 18% |
| 10306SG500 | | | | 500 gm.* | 9072 | 18% |



The Partner of Choice

Our infrastructure and easy to use website helps you to select the right product at the right time for your desired application along with comprehensive technical information on each product

www.actylis.com

Add to your favourites today

- Product search backed by a **powerful search engine**
- **Search** using partial names
- **COA & MSDS** available online
- **Comprehensive information** for enlisted products
- **Download** standard technical literature to aid your functions

